





BAHZ-401-A & BAHZ-402-A

Horizontal Autoclave

Thank you for Choosing Biolab products. Please read the "Operating Instructions" and "Warranty" before operating this unit to assure proper operation.



Index

1.	Precaution on Safe Operation	03
2.	About the Instrument	05
3.	Installation of Autoclave	08
4.	Maintain & Management	40
5.	Troubleshooting	45
6.	Electrical Wiring Diagram and Component Distribution Diagram	49
7.	Pipeline Diagram	52



01 Precaution on Safe Operation

• To ensure the safe and correct use of the instrument, please read this manual carefully before use and operate according to the instructions in the manual. If it is not used according to the method specified by the manufacturer, the protection provided by the instrument may be damaged.

• In addition to sterilization, drying and agar melting, the instrument shall not be used for other purposes; it shall not be used for sterilization of inflammable, explosive, oxide prone or strong acid, alkali, salt water and other substances, or it may cause corrosion of sterilization chamber and pipeline, or even explosion.

• During installation, it is required to connect correctly according to the power requirements on the instrument nameplate; if the voltage fluctuates too much, it is required to use a regulated power supply to ensure the best performance of the instrument; if other types of voltage are used, it is required to use a transformer, otherwise the instrument will be damaged.

• The instrument must be reliably grounded. Do not connect the ground wire of the instrument to the plastic pipe, gas pipe, telephone ground wire, lightning rod, etc.

• Do not let the object block the exhaust port on the safety valve, so as to avoid that the safety valve can not exhaust and relieve the pressure in case of abnormal situation.

• Before opening the chamber cover, make sure that the reading of the pressure gauge is "0 MPa"; when the pressure in the sterilization chamber is higher than "0 MPa", do not open the chamber cover and drain valve, otherwise it will cause high-pressure steam to spray out and hurt people.

• When adding distilled water into the sterilization chamber, do not leak the water into the control circuit, so as to avoid electric shock accident or other faults.

• When using cleaning or other bags, please put the bag in the stainless steel basket first, and then put it into the sterilization chamber, otherwise the accuracy of temperature may be affected.

• Pay attention to observe the temperature in the sterilization chamber. The temperature is high at the end of operation. When opening the cover, pay attention not to put your face and hands close to the sterilization chamber to prevent scalding caused by steam spraying. When taking out the articles from the sterilization chamber, wear heat insulation gloves. Since the liquid needs to be cooled for a certain time, when the sterilized liquid material is taken out from the sterilization chamber, it is necessary to confirm that the temperature has dropped to a sufficiently low level to avoid scalding.

• Distilled water must be used as sterilization water to avoid affecting the service life of

sterilizer. When the instrument works continuously, it shall be ensured that there is an interval of more than 15 minutes for the instrument to cool down. Otherwise, the instrument will not be able to produce enough saturated steam.

• In case of any abnormal situation (such as abnormal sound, smell, smoke), turn off the power supply immediately, pay attention to observation, and contact the local dealer or our after-sales service department after the abnormal situation no longer continues.

• It is recommended to place a pressure steam sterilization chemical prompt card (hereinafter referred to as the chemical indicator card) on the sterilized substance for each sterilization. After a sterilization cycle, when the color change of the chemical indicator card coincides with the temperature and temperature duration to be represented, it indicates that the temperature and temperature duration reached have met the requirements of the sterilization Institute If necessary, sterilization can be carried out; otherwise, sterilization requirements are not met.

• If the instrument is used in China, it shall comply with the relevant provisions of "Regulations on safety supervision of special equipment" and "TSG 21-2016 Safety Technical Supervision Specification for stationary pressure vessels".

• Precautions for instrument operation

• The user shall read and understand the contents specified in the operation manual before the first use;

• The user must receive correct and complete training related to work safety and accident prevention or complete training on operation guidance of autoclave;

- Users should follow the safety, maintenance and maintenance regulations in the manual;
- Explanation of relevant marks in the manual

igsquircless This symbol indicates the need to pay attention to safety signs during operation.



02 About the Instrument

1. Application

• This series of products are used for sterilization of scientific research institutions, laboratory utensils, culture media and unsealed liquids or preparations.

2. Types of microorganisms killed

• Using hot and humid high pressure steam as sterilization factor to kill loaded microorganisms, including spores of bacteria, spores of fungi, etc.

3. Normal working conditions

- Ambient temperature: 5 $^{\circ}$ C ~ 40 $^{\circ}$ C
- relative humidity not greater than 85%.
- Atmospheric pressure: 70kPa~106kPa.
- Suitable for power supply AC380V \pm 38V, (50 60) Hz \pm 1Hz.

4. Transportation requirement

Instruments are not allowed to stand upside down, overlap, below is not allowed to put items, avoid rain, carefully handle, there should be anti-movement measures.

5. Storage requirement

- Ambient temperature:-20 °C ~ 55 °C.
- Relative humidity not greater than 93%.
- An indoor or sheltered place free of corrosive gas and well ventilated

6. Service life

Production date:

Service life: 8 years. It is not allowed to use it after the service life, so as to avoid danger.

7. Technical specifications

Model	BAHZ-401-A	BAHZ-402-A	
Capacity(L)	180L	280L	
Dimension(L*W*H,mm)	1450*800*1500	1450*800*1500	

Chamber dimension(Dia*H,mm)	φ508*1130	φ608*1130	
Net weight	650kg	720kg	
Total power	16k	×W	
Chamber material	S30-	408	
Sterilization temperature	115°C ~	• 135°C	
Rated working pressure	0.25	MPa	
Display range of pressure gauge	-0.1 ~ 0.5MPa		
Set pressure of safety valve	0.29	MPa	
Sterilization time range	1min ~ 6	000min	
Preset range of drying time	0min ~ 300min		
Pulsation times	1-99times		
Cooling time	10min ~ 20min		
controller	Programmable intelligent logic system		
Cavity door opening mode	Automatic door		
Operation Mode	 Fabric mode: standby - preheating sterilization - Exhaust - drying - keep Device mode: standby - preheating sterilization - Exhaust - drying - keep Liquid mode: standby - preheating sterilization - Exhaust - balance preservation exhaust balance preservation exhaust balance preservation exhaust balance preservation exhaust balance preservation - end. Self-defined program: Vacuum sterilization procedure pulsation - heating - sterilization preserve - end Liquid program: standby - preserve - end Special lower arrangement procedure preserve - prese	ng - pulsation - heating - balance pressure - end ng - pulsation - heating - balance pressure - end ng - cooling - heating - bressure - end r preheating pulsation heating sure end cuum - Leak Detection - e: standby - preheating - - Exhaust - drying - balance heating - cooling - heating - bressure - end bcedure: standby - preheating -	

Safety device	Pressure safety interlock, safety valve, over-current and short- circuit protection, automatic fault detection, dry burning protection, water level detection, over temperature protection and leakage protection		
Standard spare parts	Stainless steel baskets, storage board		
Optional spare parts	Article thermometer, stainless steel cover, printer, water pressure sensor (optional), back pressure system		

Note 1: the material of sterilization chamber can be purchased s31603, and the model is equipped with - L, for example: BAHZ-401-A becomes BAHZ-401-A-I;

8. Introduction to the main structure of the product

It is mainly composed of shell (shell decoration cover), sterilization chamber, sterilization chamber door, built-in steam generator, pipeline system, temperature control system, pressure detection, safety interlock device and pulsating vacuum pumping device. It is summarized as follows:

•Sterilization chamber container part:

The design, manufacture and acceptance of pressure vessels shall be carried out in accordance with GB / t150-2011 pressure vessels, and shall be supervised by tsg21-2016 Supervision Regulation on safety technology of stationary pressure vessels. It is mainly composed of cylinder body and door body. The cylinder body forms a limited sterilization chamber, which is composed of an inner cavity, a front sealing plate and a container fixing frame;

•Steel structure part:

The utility model is composed of a chassis and a frame, the chassis supports the main body of the container, and the frame connects the outer decorative cover plate;

•Exterior decoration part (cover plate):

It is made of 304 steel plate or cold rail plate by spraying and pressing, which can not only decorate the equipment, but also prevent the influence of thermal radiation on the environment;

•Control system:

It is composed of programmable intelligent controller (PLC) / intermediate relay, pressure transmitter, temperature transmitter, travel switch and other auxiliary devices. PLC processes temperature, pressure, some switches and other input information, and outputs different

control signals to automatically complete the control process;

•Built in steam generator unit:

The electric heating generator is a device that can automatically generate steam, and has the functions of automatic control of water inlet, heating, dry burning prevention and overpressure protection;

• Piping system:

It is composed of solenoid valve, filter, check valve, vacuum pump, booster pump, pressure gauge, steam filter, air filter and other pipe fittings,

According to its function, it is mainly divided into the following branches:

- 1. Steam inlet branch 2. cooling water branch 3.water inlet branch
- 4. drainage branch 5. Booster pump water inlet branch 6. evaporator water inlet branch
- 7. drainage steam exhaust branch 8. steam outlet branch
- 9. Inner chamber steam inlet branch 10. evaporator drainage branch
- 11.general exhaust branch

Safety interlock

As the protective device of the instrument in the sterilization process, the interlock device will be triggered as long as the instrument is in the working state, and the interlocking status will be displayed on the touch screen. When the interlock fails, click "start program" to run the program.

9.Description of the working principle of the product

•Principle description:

The steam sterilizer can be used as a combination of the temperature of the steam sterilizer and the temperature of the steam sterilizer after the sterilization.

03 Installation of Autoclave

1.Placement of Autoclave

Lift the instrument to the installation site. For the convenience of heat dissipation and maintenance of the equipment, it is necessary to follow the installation drawing (as shown in the figure below). It can be fixed after determining the orientation. If forklift is used, the forklift should be between two sets of foot cups.



1.Preparation before installation

•According to the installation drawing provided by our company, under the guidance of professionals, professional construction personnel are responsible for it;

•Location and environmental requirements:

1) According to the size of the instrument's installation drawing, reserve the installation space of the instrument to meet the convenience of instrument installation, maintenance and use;

2) The selected sterilizer space needs good ventilation and heat dissipation;

3) Avoid installation in heavy dust, oil mist containing conductive particles, corrosive gas,

combustible gas environment.

4) Avoid installation where electric shock or vibration may occur.

5) Avoid installation in places with high temperature and humidity or easy to be wet by rain.

6) Avoid installation in strong magnetic field.

• Energy requirements:

A. Water source inspection:

Tap water source requirements: pressure 0.15 ~ 0.3MPa; Suggestions:

•If the water source pressure exceeds 0.3MPa, a pressure reducing valve shall be installed; otherwise, if the water source pressure is lower than 0.15Mpa, a booster pump shall be installed.

•It is suggested to install a water source pressure gauge for the tap water source, so as to monitor the pressure of the water source and ensure the normal use of the instrument. Requirements for source of demineralized water: demineralized water / softened water / pure water or ultra pure water with pressure of 0-0.3 MPa.

Suggestion: if the barreled demineralised water (without pressure water source), the outlet pipe should be installed at the bottom of the barrel, and a minimum water level should be drawn above the water outlet of the barrel to ensure that there is water at the minimum water level. At the same time, for the first start-up, the demineralization water pipe (including booster pump head) should be filled with water without air.

B. Power supply check:

1) It is required to install a power switch box on the right or right rear wall of the sterilizer. A three-phase air switch and an over-current protection device must be installed in the switch box to realize the function of power on-off and over-current protection. Make sure the switch box is close to the equipment and accessible to the operator. The switch box should be marked as special for sterilizer. In order to ensure the safety of personnel and equipment, it is necessary to ensure the reliable connection between the ground wire of the equipment and the ground wire in the power switch box. The pipelines and lines should be horizontal and vertical and effectively fixed.

2) The power supply is 380V 50 / 60Hz 40A three-phase AC. according to the relevant international electrical equipment regulations, the conductor color is yellow for phase a (u.l1), green for phase B (v.l2), red for phase C (w.l3), light blue for neutral line (n), yellow and green for grounding wire PE.

2. Arrival inspection instructions

2.1 after receiving the goods, the customer should check the appearance of the packing box immediately. If there is any damage, take a picture and mark it on the delivery note. Contact the distributor, inform the manufacturer and prepare for claim.

2.2 after the outer package of sterilizer is opened, check whether there is mechanical damage, remember the packaging method and keep the packaging materials until the instrument is accepted. The mechanical inspection includes whether there are traces of physical damage,

such as whether there are scratches on the surface of the panel, whether the knob is broken, etc. if you find any problems, please contact our company immediately.

3.Instrument installation

•After the instrument is hoisted or transported in place, adjust the height of the foot cup to make the instrument adjust to the horizontal position (based on the storage plate inside the container) and make the sealing door open and close freely.

When adjusting, prepare tools, such as a 300 mm spanner, a small jack (bearing 1000kg), and a stainless steel plate (thickness of 1 mm / 2 mm) of 80 * 80 size. The instrument cannot be fixed by casters.

•The surface of instrument installation base should be solid and flat. If installed on the second floor of the user, consider whether the floor strength is reasonable or not, even strengthen the floor.

•Nozzle installation:

According to the mark of the drain valve cover, connect the inlet branch of vacuum pump (as shown in the right figure of tap water), the main discharge branch (main drain) of booster pump inlet branch (as shown in the right figure of demineral water mark) and the evaporator drainage branch



(sewage) marks. The distance between the equipment drain outlet and the connecting pipe of drainage pipe shall not be more than 0.5m, and the drainage pipe shall not have reverse pressure.

Note: Chapter 7 piping diagram for reference.

•Pipe size:

Sewage, distilled water, tap water connection 1 / 2 joint, the total discharge is 1 inch joint.

•The power supply is installed according to the electrical wiring diagram.

Remarks: Chapter 6 electrical wiring diagram for reference.

Instrument cleaning

•Remove the protective package of the power plug and connect the power supply.

•Open the sterilization chamber door and remove the foam from the container. Clean the container and put it into the storage plate and stainless steel basket in turn.

Be careful to clean up the foam and other miscellaneous items in the sterilization cavity so as not to jam the pipes.

Instrument debugging

Some programs have been set in advance before the instrument leaves the factory. In the process of equipment debugging, users (technicians) can set the sterilization parameters according to the use of sterilization materials.

Before debugging, check whether the electrical parts are loose and fall off; whether the pipeline connection is reasonable; whether the main switch of the circuit breaker is damaged;
Open the cut-off valve of distilled water and tap water, close the blow down stop valve, turn on the power supply, and observe whether the pressure gauge meets the use requirements;
Adjust the rotation direction of vacuum pump:

Specific operation: turn on the power switch, click the "open door" button on the display screen to check whether the pressure drops, and the value becomes smaller (atmospheric pressure 100KPA). If the pressure value does not change, it indicates that the power line is connected in reverse; if the pressure value will drop, the wiring is correct.

03 Operation Instruction

I. Basic Sterilization Illustration & Operation Instruction



Turn off the power

🔼 It is strictly forbidden for non operators to touch the control screen at will

1. Turn on the power supply and water source switch

Turn on the water source and main power switch of the equipment (first close the small blue handle, then the 4P handle), and then turn on the key switch and emergency stop switch on the left side of the pressure gauge. After a period of self inspection, the instrument enters the self inspection interface, as shown in the following figure: If there is a fault display in the self-test completion screen, please remove the fault according to Chapter 5 "fault guide" and click "confirm" to exit the fault.

Note: If you want to enter the training interface, you need to enter the training interface after power on self-test. After the training interface is started, the steam generator will stop heating (method: click "Administrator" to enter the administrator interface after entering the number and password, and the administrator's authority can enter the training interface)

Evaporator Pressure 98.0 kPa Object Temperature 30.5 °C
Water Pressure 0.0 kPa safety valve 353 Day
Door Status Closed Lock Status Unlocked

2.Open the chamber door of sterilizer

2.1.In the system login interface, click the "operator" in the next line and enter the corresponding operator number and password according to the prompt to enter the operator interface (the operator number and password are provided by the administrator, and the password and authority can be set in the administrator management interface, and the password is composed of at least 6 digits and is composed of numbers)



Note: When the machine leaves the factory, the administrator number is set as 15; the password of the administrator is 1234567; the password of the operator (No. 1-10) and the technician (No. 11-14) is set by the administrator, and the user without password is the user who is not enabled.

2.2. After entering the number and password, the door operation interface as shown in the figure below will pop up under the correct condition:



If you input the wrong password, a dialog box will pop up to remind you that the password is wrong. Please re-enter it, as shown in the following figure:



Press OK to return to the operation login interface.

2.3. Click "open door" to display the following screen:

After the interlock is opened, the following interface appears:

Inner Temp. 31. 9	℃ Inner Pressure	99.3 kPa	Object Temp. 30. 5	r
Door Status Close	Door Status Unlocked	Negative Pressure Release	Open Door	
Automatic openin	g, please wait			
Door Operation	Program	Man	agement	Exit

When the inner chamber pressure returns to atmospheric pressure, the chamber door will automatically pop open, indicating the following interface:

Inner Temp. 31.	9 ℃ Inner Pressur	e 99.3 kPa	Object Temp. 30. 5	C
Door Status Close	Door Status Unlocked	Negative Pressure Release	Open Door	
Door Operation	Program	Manager	nent .	Exit

Note: If the system can not automatically restore to atmospheric pressure, please click "negative pressure release" to manually restore. The system shows absolute pressure

3.Loading of articles to be sterilized

1> Open the chamber door, take out the stainless steel basket, and put the sterilized material into the basket

2> Put the basket into the sterilization chamber again after placing the sterilized materials3> Preparation of instrument sterilization

The equipment must be cleaned and then put into the sterilizer to avoid the harm of residual substances on the equipment to the sterilizer and the equipment itself. For example, blood and other impurities.

The following cleaning plan is for your reference:

For instruments that need to be sterilized after use, the residues attached to the instruments should be cleaned in time. Cleaning machine, cleaning agent and distilled water are recommended. After cleaning, it is recommended to rinse with water again to ensure its cleanness.

When putting the instrument into the sterilizer storage board, please pay attention to: Instruments should be evenly arranged and spaced without overlapping. Otherwise, it will cause insufficient sterilization. Different types of instruments should be placed in different storage boards, such as stainless steel, carbon steel, etc. If carbon steel instruments are stored in the storage board, the storage board should be padded with several layers of disinfection paper or kapok paper to avoid direct contact between carbon steel and stainless steel. Packaging materials with good air permeability should be selected for the instruments to be wrapped, such as sterilization bags, sterilized paper, gauze fabric and other containers. The sterilization should be placed with the opening downward or side.

Preparation before sterilization of rubber tube:

Please clean the rubber tube with warm water first, and then place it on the sterilizer storage board. At the same time, make sure that the pipe is a hollow pipe with open ends, and there is no sharp turning, twisting or kinking.

Preparation before sterilization of dressing package:

Put the dressing bag vertically on the tray, and pay attention to avoid contact with the inner wall of the sterilizer, and there should be enough space between the packages.

4) Preparation before liquid sterilization

It is only limited to heat-resistant glass bottles. When loading, it should not exceed 1 / 2 of the volume, so as to avoid liquid overflow. Vent plug should be used to seal the container, and the cover should be fully relaxed. As the sterilizer is stainless steel, it is sensitive to chloride ionic liquids and easy to be corroded. Therefore, sterilize such liquids (such as normal saline, etc.) carefully

4. Close the chamber door of the sterilizer

After closing the cavity door, the following prompt interface appears, and operate according to the prompt:



After the prompt "please press the chamber door" disappears, the hand can leave the cavity door, and the prompt interface is as follows:

Inner Temp. 3	Door Status	re 99.3 kPa Negative Pressure Release	Object Temp. 30, 5	t
Automatic closing	, please wait			
Door Operation	Program	Manageme		Exit



 ${f C}$ Inner Pressure 99.3 kPa Object Temp. 31.9 30.5 C Inner Temp. Door Status Door Status Negative Pressure Release Locked Close 1 Door Management Operation

After the prompt disappears, the following screen appears:

When the door lock is locked, the prompt will disappear after a few seconds. (Note: the dooroperation cannot be performed when the sterilization program is started)

5.Select program



Click "program" on the main interface to enter the following figure:

There are three sterilization procedures, two test procedures, three self-defined programs and one new program. Sterilization procedures are: fabric program, device program, liquid program; test procedures are: BD test, vacuum test ; self defined programs are according to user's own needs.

Note: the operator account has no authority to set parameters, and the administrator account can set parameters.

5.1 Fabric program:

Sterilization process: standby \rightarrow preheating \rightarrow pulsation \rightarrow heating \rightarrow sterilization \rightarrow exhaust \rightarrow drying \rightarrow balance pressure \rightarrow end



Sterilization purpose: fabric sterilization, drying after sterilization Click the fabric program to pop up the following screen:

	Inner Temp. 31. 9	℃ Inner Pres	ssure	99.8 kPa Object	Temp. 30. 5	τ
Fabrics	Sterilizing Temp. Setting	134. 0	°C	Negative Pressure Pulse	3 times	
- U - U - U - U - U - U - U - U - U - U	Sterilizing Time Setting	4	min	Cross Pressure Pulse	1 times	
Preneating 	Sterilizing Pressure Setting	304. 0	kPa	Positive Pressure Pulse	3 times	
U Heating Un	Chamber Door Open Temp.	99. 0	°C	Drying Time Setting	10 min	
Sterilizing	Exhaust Speed Setting	5		Drying Press. Setting	330. 0 kPa	1
				Set Printin	g Enable Printing	
A Balance Press.						
	Start		Dig	sl Query	Jure	
Administrator	Door Operation	Program		Management	Flow Char	

Name	Default parameters	Program parameter range	
Sterilizing Temp.	134°C	115°C ~ 135°C	
Sterilizing Time	4min	1min ~ 6000min	
Sterilization pressure	It changes automatically with the	No modification is	
setting	change of sterilization	recommended	
	temperature		
Chamber cover opening	99°C	40°C ~ 99°C	
temperature			
Exhaust speed setting	Level 5	Level 0 ~ 5	
Negative pressure	3times	2 ~ 99times	
pulsation times			
Times of cross pressure	1time	1 ~ 99times	
pulsation			
Positive pressure	3time	1 ~ 99times	
pulsation times			
Drying time setting	10min	0~300min(0 means to cancel	
		the drying function)	
Drying pressure setting	330kPa	200-330kPa	

Note: 1) exhaust speed setting Description: stage 0 means exhaust closing, stage 1, 2, 3 and 4 exhaust time is lengthened in turn, and stage 5 is exhaust normally open.

2) When the administrator sets "print setting" to be controlled by "operator", the operator can start and stop printing

5.2 Device program:

Sterilization process: standby \rightarrow preheating \rightarrow pulsation \rightarrow heating \rightarrow sterilization \rightarrow exhaust \rightarrow drying \rightarrow balance pressure \rightarrow end

Sterilization purpose: instrument sterilization and drying after sterilization Click the device program to pop up the following screen:



Name	Default parameters	Program parameter range	
Sterilizing Temp.	126°C	115°C ~ 135°C	
Sterilizing Time	15min	1min ~ 6000min	
Sterilization pressure setting	It changes automatically with the change of sterilization temperature	No modification is recommended	
Chamber cover opening temperature	99°C	40°C ~ 99°C	
Exhaust speed setting	Level 5	Level 0 ~ 5	
Negative pressure pulsation times	3times	2 ~ 99times	
Times of cross pressure pulsation	1time	1 ~ 99times	
Positive pressure pulsation times	3time	1 ~ 99times	
Drying time setting	10min	0~300min(0 means to cancel the drying function)	
Drying pressure setting	330kPa	200-330kPa	

Note: 1) exhaust speed setting Description: stage 0 means exhaust closing, stage 1, 2, 3 and 4 exhaust time is lengthened in turn, and stage 5 is exhaust normally open.

3) When the administrator sets "print setting" to be controlled by "operator", the operator can start and stop printing

5.3 Liquid procedure:

Sterilization process: standby \rightarrow preheating \rightarrow cooling \rightarrow heating \rightarrow sterilization \rightarrow steam exhaust \rightarrow balance pressure \rightarrow end

Purpose: liquid sterilization

Click the liquid program to pop up the following screen:

	Inner Temp. 31. 9	℃ Inner Press	ure 99.3 kPa O	bject Temp. 30, 5 °C
Liquid	Sterilizing Temp. Setting	121. 0 °C	Exhaust Speed Setting	3 0
Standby ⊕	Sterilizing Time Setting	20 mir	Cold Air Purging Time Settin	g 10 min
Preheating &	Sterilizing Pressure Setting	205. 0 kl	Pa	
Air Purging	Chamber Door Open Temp.	80.0 [.] C		
Heating Up				
Sterilizing				
- Exhaust			Set P	rinting Enable Printing
⊕ Balance Press.				
End	Start		Dig tal Query	() yes
Administrator	Door Operation	Program	Management	Flow Chart

Name	Default parameters	Program parameter range
Sterilizing Temp.	121°C	115°C ~ 135°C
Sterilizing Time	20min	1min ~ 6000min
Sterilization pressure setting	It changes automatically with the change of sterilization temperature	No modification is recommended
Chamber cover opening temperature	80°C	40°C ~ 99°C
Exhaust speed setting	Level 0	Level 0 ~ 5
Cold Air Purging Time	10min	10min~20min

Note: 1) exhaust speed setting Description: stage 0 means exhaust closing, stage 1, 2, 3 and 4 exhaust time is lengthened in turn, and stage 5 is exhaust normally open.

2) When the administrator sets "print setting" to be controlled by "operator", the operator can start and stop printing

5.4 BD Program

It is suggested that the verification should be carried out once a week according to the provisions of 6.8.1.1 in gb8599-2008 standard.

BD	Inner Temp. 32. 0	℃ InnerI	Pressure	99.6 kPa Ob	ject Temp.	30. 6	r
BD Test	Sterilizing Temp. Setting	134. 0	°C	Negative Pressure Pulse	3	times	
Standby ⊕	Sterilizing Time Setting	4	min	Cross Pressure Pulse	1	times	
Preheating	Sterilizing Pressure Setting	304.0	k Pa.	Positive Pressure Pulse	3	times	
Pulse	Chamber Door Open Temp.	99. 0	℃				
Heating Up	Exhaust Speed Setting	5					
Sterilizing							
Exhaust				Set Pri	nting Enable	Printing	
Balance Press.							
End	Start		Dig	lQuery	Q	arve	
I	Door Operation	Program		Management		Flow Chart	

Name	Default parameters	Program parameter range
Sterilizing Temp.	134°C	115°C ~ 135°C
Sterilizing Time	4min	1min ~ 6000min
Sterilization pressure	It changes automatically with the	No modification is
setting	change of sterilization	recommended
	temperature	
Chamber cover opening	99°C	40°C ~ 99°C
temperature		
Exhaust speed setting	Level 5	Level 0 ~ 5
Negative pressure	3times	2 ~ 99times
pulsation times		
Times of cross pressure	1time	1 ~ 99times
pulsation		
Positive pressure	3time	1 ~ 99times
pulsation times		



Note: 1) description of exhaust speed setting: level 0 means no exhaust, stage 1, 2, 3 and 4 exhaust speed increases in turn, and level 5 means exhaust normally open. 2) When the administrator sets "print setting" to be controlled by "operator", the operator can start and stop printing

5.5 vacuum test program

3) Vacuum test process: standby \rightarrow vacuum pumping \rightarrow leak detection \rightarrow pressure balance \rightarrow end Purpose: sealing test

	Inner Temp. 32. 0	℃ Inner Pre	ssure	99.8 kPa Object I	emp. 30.6	r
Vacuum Test	Vacuum Time	0	min			
Standby	Test Time Setting	20	min	Test Countdown	20 min	
Vacuum						
Leak Detection				Set Printing	Enable Printing	
Balance Press.						-
End						
			Ste	art		
Administrator	Door Operation	Program		Management	Flow Chart	

Parameters do not need to be set.

Note: 1) when the administrator sets "print setting" to be controlled by "operator", the operator can start and stop printing

5.6 Self-defined program

Click the user-defined program (user defined 1, user defined 2, user defined 3) to enter the program. There are four modes to select, which are as follows:

Sterilization process:

5.6.1, Vacuum sterilization procedure: standby \rightarrow preheating \rightarrow pulsation \rightarrow heating \rightarrow sterilization \rightarrow steam exhaust \rightarrow drying \rightarrow pressure balance \rightarrow end

(application: solid sterilization and drying, also applicable without drying)

(Inner Temp. 32. 1	°C Inner	Pressure	99. 3 kPa Object	Temp. 30, 6	°C
Self-defined 1	Ste	erilizing Mode	Vacuum Ste	nlizing 💌		
3 Brokenting	Sterilizing Temp. Setting	134. 0	'C	Negative Pressure Pulse	3 times	
1 retirating	Sterilizing Time Setting	4	min	Cross Pressure Pulse	1 times	
Heating IIn	Sterilizing Pressure Setting	304.0	kPa.	Positive Pressure Pulse	3 times	
	Chamber Door Open Temp.	99. 0	ю	Drying Time Setting	10 min	
	Exhaust Speed Setting	5		Drying press. Setting	330.0 kPa	
Drying						
Balance Press.				Set Printing	Enable Printing	
8 End	Start		Digot	Query	 Image: Second sec	
Administrator	Door Operation	Pro	gram	Management	Flow Chart	

Parameter setting range:

Name	Default parameters	Program parameter range
Sterilizing Temp.	134°C	115°C ~ 135°C
Sterilizing Time	4min	1min ~ 6000min
Sterilization pressure setting	It changes automatically with the change of sterilization temperature	No modification is recommended
Chamber cover opening temperature	99°C	40°C ~ 99°C
Exhaust speed setting	Level 5	Level 0 ~ 5
Negative pressure pulsation times	3times	2 ~ 99times
Times of cross pressure pulsation	1time	1 ~ 99times
Positive pressure pulsation times	3time	1 ~ 99times
Drying time setting	10min	0~300min(0 means to cancel the drying function)
Drying pressure setting	330kPa	200-330kPa

Note: exhaust speed setting Description: stage 0 means exhaust closing, stage 1, 2, 3 and 4 exhaust time is lengthened in turn, and stage 5 is exhaust normally open.

5.6.2.Liquid program: standby \rightarrow preheating \rightarrow cooling \rightarrow heating \rightarrow sterilization \rightarrow steam discharging \rightarrow pressure balancing \rightarrow end Purpose: liquid sterilization



Name	Default parameters	Program parameter range
Sterilizing Temp.	134°C	115°C ~ 135°C
Sterilizing Time	4min	1min ~ 6000min
Sterilization pressure setting	It changes automatically with the change of sterilization temperature	No modification is recommended
Chamber cover opening temperature	99°C	40°C ~ 99°C
Exhaust speed setting	Level 5	Level 0 ~ 5
Cold Air Purging Time	10min	10min~20min

Parameter setting range:

Note: exhaust speed setting description: stage 0 means exhaust closing, stage 1, 2, 3 and 4 exhaust time is lengthened in turn, and stage 5 is exhaust normally open.

5.6.3.Special lower arrangement procedure: standby \rightarrow preheating \rightarrow cooling \rightarrow heating \rightarrow sterilization \rightarrow steam discharging \rightarrow drying \rightarrow pressure balancing \rightarrow end Application: solid sterilization, drying after sterilization

(Inner Temp. 32. 1	℃ Inne	r Pressure	99.6 kPa Object Te	mp. 30.6
elf-defined 1	Ster	ilizing Mode	Special low dr	ainage 🔻	
J.	Sterilizing Temp. Setting	134. 0	°C		
Treneating T	Sterilizing Time Setting	4	min		
A Heating IIn	Sterilizing Pressure Setting	304.0	kPa		
	Chamber Door Open Temp.	99.0	℃	Drying Time Setting	10 min
J. Frhaust	Exhaust Speed Setting	5		Drying press. Setting	30.0 kPa
Drving	Cold Air Purging Time Setting	10	min		
Uance Press.				Set Printing	Enable Printing
0 End	Start		DigtalQ	iery	Surve
ň	Door Operation	Pn	ogram	Management	Flow Chart

Name	Default parameters	Program parameter range
Sterilizing Temp.	134°C	115°C ~ 135°C
Sterilizing Time	4min	1min ~ 6000min
Sterilization pressure	It changes automatically with the	No modification is
setting	change of sterilization	recommended
	temperature	
Chamber cover opening	99°C	40°C ~ 99°C
temperature		
Exhaust speed setting	Level 5	Level 0 ~ 5
Cold Air Purging Time	10min	10min~20min
Drying time setting	10min	0~300min(0 means to cancel
		the drying function)
Drying pressure setting	330kPa	200-330kPa

Note: description of exhaust speed setting: level 0 means no exhaust, stage 1, 2, 3 and 4 exhaust speed increases in turn, and level 5 means exhaust normally open.

Positive pressure procedure (applicable to items that cannot be pre vacuum and need to be dried, and need multiple positive pressure)

(Realization function: articles that cannot be used for pre vacuum sterilization and need to be dried. If the vacuum drying is not possible, this mode can also be used when encountering the liquid that needs positive pressure. Only when the drying time is set to 0, the program will skip the drying stage and realize atmospheric pressure drying; you can also set the drying time not to 0 and run the drying stage until the end of the program.)

5.6.4.Auxiliary drying: standby \rightarrow preheating \rightarrow drying \rightarrow ending Application: suitable for fabrics or instruments which is no need of sterilization



Note: In the program, if the item is 0, such as 0 times of positive pressure pulsation, the next function will be automatically entered in order not to perform this function. Drying is atmospheric pressure drying, that is, drying under the same internal chamber pressure and external pressure.

Parameter setting range:

Name	Default parameters	Program parameter range
Chamber cover opening	99°C	40°C ~ 99°C
temperature		
Exhaust speed setting	Level 5	Level 0 ~ 5
Drying time setting	10min	0~300min(0 means to cancel the drying function)
Drying pressure setting	330kPa	200-330kPa

6.Start sterilization

Click the "start program" button, the system will pop up with the following prompt:



6.1 According to the prompt, click "confirm" to exit the prompt and start the program directly.

6.2 Click "Cancel" to stop, the system will return to the parameter interface

Note: please confirm the current inner chamber pressure before starting the program. If the inner chamber pressure is greater than 110kpa, there will be a "cavity pressure is too high" prompt, and the program cannot be started directly.

6.3 Click program to return to the program interface

6.4 After the sterilization is started, the interface will jump to the digital query interface (after the sterilization is started, the default is the digital query screen)

Evaporator Press.	97.9 kPa Inner Temp.	32. 1 °	C Inner Pres	s. 99.3 kP	a Object Ten	np. 30, 6	r
Fabrics	Sterilizing F0 Value	0.0		Total F0 Value	0.0		
Standby	Negative Pressure Pulse Times	3	times	Surplus	3	times	
Pulse	Cross Pressure Pulse Times	1	times	Surplus	1	times	
J. Heating Up	Positive Pressure Pulse Times	3	times	Surplus	3	times	
Sterilizing	Sterilizing Temp. Setting	134.0	۳C				
Exhaust	Sterilizing Time Setting	4	min	Surplus	4	min	
Drying U Balance Press.	Drying Time Setting	10	min	Surplus	10	min	
0 End							
	Door Operation Pro	ogram	Manager	nent Flow	7 Chart	Return	

6.5 Click the curve to pop up the following interface:



6.6 click the working condition to open the following interface:



6.7 Click stop program to pop up the following interface:

Selection confirmation: stop sterilization, return to program interface (if the original fabric program is returned to the fabric program interface), select Cancel and maintain the original.

K	Inner Temp. 27.	.1 °C :	Inner Pre. 100. 4	kPa Object Temp.	C
Fabrics	5			e <u>3</u>	times
Proheating	S			e <u>1</u>	times
Pulse	Ste	To suspend	l sterilizing?	е З	times
theating Un	Char			ng 10	min
Sterilizing				ng <u>33</u> 0	.0 kPa
Exhaust			·		
		YES	NO	inting St	p Printing
Balance Press.					
U End	Stop		Digital Query		Curve
Technologist	Dootoperation	Fogan	Man	egement	Flow Chart

7. Program end and cover opening



7.1 After the end of the program, it will automatically jump to the "door operation" interface. See operation in "2". The door opening mode is "open sterilization room door".

8. Take out the sterilized articles

1) Heat insulating gloves should be worn when taking out sterilized articles from the sterilization chamber, and hands can be put into the sterilization chamber after the steam disperses.

2) When sterilizing liquid articles, due to the slow cooling speed of liquid, it is necessary to confirm that the temperature has dropped to low enough before taking out, otherwise scalding may occur.

3) When taking out the basket with the thermometer, it should be noted that the thermometer can only be taken out if it is not jammed.

9.Turn off the Power

If the end of the day's sterilization work or not used for a long time, the power supply and water source switch should be turned off.

Description of each function management interface:

1.Operator management interface:

In operator mode, click manage to open the operator management interface, as shown in the following figure:

Note: the operator has no authority to modify parameters or set new program functions, but allows to view program parameters.



1.1.Click password setting to open the following interface:

	Modify passwo	rd	
	Old Password	0	
	Please input new password	0	
Switch to Operato	Confirm to mothfy	Return	witch to Engineer
Door Operation	Program Managem	mt Flow Chart	Help

The operator modifies his / her password. After the modification, click confirm to modify and click back to return to the operator management interface

1.2 .Click the alarm information to pop up the following interface:

It can query the fault information that has occurred up to the current time



1.3.Click switch to technologist, administrator and engineer in the following line to switch to the corresponding function interface.

2. Technician management interface:

In the case of technologist, click manage to open the technician management interface, as shown in the following figure:



2.1. Click parameter setting to open the following interface:

Inner chamber quick re	lease open 🚺 🧎	. 0 S		Inner chamber quick	release interval	50. 0	S
Upper limit in positive press	ure drying 11	5.0 kF	Pa Lov	ver limit in positive ;	pressure drying	102. 0	k
Upper limit in negative pres	sure pulse 🛛 8(). () kF	P a _{Lo}	wer limit in negative	pressure pulse	10. 0	k]
Upper limit in cross pres	sure pulse 12	0. 0 kF	Pa	Lower limit in cross	pressure pulse	60. 0	k]
Upper limit in positive pres	sure pulse 🛛 🚺 🕹	0.0 kF	Pa L	ower limit in positiv	e pressure pulse	110.0	kl
Water tank cooling ter	mperature 48	5. 0 °C					

Set the parameters related to sterilization (it is recommended not to set this parameter). The setting range is as follows:

Name	Default parameters	New program parameter
		range
Inner chamber fast exhaust switch on	1 second	1 ~ 9seconds
Upper limit of positive pressure pulse	115kPa	100 ~ 145kPa
drying		
Upper limit value of negative pressure	80kPa	3 ~ 99kPa
fluctuation		
Upper limit value of pressure	120kPa	100 ~ 170kPa

fluctuation		
Upper limit value of positive pressure	140kPa	100 ~ 170kPa
fluctuation		
Water tank temperature on	45°C	45 ~ 65°C
Rapid ventricular septum	50 second	30 ~ 90 seconds
Lower limit of positive pressure pulse	102kPa	100 ~ 110kPa
drying		
Lower limit value of negative pressure	10kPa	5 ~ 90kPa
fluctuation		
Lower limit value of cross pressure	60kPa	5 ~ 100kPa
fluctuation		
Lower limit value of positive pressure	110kPa	100 ~ 170kPa
fluctuation		

2.2. When you click the password setting, it is the same as the password setting in the operator interface.

2.3. When clicking the alarm information, it is the same as the alarm information in the operator interface.

2.4. When you click time setting, the following interface will pop up:



Current Time (hh:mm:ss): 11:49:37 Current Date (yy/mm/dd): 20/08/10 + Current Day of Week: Mon.	
Current Date (yy/mm/dd): 20/08/10 + + Current Day of Week: Mon.	
Current Day of Week: Mon.	
Hour: 11	
Minute: 48	
Second: 53	
Year: 20	
Month: 8	
Day: 10 Set & Exit	
Day of Week: Mon.	
RTC Adjustment (sec./day): 0.00 Cancel	
Time Zone: UTC + 08:00	
Retu	n

The time can be changed to the current time

2.5. Click the following line to switch to technologist, administrator and engineer to switch to the corresponding function interface.

3. Administrator management interface:

In administrator, click manage to open the administrator management interface, as shown in the following figure:



3.1. When you click parameter setting, it is the same as the parameter setting in the technician interface.

3.2. When you click password management, the following interface will pop up:

Parameter		Manage Password		Time
Mannual	Modify password	Set password	Return	
Swith to	Operator	Switch to Technolo	gist	Switch to Engineer
Door Operation	Program	Management	Flow Chart	Help

3.2.1. When you click modify password, the following interface will pop up (only change the password of the administrator himself)

	Modify password
Parameter	Old Password 0 Time
	Please input new password
Mannual	Confirm to modify Return
Door Operation	Program Management Flow Chart Help

If the number of digits is wrong or the password is wrong, there will be a corresponding prompt;

3.2.2. Click back to return to the management interface;

3.2.3. Click set password to pop up the following interface (user number and password of technologist and operator can be modified);

Set passwo		
User Account	1	_
Please input new password	0	
Confirm to modify	Return	

3.3 when the alarm information is clicked, it is the same as the alarm information in the operator interface.

3.4. When you click time setting, it is the same as the time setting in the technician interface.

3.5. When you click system manual, the following interface will pop up:

	Input	Data	Mannu	al
	Water Tank Temp. <mark></mark> [®] C	Door in place OFF	F2 Inner chamber steam inlet valve close	F9 Steam inlet val close
Develop	Water Pressure 428. 3 kPa	Door in open place ON	F3 Return air valve close	F10 Steam outlet va close
running	Vacuum Status Normal	Door in close place OFF	F4 Vacuum water inlet valve close	F11 Cooling down valve close
Status		Evaporator High Water Level Not full	F5 Air pump out valve close	Inflater Pump clos
		Evaporator Low Water Level <mark>Unlack</mark>	F6 Cooling down valve close	Evaporator heating close
		Electromagnet interlocking Not in place	F7 Water outlet valve close	Vacuum pump OF
			F8 Quick release valve close	Buzzer OFF

The corresponding solenoid valve can be opened and closed manually **3.6. When you click system settings, the following interface will pop up:**

Language	中文 English		D	Authority Tum on 💌 Control
Screen Touch Sound	OFF ON		Printer	Operator Control 60.0 S
Load Thermometer Option	Optional		Pressure Unit	kPa 🔽
Water Supply Pressure Option	Not optional		Pressure Unit	YYMM/DD 🔽
Printing Temperaure Type	Chamber Temp.		Sofety Volue Test	Test Cycle 353 Day
Printing Serial No.	YES		Salely Valve Test	Last Test Time <mark>2019 – 8 –22</mark>
Drying after liquid sterilizing	No drying		PLC Software Version	P1. 2. 0. 3. 3
Local Alltitude Setting	100	m	HMI Software Version	H1. 2. 0. 3. 3
				Return

It can modify the language, touch prompt tone, article temperature selection, water source pressure selection, print temperature type, print batch number, printer function, pressure unit, date format and altitude

3.6.1 when setting the printer function, the authority control includes closing, opening and operator control; the printing cycle is added and subtracted once for 30 seconds, and the minimum is 30 seconds

3.7. Click switch to technologist, administrator and engineer in the following line to switch to the corresponding function interface.

3.8. Help interface

Click the "help" button in the system login interface to enter the help interface:

Contact u	E	
	Zealway Instrument Inc.	
	E-mail: sales@zealway.us	
	Veb:www.zealway.us	
		Return

This interface has the company's after-sales contact information, you can seek help from the company through the above methods

Welcome to scan the two-dimensional code, pay attention to the official account of the micrometer.

3.9 When you click the deviation compensation setting, the following interface will pop up:

Item	Test Va	alue	Correction	Value	0	Calibration V	alue
Evaporator Press.	98. 9	kPa	-1.0	kPa		97. 9	kPa
Inner Pressure	99. 8	kPa	-0. 5	kPa		99. 3	kPa
Inner Temperature	33. 0	°C	-0. 7	°C		32. 3	°C
Object Temperature	32.1	۳C	-1.3	•C		30.8	°C
Water Tank Temp.	33. 3	۹C	0. 0	۰C		33. 3	۳C
Factory Reset							
a Return							

3.9.1 Click the correction value in the middle of the above figure to modify, with the correction range of - $5 \sim + 5$. For example, if the measured value is 320.7kpa and the correction value is -0.7kpa, the corrected value will display 320kpa and take effect. Press to return to exit the interface.

Options

Note: only administrators have permission to set this item.

1) Selection of water source pressure sensor or pressure gauge:

•If the machine needs to monitor the water source pressure in real time, the water source pressure sensor can be selected and the "management" in the "Administrator" interface can be selected

In the parameter setting of the interface, the water source pressure sensor is enabled and the corresponding pressure unit (kPa, bar or PSI) is set.

•If it is necessary to observe the water source pressure, it is recommended to install a water source pressure gauge for the tap water source, so as to observe the water source pressure and ensure the normal use of the instrument.

2) Load thermometer (object temperature sensor)

•When the cavity temperature reaches the preset sterilization temperature, the temperature of the articles has not reached the preset sterilization temperature. In order to obtain the ideal sterilization effect, it is suggested to use the article thermometer to monitor the actual temperature of the articles in real time. In this case, the sterilization timer starts only when the actual temperature of the article reaches the preset sterilization temperature.

 \star Note: 1. After the item thermometer is installed, enter the administrator settings menu and

set the item thermometer to enable. After use, the article thermometer should be placed on the storage board and properly kept. If it is not used for a long time, the article thermometer can be set to not be enabled in the user administrator settings menu.

•It is recommended to test once a month. The test method is as follows: fix the stainless steel tube of the article thermometer in the grid close to the temperature sensor (filter screen position), set the sterilization temperature at 121 °C and the sterilization time for 20 minutes. When the sterilization reaches 5 minutes, confirm that the error between the display temperature of the article thermometer and the temperature sensor is 0 $^{\circ}C \sim + 0.2 ^{\circ}C$. If it is no longer in this range, modify the thermometer compensation value of the article (enter the "deviation compensation" of "Administrator" to modify).

Note: after installing the item thermometer, you should enter the "system settings" interface of "Administrator" to set the item thermometer to

Enable. After use, the thermometer should be placed on the water level board and kept properly. If it is not used for a long time, it should not be turned on Use.

3) Printer

•If you need to record the sterilization process parameters, you can choose a printer. To enable, you need to enter the "system settings" interface of "Administrator", set the printer to enable, and you can select the printing language (English or Chinese) and the display mode of date (day / month / year or year / month / day) according to the needs.

Install the roll

OPut the paper roll into the paper inlet, press the paper in and out button to confirm that the paper roll comes out of the paper outlet, and then press the paper in and out button, and the paper roll will stop automatically.



Paper outlet

Paper inlet

Paper roll

Paper in and out button

•The printer prints the F0 value

A total of 2 F0 values are printed, one of which is in sterilization state, and the other is F0 value of heating process, sterilization process and cooling process in the program process. The contents of the printer are as follows:

F0 value in sterilization stage: XX

Total F0 value: XX (excluding F0 value of steam stage after sterilization)

The F0 value is displayed with 1 decimal point. If an article thermometer is used, the calculation shall be based on the value recorded by the article thermometer, otherwise, it shall be calculated

According to the cavity temperature value.

4) Selection of back pressure system:

Ensure that there is no negative pressure in the chamber during the rapid cooling process after sterilization in liquid mode. Even if the liquid loading is more than 80% of the container, the sterilization liquid will not overflow during the sterilization process

•Back pressure system: pipeline diagram of air compressor

04 Maintain&Management

1. The correct use and routine maintenance of sterilizer equipment are necessary conditions for prolonging the service life of the equipment and reducing failures. This chapter will briefly introduce the precautions and maintenance work in use.

2. Before maintenance, please cut off the power supply and confirm that the sterilization chamber is cooled before maintenance or maintenance.

Instrument maintenance

1. Daily maintenance of main parts and components

•Sterilization room and storage board: when the temperature of sterilization room and storage board drops to close to room temperature, scrub with neutral detergent, then rinse with tap water, and finally dry it with soft cloth.

When cleaning the inner room, effective measures should be taken to prevent the door from closing automatically, so as to avoid danger from harming people's life.

•Filter screen: clean the filter attached to the filter screen at the front of the sterilization room to ensure the vacuum rate and condensate flow, and the temperature indication is consistent with the pressure.



 Vacuum pump: If the working fluid used is harmful to the human body or the environment, for the sake of safety, clean the pump with clean water before emptying.

When emptying, open the screw plug (068a) on the pump cover (061a) to let the liquid flow out. Rotate the pump by hand until no liquid flows out.

By tilting the pump cover 45 degrees, the pump can be basically emptied, so that even if the pump is stopped for a long time or in the mist, the pump will not be damaged.

1) Treatment method for long-term pump shutdown:

2) If the pump stops conveying for about 4 weeks, the liquid in the pump should be drained first. For cast iron pump, the pump should be preserved with 1 / 2 I anticorrosive oil, and the oil should be poured into the suction port and exhaust port. If the impeller is stuck after stopping the pump for a long time due to the use of hard water, the pump cavity should be filled with 10% oxalic acid for about 30 minutes.

3) If the pump is stopped for about 4 weeks, remove the inlet and outlet of the pump and blow the water inside with an air gun. If the water in the inlet and outlet pipes is not drained, it is easy to stink or scale, causing impurities to affect the normal use of the pump chamber.

• Drain valve (water vapor separator): the normal operation of a drain valve in the drain pipe directly affects the sterilization effect. If the accumulated water can not be discharged normally, the drain valve should be opened for cleaning. During operation, it is normal to have a small amount of steam discharged, but the exhaust volume should not be too much. The drain valve should be cleaned once every three months.

•Safety valve: it is installed on the sterilizer to protect equipment and personal safety. It is installed on the top of the container. Its opening pressure has been adjusted before delivery. It is forbidden to adjust it at will to avoid accidents due to improper adjustment. However, it is necessary to lift the handle several times every half a year and wash it with steam to avoid its action failure.



Please do not touch the handle of safety valve directly by hand to avoid steam scald. Screw driver or other strip-shaped objects can be inserted into the pull ring to lift the handle.

•Sterilization filter: the filter works when the air enters into the room. After the gas is filtered, it can effectively prevent the air entering the inner room from polluting the sterilized articles. The filter is required to be replaced regularly. For details, see the list of consumables in Chapter 4, Part 5.



•**Filter**: a filter is installed on the water inlet and exhaust pipe, which is used to filter out the impurities in the exhaust pipe and water inlet pipe. Ensure the smooth pipeline and reliable action of each valve. Therefore, regular cleaning is required to prevent blockage. When cleaning, screw out the screw below, clean the filter screen, and then tighten the screw plug.





•Check valve: the vacuum pipeline and water inlet pipeline are equipped with check valve, which should be checked regularly to avoid foreign matters affecting its one-way sealing.

•Electric heating tube: the electric heating tube is located at the back end of the steam generator. Once it is damaged, it will affect the normal use of the whole equipment. Therefore, the electric heating pipe should be inspected regularly. If it is damaged, it should be replaced in time. After replacement, it should be ensured that there is no steam leakage and water leakage under the working pressure, and the grounding should be firm and reliable. b: The parts of pressure gauge shall be assembled firmly without looseness;

c: The pressure gauge indexing plate shall have the following marks: manufacturer, product

name, measuring unit and number of digital measuring instruments; vacuum gauge shall have "negative" sign, accuracy grade and factory number.

d: The pressure gauge shall be equipped with safety hole, and dust-proof device shall be installed on the safety hole.

e: The pointer indicating end should cover 1 / 3 of the shortest dividing line length, and the pointer width should not be greater than the dividing line width.

f: In the reading part, the pressure gauge glass shall be colorless and transparent, and there shall be no defect hindering the reading.

In case of any discrepancy, please stop sterilization immediately and contact the distributor.

2. Maintenance of electrical components:

•Electrical components and connecting wires are strictly prohibited from contacting with water. If water is accidentally stained, the power supply can be connected only after treatment.

•Electrical components should be dustproof. It is recommended to remove dust once a season

•All connecting wires, plugs and sockets, and components should be checked for looseness frequently. If they are loose, they should be plugged and locked immediately.

Instrument maintenance

1. Inspection of leakage protection circuit breaker (main power switch)

•Press the test button of the leakage protection circuit breaker. If the leakage protection circuit breaker trips, it indicates that it is normal. Otherwise, please close the leakage protection circuit breaker and contact the dealer.

2.Safety valve working condition test

•Enter the technologist interface, enter the "process parameter setting", click the safety valve test, press "start safety valve", the instrument will automatically jump to the safety valve interface and start the safety valve test, and press "stop" to terminate the safety valve test. •When the temperature rises to 142 ~ 145 °C, the pressure rises to 0.275 MPa ~ 0.30 MPa. The safety valve will trip to exhaust steam. At this time, the pressure will no longer rise or slowly fall back, indicating that the safety valve is normal.•If the temperature exceeds 145 °C, the safety valve does not trip, then the safety valve is abnormal. Must stop immediately and contact the dealer.

3.Replace the sealing ring

•Open the chamber door \rightarrow Click to open the door. The vacuum pump starts to extract steam. When the air pressure is less than 100KPA, the valve will open automatically.

OPull the corner of the seal ring to take it off directly

OClean the dirt on the sealing ring and the contact part between the cavity and the sealing ring with a cloth.

OWet the seal ring with water, the thick side of the seal ring faces inward, and then press the seal ring into the groove in an average section

OOperate the instrument according to the conventional method and observe the air tightness of the sealing ring.

4.Parts list that must be checked or provided by the manufacturer or its agent

Number	Part name	Number	Part name
1	Temperature switch (for heating	14	Air filter
	pipe)		
2	Cavity temperature sensor	15	Water level sensor
3	Article temperature sensor	16	Interlocking microswitch
4	Pressure sensor	17	Power switch
5	Electric heating tube	18	Vacuum pump
6	Seal ring	19	Booster pump
7	Pressure gauge	20	Safety valve
8	Stainless steel basket	21	Leakage protector
9	AC contactor (for heating	22	Touch screen
10	Solid state relay	23	Key switch
11	PLC	24	Emergency stop switch
12	Solenoid valve	25	Small relay of distribution box
13	Printer		

5.List of consumables:

Number	Part name			
1	Electric heating tube			
2	Seal ring			
3	Sterilizing filter			
4	Article temperature sensor (optional)			

05 Troubleshooting

Equipment alarm information table:

Number	Alarm information	Possible causes	Handling suggestions
1	E01	Sterilization chamber door not locked or chamber door interlock failure	Check whether the switch handle is in place; check whether the door is closed in place; check whether the electromagnet mandrel top is in place; check the 24 V and 15 V power supply
2	E02	Heating pipe damaged	Check the working state of heating controller; Check the working indicator of solid state relay
3	E03	Temperature switch or pressure switch failure	Turn off the temperature control switch and the pressure switch. Check whether the tap water is sufficient
4	E04	Abnormal fluctuation of	The chamber temperature is

		power supply	higher than the boiling point temperature during start-up operation; The temperature of sterilization is higher than that of sterilization
5	E06	Inner chamber temperature control probe failure	Plug in the temperature sensor again
6	E10	Sterilization is not normally interrupted	Re sterilization
7	E11	Abnormal water inflow	Confirm whether the stop valve of demineralised water of drain valve seat is open; Whether the booster pump is overloaded; Check whether the tap water is sufficient.
8	E12	Micro pressure switch failure	Contact local dealer
9	E16	Temperature control probe failure	Plug in the temperature sensor again
10	E19	Dirty level sensor	Check and clean the liquid level sensor in the cavity for dirt
11	E21	The leakage of the instrument exceeds the standard	Check whether there are foreign matters in the vacuum pump; Replacement of sealing ring
12	E22	The leakage of the instrument exceeds the standard or there is water	Remove the water in the inner chamber; Replace the door sealing ring
13	E23	Temperature and pressure do not reach the set value	Contact local dealer
14	E24	The safety valve is unqualified	Replace the safety valve
15	E26	Inner chamber pressure sensor failure	Replace the pressure sensor
16	E36	Water source pressure sensor failure	When selecting the water source pressure sensor, replace the water source
17	E46	Interlayer pressure sensor failure Evaporator pressure sensor failure	Replace the sandwich pressure sensor Replace the evaporator pressure sensor

Instrument **Possible causes** Handling suggestions failure 1. After the power 1. Check whether 380V power 1. Check touch screen power is turned on, the supply is normal supply 2. Check the 24 V power supply touch screen does 2. Touch screen power not not light up connected 3. No 24 V power supply 1. The door is not closed 2. The program 1. Please close the cavity door does not start 2. In detection state 2. Please switch to the running state, which can be switched in the detection interface 3. Communication 1. The interface in communication 1. Replace the communication line and check whether the interface is interruption is burnt out by hot plug 2. Poor contact burnt out 4. Touch screen 1. Enter screen saver 2. Reconnection after power failure black screen 2. Touch screen failure 5.Pumping 1. Whether the water source 1. Check water source pressure vacuum is too pressure is satisfied 2. Check the connecting parts of slow, negative 2. There is a leak in the line each pipeline and conduct 3. Is there dirt on the surface of pressure maintaining test pressure can not reach the standard 3. Check for water supply the sealing ring 4. Replace the pressure feeder 4. No water source 5. Check whether the vacuum 5. Pressure transmitter failure valve is powered on 6. Vacuum valve not open 6. Check the sealing ring 6. The vacuum 1. Check the water source 1. No water supply pump is noisy 2. Reversal of vacuum pump 2. Adjust any two power lines 3. Serious scaling of vacuum 3. Descaling for vacuum pump and vacuum pipeline pump 1. Solenoid valve coil damaged 7. Solenoid valve 1. Replace the solenoid valve 2. The solenoid valve is not 2. Check the power supply of does not operate powered solenoid valve 1. Cleaning the surface dirt of float 1. Poor insulation or short circuit 8. There is no water in the steam of floating ball level gauge ball level gauge 2. The water inlet solenoid valve 2. Repair or replace the water inlet generator

2. Analysis and elimination of common faults

	does not work when it is powered	solenoid valve	
	on	3. Repair or replace the water inlet	
	3. Water pump damaged	pump	
9. The water in the	1. Floating ball level gauge failure	1. Replace or repair the float level	
steam generator is	2. Water inlet solenoid valve not	gauge	
not stopped	closed tightly	2. Repair or replace the water inlet	
		solenoid valve	
10. Heating tube	1. Heating pipe damaged	1. Replace the heating pipe	
not heated	2. Solid state relay failure	2. Repair or replace the solid state	
	3. Protect AC contactor from	relay	
	power failure	3. Check the power supply of AC	
		contactor	



06 Electrical Wiring Diagram and Component Distribution Diagram

1. Electrical wiring diagram



Symbol description						
Number	Symbol	Name	Number	Symbol	Name	
1	QF0	Total power supply	15	KA4	Electromagnet pull in relay	
2	QF1	Motor and control power supply	16	KA5	Electromagnet protection relay	
3	QF2	Control power supply	17	HMI	Touch screen	
4	SPD1	Surge protective device	18	PLC1-1	PLC IVC1-1614MAR6	
5	PBU1	Buzzer	19	PLC1-2	PLC IVC1-0016ENR	
6	KM1	Booster pump	20	PLC1-3	PLC IVC1-4PT	
7	KM2	Vacuum pump	21	PLC1-4	PLC IVC1-4AD	
8	KM3	UC heating contactor	22	U1	DC24 V power supply	
9	FR1	Booster pump thermal relay	23	U2	DC15V power supply (printer power supply)	
10	FR2	Vacuum pump thermal contactor	24	U3	DC5V power supply (electromagnet holding power supply)	
11	SSR1	UC heating solid state relay	25	FU1	24 V power supply protection fuse	
12	KA1	Automatic door closing relay	26	FU2	5V power supply protection fuse	
13	KA2	Automatic door opening relay	27	FU3	15V power supply protection melting	
14	КАЗ	Safety valve switch relay	28	FU4	AC220 protective fuse	

2. Component distribution diagram



www.biolabscientific.com 51







BIOLAB SCIENTIFIC LTD.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8 Canada Email: contact@biolabscientific.com Tel: +1 707 533 1445 Website: www.biolabscientific.com