



## SAMPLE CONCENTRATOR

BFA1CE1

# INDEX

1. Introduction	2
2. Specifications	3
2.1 The normal working condition.	3
2.2 Basic Parameters and Performance.	4
3. Preparations	5
3.1 Structure Description	5
3.2 Instrument Installation	6
3.3 Operation panel	8
4. Operation Guide	9
4.1 Temperature and time set	9
4.2 Operation and stop:	11
4.3 Temperature error adjustment	12
4.4 The exchange of the block	14
5. Error Analysis and Recovery Processing	16
6. Wiring Diagram	17

# 1. Introduction

Sample Concentrator (visible) is mainly used for concentrating or preparing sample in batch. Such as drug screening, hormone analysis, liquid phase and mass spectrometry in the analysis of sample preparation. It works by blowing nitrogen, in the surface of sample which is being heated to accelerate evaporating and separating the solvent in the samples without oxygen.

Features:

1. The heater causes the sample to be rapidly heated to the evaporation temperature, and at the same time, the gas is blown to the surface of the solution through the gas needle, which promotes the rapid evaporation of the solution and the concentration of the sample.
2. The height of the air chamber assembly is adjustable, and the standard air needle length is 150mm.
3. When the sample is concentrated, the liquid level position of the concentrated sample can be observed at any time.
4. Gas needle is controlled independently. Separately blow of each needle and flow regulating of each needle are available to avoid gas waste.
5. The entire equipment can be put into ventilation cabinet when the concentration sample in toxic solvents.
6. Built in overheat protection, automatic fault detection and fault beep alarm devices.
7. LED display immediate temperature and diminishing time. Operation is simple and convenient.

## 2. Specifications

### 2.1 The normal working condition.

The room temperature: 4C 45C

The relative humidity: ≤70%

The using power: AC220V / AC110V ~ 50-60Hz

### 2.2 Basic Parameters and Performance.

Model Parameters	BFA1CE1
Temperature Range	R.T.+5°C ~ 150°C
Temp. Setting Range	5°C ~ 150°C
Temperature Accuracy (40~100°C)	±0.5°C
Temperature Accuracy (100~150°C)	±1°C
Temperature Uniformity (100°C)	±0.5°C

Blocks Quantity	1 Standard Visible block
Time Range	1min~99h59min
Needle Plate Max. Lift Stroke	285mm
Gas-in Joint Outer Diameter	Φ7mm
Nitrogen Pressure	≤0.1MPa
Nitrogen Flow Rate	0~10L/min
Needle Length	150mm
Power	400W
Dimensions (mm)	260 x 220 x 525
Weight(kg)	5.8kgs

Table 1

## 3. Preparations

This chapter introduces the instrument's mechanical structure, keyboard and each key's function and some preparations before power on. You should be familiar with this chapter before the instrument is first operated.

### 3.1 Structure Description



Figure 1

## 3.2 Instrument Installation

### 3.2.1 Assembly




Figure 2

Installation steps:

- 1) Put the dry bath on the bench smoothly and solidly
- 2) Fully insert the column into a square hole of the column holder behind the host, then screw knob into the column mounts, and finally tighten the knob to secure the uprights
- 3) Screw fixed knob of the mobile seat into the threaded hole on the left side of the mobile seat, then hold knob of mobile seat and counterclockwise rotate the hand wheel adjustment. When reach the desired position, we can tighten the knob to fix mobile seat. Clockwise rotating the hand wheel adjustment, the mobile seat will decline (as shown in the picture).

Note: When move the seat, two hands must be used in conjunction with operating mobile seat fixed knob and hand wheel adjustment

- 4) Screw fixed knob of the gas chamber into the threaded hole on the left side of the movable block, then insert column of gas chamber into the corresponding hole of the front movable block. Fully inserted column of gas chamber, then tighten fixed knob. (as shown in the picture)

3.2.2 Gas Control Valve (optional purchasing auxiliary) gas-out port

Structure




Figure 3

Pull out the gas pressure adjust knob, clockwise rotate the knob to increase pressure, anticlockwise rotate the knob to decrease pressure. Push back the knob to lock pressure.

Clockwise rotate the gas flow adjust knob to increase flow rate, anticlockwise rotate the knob to decrease flow rate till shut off it.

### 3.2.3 Installation of nitrogen flow rate adjusting valve

The composition of nitrogen flow rate adjusting valve.




Figure 4

As the above picture, connect the outlet of the nitrogen flow regulator valve with outlet port of the host with a short hose (about 1.5m). See page 8 the no. 6, Connect the inlet of nitrogen flow regulator valve with outlet port of the host with a long hose (about 3m)

Note: Slowly open the valve of the nitrogen bottle to control the outlet pressure between 0.1MPa to 0.2MPa, and then open the valve adjusting knob to keep pressure value at 0.02MPa. According to the number of holes, we can adjust the pressure appropriately (usually between 0.02MPa and 0.05MPa optional).

## 3.3 Operation panel




Figure 5

### 3.3.1 Button instruction

▼ Decreasing button: decrease the figure set;

▲ Increasing button: increase the figure set;

START/STOP : Run/stop button: Run when you press the button once after setting the temperature and the time. Stop when you press this button for 2s.

## 4. Operation Guide

### 4.1 Temperature and time set

a). The LED will display "8" as the chart when the Instrument powers on and the Instrument goes into the initial state with the sound of "du...".




b). About 2s later, the figure 28.5 is the block's current temperature; 00:35 in the time display is the last set time.






Press ▲ or ▼ of the "Temp" key and keep off at once. Now, the value in the temperature display is the former setting temperature. As shown in the left drawing, the last digital of the setting temperature is flickering. If you want to set the temperature to 55.5°C, do as follows:






Keep pressing  , when it reach 55.5°C, stop Pressing, it will auto save as 55.5°C  
Keep pressing the above keys for 2 seconds. We can modify the value fleetly. It is very convenient.



c). Press  or  of the "Time" key and keep off at once, the value in the time display is the former setting time. Shown in the left drawing is 00:35 (35minutes). At same time, the last digital is flickering.  
If you want to change the time to 01:20, do as follows:



Keep pressing  , when it reach 01:20, stop Pressing, it will auto save as 01:20.

Notes: If the time setting is 00:00, meaning the time of operation is forever, the Instrument run constantly in setting temperature.

## 4.2 Operation and stop:

a) After accomplishment of the temperature and time setting, press the button of START/STOP and keep off at once, the Instrument begins operating. The temperature rises with the sound of the "du...".



Meanwhile, there is current temperature in the temperature display, and the radix point is flickering regularly during the rise process.

When the temperature stops to rise, the flickering also stops. " : " in the time display begins flickering, then the counting down starts.



When the timing is stopped, the operation is also stopped. The buzzer alarms. Current temperature in the temperature display and "over" in time display mean "over", the accomplishment of the operation.



- b). After accomplishment of the operation, the instrument is in standby mode waiting for new declaration. Then pressing the keys near "Temp" or "Time" can reset the temperature or the time. Pressing START/STOP, it will operate according to last setting.
- c). Press START/STOP for 2s during the operation, it will stop running. Press this button again, operation will continue.

## 4.3 Temperature error adjustment

The temperature of the Instrument has been adjusted before it is sold out. In case there is deviation between the actual temperature and the displayed temperature, you can do as follow for temperature calibration.

### NOTICE:

- 1) The instrument has 2 calibration temperature points to ensure the veracity. It is linearly adjusted on 40C and 100C. The temperature accuracy will be within  $\pm 0.5C$  after temperature calibration.
- 2) Both the environmental and the block temperature should be lower than 35C when calibration.

### Adjustment Methods:

- a) Power on the instrument, it enters into waiting interface. Make sure the current temperature in display is below 35C. If the temperature is higher than 35C, please wait until it down below 35C.
- b) Inject olefin oil into one of block well, and then put a thermometer into this well (the precision of the thermometer should be 0.1C and the temperature ball should be absolutely immersed into the olefin oil in the block well). Adiabatic material is needed on the block to separate it from the circumstance. (refer to below Figure 6)







Figure 6


NOTICE: To ensure the calibration precision, read the actual temperature value after the temperature reaches calibration points for at least 20 minutes.

c). Press  and  simultaneously to enter the temperature calibration interface as the right chart.

LED display "  " in the time display, the temperature display the current temperature and the program auto control the temperature to 40C.




When the temperature reaches 40C, the decimal digit begins to flicker. Wait for at least 20 minutes, after 20 minutes of constant temperature, read the measured temperature of the thermometer.



Notice! In order to ensure the temperature calibration accuracy, it is recommended that users read the temperature after 20 minutes of constant temperature  
Measured temperature!

If the thermometer reads 39.6 °C, the actual temperature of thermometer is 39.6C, press ▲ or ▼ of temperature to amend the temperature value to 39.6, Press "Start/Stop" to confirm.



d). Then the instrument automatically heats up to 100°C, and waits for the temperature calibration value to be input after a constant temperature of 100°C. It is also required to read the measured temperature of the thermometer after 20 minutes of constant temperature.

Notice! In order to ensure the temperature calibration accuracy, it is recommended that users read the temperature after 20 minutes of constant temperature  
Measured temperature!

If the actual temperature of thermometer is 101.5C, press ▲ or ▼ of temperature to amend the temperature value to 101.5, Press "Start/Stop" to confirm.



After temperature have been adjusted, press ▲ and ▼ simultaneously to exit the temperature calibration and return to waiting interface.




Notice:

After temperature calibration, the temperature display is the same as the practical temperature of block.


Press ▲ and ▼ simultaneously during the calibration to exit the adjustment procedure, the calibration is invalid.

## 4.4 The exchange of the block


a) Pull out the two screws which fix the block to the heating board with the screwdriver.





b) Fix the raiser in the center well of the block.



c) Pull out the raiser with the block.



d) Screw the raiser out from the block, then fix the raiser to another block needed. Put the block in the instrument and fix it with screws.



Notice: It is a little complicated to exchange the block, please be a bit patient.

## 5. Error Analysis and Recovery Processing

NO	Fault phenomenon	Cause analysis	Recovery processing
1	No signals on the display when the instruction is powered on.	No power	Check the connection of power
		Broken fuse	Exchange fuse(250V 4A Φ5x20)
		Broken switch	Exchange the switch
		Others	Contact to the seller
2	The actual and displayed temperatures are quite different.	Broken sensor or loose contact of the block	Contact to the seller
3	" ERR " in the display with the alarm of "du..."	Broken sensor or the environmental temperature is below zero.	Contact to the seller
4	no heating of the block	Broken sensor.	Contact to the seller
		Solid state relay damage	
		Broken heater	
5	Key doesn't work	Key broken	Contact to the seller

Table 2

## 6. Wiring Diagram

(Below diagram is just for reference. It is subject to change without prior notice.)




Figure 7



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
 Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)