



DUCTED FUME HOOD

BFSD-401 BFSD-402 BFSD-403 BFSD-404

INDEX

1. Unpacking, Installation and Debugging	2
1.1. Unpacking of Main Body	2
1.2. Unpacking of Base Cabinet	4
1.3. Accessories Checking	5
1.4. Installation Conditions and Operating Environment	6
1.5. Installation	6
1.6. Inspection after Installation	8
2. User Instructions	9
2.1 Functions	9
2.2. Product Structure	10
2.3. Control Panel	13
2.4. Regular Maintenance	16
2.5 Wiring Diagram	17
3. Trouble Shooting and Labels	19
3.1 Common Failures and Solutions	19
3.2 Label Description	21
4. Warranty	22

1. Unpacking, Installation and Debugging

Please firstly check whether the packing box is in good condition. If the packing box is damaged, please take photos and contact the freight carrier. Biolab and its dealers are not responsible for shipping damages.

1.1. Unpacking of Main Body

Choose a proper unpacking method according to the actual situation.

1.1.1. For wooden box

a) Method 1 Necessary tools for unpacking: Electric drill with hexagon dead M8



Figure 1

b) Method 2 Use M8 wrench to unpack



Figure 2

The following diagram demonstrates quick unpacking procedures (Figure 3). Remove the screws shown in the diagram below, then move the wooden pieces to right and left.




Figure 3

1.2. Unpacking of Base Cabinet

Choose a proper unpacking method according to the actual situation.

Please refer to the main body unpacking method (of the carton box) in the previous pages, use a scissor to cut the packing tape, take off the package cover and move up the box body.




Figure 4

1.3. Accessories Checking

Accessories are placed inside the working area and base cabinet. Please refer to 1.5.6. to take out all accessories and check the completeness referring to this packing list.

Packing list (BFSD 400 series).

Main body box:

No.	Items	Quantity	Remarks
1	Upper cabinet	1	
2	Lower cabinet	1	
3	Certificate	1	
4	Warranty card	1	
5	Power cord	1	
6	Testing report	1	
7	Exhaust pipe (with two clamps)	1	BFSD-404 has two exhaust pipes (4 with clamps)
8	PP centrifugal fan	1	BFSD-404 is two
9	PP small sink (with drain hose)	1	Optional
10	Faucet	1	Optional
11	Faucet control valve	1	Optional
12	Air faucet	1	Optional
13	Air faucet control valve	1	Optional

14	Countertop	1	
15	Stainless steel bolt	2	

Table 1

1.4. Installation Conditions and Operating Environment

1.4.1. Location requirements

The installation position of the fume hood should avoid facing the road and places where people pass frequently, avoid blocking windows and places with unfavorable ventilation and lighting, avoid blocking entrances and exits that affect the switch of the door, and avoid facing and placing them on corner walls.

1.4.2. Environment requirements

a) Only applicable to indoor operation

b) Ambient temperature: 15°C~35°C

c) Relative humidity: ≤75%

d) Atmospheric pressure range: 70 kPa~106 kPa

e) Power supply: the same as the rated voltage and rated frequency of the cabinet, see 2.1.4 technical parameters and performance indicators table;

f) The power supply needs to be reliably grounded; (judgment method: use a multimeter to test the live wire-to-ground voltage of the power supply, which should be the grid voltage, and the neutral wire to ground should be 0, otherwise the power supply grounding is not good).

1.5. Installation

1.5.1. Remove all the package materials

1.5.2. Check the surface of the base cabinet and the main body to make sure there is no scratch, deformation or foreign bodies

1.5.3. Carefully check the attachments and materials against the packing list in the instruction manual;

1.5.4. Move the entire equipment to a site that is as close as possible to the final location and is convenient for installation;

1.5.5. Check whether the ambient voltage frequency is consistent with the voltage frequency displayed on the label, and then power on.



Figure 5

1.5.6. Connect base cabinet with main body

Put the bottom cabinet in a suitable installation position, then put the fume hood cabinet on top of the bottom cabinet, adjust the upper cabinet to be flush with the back and both sides of the bottom cabinet, and fix it with bolts, and the fixing points are behind the two sides of the cabinet.




Figure 6

Put the pipe clip on the outer pipe first, then put the outer pipe on the air outlet of the cabinet, and fasten the pipe clip with a flat screwdriver; BFSD-404 has two air outlets, and the main pipe needs to reserve two air inlets for the installation of the exhaust pipe of the fume hood; BFSD-401, BFSD-402, and BFSD-403 all have one air outlet;

(4) Put on the countertop, install the sink on the countertop, and then connect the drain pipe of the sink in the bottom cabinet to the drain outside the cabinet;

1.6. Inspection after Installation

Refer to this table and follow the instruction in 2.3.2, check the following items after powering on the Fume Hood.

Checking Items	Normal working status
Power status	Equipment could be powered on/switched off when press the power button
Fan	Runs normally after pressing the Fan button; speed could be adjusted by pressing the adjusting button
Fluorescent lamp	Lamp lights up after pressing button
Socket	Use multimeter to test voltage output after pressing the socket button

Table 2



NOTE: Please contact Biolab technical department or agent for inspection or trouble shooting when problems could not be solved.

2. User Instructions

2.1 Functions

2.1.1. Product concept

In the chemical laboratory, various harmful and odorous gases, odors, moisture, and corrosive substances will be produced during the experimental operation. In order to protect the safety of users and prevent the pollutants in the experiment from spreading to the laboratory, a fume hood should be used near the pollution source.

2.1.2 Operating principle/air flow pattern




Figure 7

2.1.3. Protected object

Safety is the biggest mission pursued by fume hoods. The use of fume hoods in laboratories is to ensure the safety of users and prevent pollution to the surrounding environment.

2.1.4 Technical parameters

Model Parameter	BFSD-401	BFSD-402	BFSD-403	BFSD-404
Rated Voltage AC	220V±10% 110V±10%			
Rated Frequency	50 Hz 60Hz			
External Dimension (W*D*H)	1047*800*2450 mm	1247*800*2450 mm	1547*800*2450 mm	1847*800*2450mm
Working Zone Dimension (W*D*H)	787*560*700mm	987*560*700 mm	1287*560*700 mm	1587*560*700 mm
Power Supply Consumption	400 W	400 W	500 W	700 W
Inflow Velocity	0.3~0.8m/s			
LED fluorescent light power	24W	24W	30W	36W
Noise	≤70dB(A)			

Table 3



NOTE: a) Power supply consumption includes the consumption loaded at working zone, which should not exceed 500W.

b) Biolab reserves the right to make changes in future product design, without reservation and without

notification to its users.

2.1.5. Performance index

a) Vibration amplitude

The net vibration amplitude, at a range of frequency from 10 Hz to 10 KHz, would not exceed 5 μm (rms)

b) Electrical performance

Dielectric Withstand Test: the Fume Hood would not breakdown in 5s when the voltage increases by 1390V (AC) within 5s

c) Ground resistance $\leq 0.1\Omega$

2.2. Product Structure

2.2.1. Structural composition of Fume Hood




Figure 8




Figure 9

1. fan and the wind barrels 2. control panel 3. Socket 4. Water valve 5. Gas valve 6. Gas faucet 7. fuller faucet. 8. Base cabinet door 9. Table panel 10. Power Supply 11. Clump weight 12. Access wind

Note: BFSD-404 is a double fan.

2.2.2. Structure introduction

Lighting source

LED fluorescent lamps are used for lighting to ensure that the average illuminance in the operating area meets the standard requirements.

Control panel

The operation panel of this equipment has lighting, fan, fan speed control, socket, power touch key and windshield display unit, and various function work display lights. (See below)




Figure 10

2.3. Control Panel

2.3.1. Soft touch buttons

The key to human-computer interaction is that the main operations of the device can be performed with the touch of a button. Through the function of the buttons, the device can fully realize the most basic functions;



Power button, the main switch of the Fume Hood



Fan key: , is the control key of the working state of the fan (Figure 10). Each time it is pressed, the working state and corresponding indicating state change once.




Fan speed control key: , is the control key of the fan speed (Figure 10). Each time you press it, the working status and the corresponding indication status on the windshield display unit will change once, from F1 gear to F9 gear to adjust 9 gears in turn, when in F9 gear, press the fan speed adjustment button again, then return to F1 gear.

The current value is saved every time the fan gear is adjusted, that is, every time the fan is turned on, it

is the working gear when the fan was turned off last time.



Lighting key: , is the control key of the light (Figure 10). Each time you press it, the state of the lighting tube and the indication state of the corresponding indicator light change once, that is, from on to off, or from off to on;



Socket key: , is the control key for the power-on and power-off working state of the socket. Every time you press it, the on-off state and the corresponding indication state will change once; Socket

There are sockets on the left and right sides outside the operating area, which can supply power to the equipment in the operating area.



The power of the equipment used on the socket does not exceed 500W.

Structure

- a) PP board for the shell of the fume hood;
- b) The operating area of the fume hood is made of PP board, which has high acid and alkali resistance; the countertop is made of solid psychological board, which can be disassembled for easy cleaning.
- c) The window of the fume hood is made of 5mm thick tempered glass.
- d) The control panel adopts light-touch switches, which makes the machine beautiful in appearance and easy to operate.
- e) The electric control system has the functions of anti-overload, anti-electric shock, etc., with stable performance and long service life.
- f) The socket adopts special safety products for laboratories, and the material is PC material, which is fire-resistant and flame-retardant.

2.4 Points that should be paid attention to during normal operation and normal operation and use process

2.3.1 Things to pay attention to during operation

Precautions

- 1. The fume hood is an important laboratory safety equipment. In order to ensure safe and normal use, please read this instruction manual and precautions carefully, and participate in appropriate laboratory safety and skill training;
- 2. Please read this manual before use.
- 3. Please keep this manual for future reference.
- 4. The company will not be responsible for any damage caused by improper use or self-changing product structure.
- 5. The equipment should not be placed near doors/windows or corridors where people flow frequently.
- 6. The power socket has good grounding.
- 7. Please turn off the power and unplug the plug when replacing the lighting lamp.
- 8. The packaged cabinet should be stored in a well-ventilated environment with no corrosive gas and no relative humidity of more than 75%.
- 9. The perspective window of the fume hood is made of explosion-proof tempered glass. In order to keep the window clean, it should be wiped with a damp soft cloth and kept away from hydrofluoric acid;
- 10. Lining accessories such as the deflector of the fume hood should be cleaned regularly according to the usage;

11. The air ducts and fans of the fume hood should be regularly cleaned and maintained by special personnel;
 12. Do not place any equipment within 150mm from the glass window. Large-scale experimental equipment must have sufficient space and should not affect the flow of air.
 13. During the use of the equipment, do not put soft and delicate items (such as soft paper towels) on the table, so as not to be sucked into the negative pressure air duct and fan by the suction port and affect the operation of the equipment;
 14. The product shape and color of the company's pictures are subject to the actual product, and the product model is subject to change without prior notice.
 15. The temperature of the contact plate in the operation area cannot be higher than 100°C.
- Solemnly declare: The company will not be responsible for the risks caused by the operation not in accordance with the regulations!

2.3.2 Normal operation process

- a. Turn on the power, the screen lights up and displays three horizontal bars  .
- b. Press the power button on  the control panel, and the screen will display the working time of the fan (that is, the working time of the filter, the displayed value multiplied by 10 is the actual time, in hours), Open the glass door, press the fan button  , the screen will display the fan gear at the last working time, the fan indicator light above the fan button will light up, and the experiment operation in the cabinet can be carried out normally after 5 minutes of operation.
- c. Press the lighting button  , the lighting indicator above the button will light up, press the lighting button  again to turn off the lighting (you can choose to turn on or off the lighting according to the environment in the laboratory).
- d. After use, turn off the fan and lighting, and place the glass door to the bottom.
- e. If the power is cut off during use or the power plug falls off under abnormal circumstances, the device will maintain the last working state when the power is turned on again, that is, the working state before power-off when the power is turned on again , with power-down memory function.

2.4. Regular Maintenance

A detailed daily record of operating time is recommended, as the accumulated using time will directly affect the plan of maintenance.

Maintenance should be performed annually or every 1000 operating hours, and every restart.



NOTE: a) To avoid electric shock, please cut off ALL power before applying maintenance for the equipment!

b) The blower and the exhaust duct should be inspected and maintained regularly.

c) The accumulated operating time is a vital factor of deciding when the maintenance is needed. A comprehensive record of operation is highly recommended to be taken down after each time of

operating.

Maintenance methods

a) Weekly and monthly maintenance

- i. Clean the external surface and front window (refer to 2.4.1.)
- ii. Check the various functions of the Fume Hood
- iii. Record down the maintenance result

b) Annual maintenance

- i. Check the two lifting belt (sash) of the front window tubular motor, make sure both of them are well connected to the motor with same tightness
- ii. Check the fluorescent lamp, replace it if needed
- iii. Apply for overall performance test of the cabinet annually to ensure that the safety performance has met the requirements. User is responsible for testing costs
- iv. Record down the maintenance result

2.4.1 Surface cleaning

In order to keep the cabinet clean, please clean it regularly (recommended at least once a week). When wiping, please use a soft cloth dipped in soapy water and wring it dry. Please do not spray any chemical reagents on the operation panel or other labels to prevent discoloration of the label film or illegible writing. Clean the outer surface of the cabinet and the glass with a soft cleaner or a special cleaner for glass.

2.4.2 Storage conditions

Fume Hood should be stored in a warehouse with relative humidity no more than 75% and temperature lower than 40°C. The warehouse should have good ventilation performance without acid, alkali or other corrosive gases. Storage period shall not exceed one year. Fume Hood stored for more than one year needs to be unpacked and checked before selling and using. Only the tested and qualified safety cabinet could be sold.

2.5 Wiring Diagram




Figure 11

3. Trouble Shooting and Labels

3.1 Common Failures and Solutions

Before diagnosing the fault, please confirm whether the power supply is connected well, whether the power cord is obviously damaged, and whether the fuse is in good condition.

1. Check whether the equipment has implemented reliable grounding measures according to the requirements in the manual to ensure the safety of maintenance and use. Check whether the electrical wiring of the equipment is disconnected, disconnected or short-circuited. If there are similar situations, they should be eliminated in turn;

2. Judgment and solution of common faults

Failures	Checking Part	Suggestion
Fluorescent lamp fail to work	Fluorescent lamp tube	Replace the fluorescent lamp tube
	Circuit	Check the circuit
	Control panel	Replace the control panel
Button fail to work	Control panel	Make sure the power is well connected and the fuse is in good condition
		Check if the button is broken
		Make sure the connecting wire is well connected
		Replace the control panel
Fan fail to work	Fan	Replace the fan if it is defective
	Circuit	Check the circuit
	Control panel	Replace the control panel
No electricity in socket	Socket fuse	Check if the socket fuse is broken
	Socket	Check if the socket is broken
	Circuit	Check the circuit
	Control panel	Replace the control panel
No electricity in equipment	Power supply	Check whether the power supply is well connected
	Power cord	Check whether the power cord is in good condition
	Fuse	Check if the fuse is damaged
	Potential transformer	Check whether the transformer works normally
	Control panel	Replace the control panel
Display fail to work	Signal transmission line	Check whether the signal transmission line is well connected
	Display screen	Check whether the screen is in good condition
	Control panel	Replace the control panel

Table 4

3.1.1.Simple accessories replacement

a) Replace the fuse

The round safety tube on the top of the right panel of the fume hood cabinet is determined according to the label, which is F10A ϕ 5x20 mm. When it is necessary to replace the FireWire fuse, use a flat-head screwdriver to lift out the fuse holder and replace it, and then press it back.



Figure 12

b) Replace fluorescent light

When the lighting lamp of the fume hood is broken and needs to be replaced, press the glass door up button to raise the glass door to the highest height, then disconnect the power supply, remove the power cord of the LED fluorescent lamp, Use a screwdriver to remove the lampshade fixing device, remove the old LED fluorescent tube, replace it with a new one, and then connect the power cord of the LED fluorescent lamp, and the power-on test is correct.



Figure 11



Note:

- (1) The operation of the above electrical components must be carried out by a qualified electrician under safe conditions (cut off the power supply). All other parts are not allowed to be disassembled, otherwise the consequences shall be borne by the user;
- (2) When the equipment does not have the above failures and the operator cannot eliminate them immediately, please notify our company's maintenance department immediately. For your safety, please do not repair the equipment by yourself;

(3) The maintenance work of this equipment can only be undertaken by trained and approved technical personnel;

(4) If you need to order parts, you can find our technical service department, please indicate the model and serial number of the fume hood you purchased.

3.2 Label Description

3.2.1 Fuse label




Figure 13

Note: The 10A fuse indicator label is located under the live wire fuse

3.2.2. Ground label



Figure 14

4. Warranty

4.1. Warranty is 12 months from EX-factory date (excluding consumable accessories, Fluorescent lamp, fuse)

4.2. Biolab would not be liable for any repair of damage caused by improper operation

4.3. If the warranty has been expired, Biolab would still responsible for repair with relative charges

4.4. Life time of fume hood is 8 years from production date on the label

4.5. Biolab would provide equipment drawings and necessary technical data for maintenance companies or personnel trained by Biolab engineers.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com