



DESKTOP LAMINAR AIRFLOW

BJU-2301

INDEX

1. Description	2
2. Technical Parameters	3
3. Working Requirements	4
4. Panel Instruction	5
5. Wiring Diagram	8
6. Installation Conditions and Operating Environment	9
7. Installation	9
8. Common failure and solution	10

1. Description

Vertical Laminar Flow Hood / Cabinet (also called clean bench) provides sample protection as purified air travels across the working zone of the cabinet in a vertical, unidirectional stream and leaves the main work chamber across the entire open front of the cabinet.

Most Vertical Laminar Flow Hoods also have perforations at the back wall of the work zone which are designed to eliminate air turbulence and the possibility of dead air corners in the work zone.

Humanized Features :

Ergonomic Control Panel Location, centered and angled down for an easy reach & viewing.

Brushed stainless steel work surface is for easy cleaning.

The non-direct light source design can effectively relieve the operator's visual fatigue under the condition of ensuring the illumination.

Intuitive and easy-to-use microprocessor controller with LCD display, graphical and digital display of the operating status of each function, such as filter life, UV life etc.

Safety Features :

Front sliding sash interlock with fan and UV lamp.

UV sterilization with UV lamp timing function.

Front sash and side windows are made of UV-resistant tempered glass.

The host and sockets have independent capacitors to prevent current overload.

Convenience Features :

Memory function in case of power-failure.

HEPA filter and UV life indicators facilitate a timely replacement plan.

Two side glass windows installed for easy experiment observation.

Easy-to-change pre-filter installed could extend the HEPA filter's lifetime.

Equipped with a temperature sensor to show the real temperature during operation.

Waterproof sockets provide convenience for small instruments to be used in the operating area.

2. Technical Parameters

Parameter	Specification
Model	BJU-2301
External Size (WxDxH) (mm)	590x460x700
Internal Size (WxDxH) (mm)	480x340x370
Max Opening	310 mm
Control System	Microprocessor control with LCD display, display the operating status of each function

Airflow System	Airflow Velocity	0.2-0.6 m/s, intelligent multi-level adjustable
	Pre-Filter	1 pcs, Polyester fiber, washable
	HEPA Filter	1 pcs, Efficiency: >99.995% at particle size of 0.3 μ m
	Cleanliness Grade	Class 100 (ISO 14644-1 Class 5)
Cabinet Construction	Partial Air Circulation	Perforations at the back wall of the work zone which are designed to eliminate air turbulence and the possibility of dead air corners in the work zone.
	Main Body	Electro-galvanized steel with antimicrobial powder coating
	Work Surface	304 Stainless steel
	Front Window	Manual, 5 mm UV resistant tempered glass
Electrical System	UV Lamp	8Wx1
	LED Lamp	8Wx1
	Consumption	110W
	Power Supply	AC220V $\pm$ 10%, 50/60Hz; 110V $\pm$ 10%, 60Hz
Sound Emission		$\leq$ 62 dB
Colony Count		$\leq$ 0.5 pcs/vessel (petri dish is ϕ 90 mm)
Vibration		$\leq$ 3 μ m
Work Surface Bearing		10 kg
Standard Accessory		1 LED Lampx1; 2 UV Lampx1; 3 Pre-filterx1; 4 HEPA Filterx1
Net Weight		36 kg
Gross Weight		41 kg
Shipping Size (WxDxH) (mm)		710x590x850
Shipping Volume		0.36 m ³

Table 1

3. Working Requirements

1. Electricity: 220VAC $\pm$ 10% (or 110VAC $\pm$ 10% optional)
2. Power frequency: 50/60Hz;
3. Ambient temperature: - 10~60°C;
4. Relative humidity: < 90%RH (Non-condensing)

4. Panel Instruction




Figure 1

(1) Door Status Icons

When the door is open, the icon goes out in the middle and the outer frame lights up; When the door is closed, this icon is all lit.

(2) Ambient temperature

In the running state, the current ambient temperature measurement value is displayed; When the password is entered, it will be displayed as "PA";

In the parameter setting state, it is displayed as the English abbreviation of the corresponding parameter;

(3) Fan icons

When the fan is running, the icon lights up and turns;

(4) Fan Gear Icon

Displayed according to the fan setting gear;

(5) Time remaining icon

(6) Run time icon

In the running state, the cumulative running time of the fan is displayed; In UV sterilization state, the remaining time of sterilization is displayed; When UV sterilization is completed, it will display as "End"; When the password is entered, it will be displayed;

(7) Power Icon

When the power is off, the icon lights up;

(8) UV Sterilization Icon

When the system is in UV sterilization state, the icon flashes;

(9) Lighting icon

When the lighting is on, the icon lights up;

(10) Socket icon

When the power supply is turned on, the icon lights up;

(11) Filter lifespan scale bar

When the life is 0, it goes out;

When the life is not 0, the corresponding ratio (tenth ratio) is displayed;

(12) Sterilization lamp lifespan scale bar

When the lifespan is 0, it turns off;

Button Definition (two styles exist) :




Figure 2

[Power] key: The power switch key, click to turn the power on or off.

[Lighting] key: Lighting on/off key, click to turn on or off the lighting.

[Socket] key: Socket switch key, click to open or close the socket.

[Sterilize] key: Sterilization lamp switch key (when the door is closed, the lighting is closed, the socket is closed, and the fan is closed, press and hold for 3 seconds to turn on UV sterilization).

[Fan] key: Fan switch key (when the fan is off, the fan icon is off, and the fan gear bar is off; when the fan is on, the fan icon is displayed, and the fan gear bar shows the current gear).

[▲] [▼] Increase & Decrease keys(When the fan is off, increase or decrease the set value of the gear; when the fan is on, directly increase or decrease the gear; when the parameter is set or the password is entered, the parameter value or password is increased or decreased).

[Set] key: Click the setting button to enter the password input interface, modify the input password through the increase & decrease buttons, and then click the setting button to enter the corresponding parameter setting interface ; In the parameter setting interface, click the setting button to cycle the settings, and modify the parameters through the increase & decrease buttons. Press and hold the setting button for 3 seconds to exit and save the settings.

Interlock function description :

The fan can only be turned on when the door is open;

When the door is closed, turn off the lighting, turn off the socket, and turn off the fan at the same time to start the sterilization;

Sterilization time setting :

Click the setting button to enter the password input interface, which is displayed as "PA" ; Change the input password to " 1" by the increase & decrease keys, click the setting key to confirm, and it will display as "UT" at this time ; Modify the value to the value specified by the customer through the increase & decrease keys, and click the set key to exit.

5. Wiring Diagram




Figure 3

6. Installation Conditions and Operating Environment

Laminar flow cabinet shall be placed in a protection area, workspace can not be right opposite to the door or window, and away from the outlet of air conditioning, prevent from the ventilation system, air conditioning, door, window, and personnel movement caused by the airflow impact on Laminar flow cabinet.

At least 300mm gap must be kept in the side and back side of the Laminar Air Flow for clean operation and for inspection.

Working environment :

- 1) Only applicable for indoor operation.
- 2) Ambient temperature: 15°C ~ 35°C.
- 3) Relative Humidity: ≤ 75%.
- 4) Atmospheric pressure range: 70KPa ~ 106KPa.
- 5) Electrical parameters: Adequate power supply to the laminar flow cabinet. Power supply needs to be grounded; (Judging method: test the live wire and the neutral wire of the socket with multimeter. The voltage between live and ground should equal the voltage of the local electrical grid, and the voltage between neutral and ground should equal 0. Otherwise the power supply is not grounded correctly).

7. Installation

- a. Remove all the package materials
- b. Check the surface of the base cabinet and the main body to make sure there is no scratch, deformation or foreign bodies.
- c. According to the packing list of instructions, count carefully the accessories and materials.
- d. Position the base cabinet to the final location where an appropriate power supply is nearby. Move the whole equipment as close as possible to the final installation place.
- e. Connect the base stand and main body.

Inspection after Installation :

Check whether the voltage on the nameplate of the equipment is consistent with the actual supplied voltage, power on the equipment, turn on the display switch, and check the following items.

Checking Items	Normal working status
Power status	Equipment could be powered on/switched off when press the power button
Fan	Runs normally after pressing the Fan button; speed could be adjusted by pressing the adjusting button
LED lamp	Lamp lights up after pressing button
UV Lamp	Lamp lights up after pressing button
Socket	Use multimeter to test voltage output after pressing the socket button

Table 2

NOTE: Please contact the technical department or agent for inspection or trouble shooting when problems could not be solved. Methodology of trouble shooting is stated in the After-sale Service Manual.

8. Common failure and solution

Please confirm that the power is well connected, the cord is in good condition (without any damage) and the power lock is unlocked before trouble shooting the following problems.

Faults	Check Parts	Suggestion
LED light does not light up or abnormal	Lamp holder	Check that the lamp tube and the lamp holder are firmly connected
	Light tube	Replace the lamp
	Line	Check the line
	Control panel	Replace the control board
UV lamp does not turn on or abnormal	Interlock	Check the LED lamp and the blower is open or not
	Lamp holder	Tube and lamp holder is connected securely
	Circuit	Check circuit
	Ballast	Change it
	UV lamp	Change it
	Micro Switch	Check if Micro Switch is broken
	Control panel	Change it
Buttons do not work	Control panel	Make sure the power connects well and the fuse is well
	Control panel	Check if the button is broken
	Control panel	Make sure the connecting wire is connected well
	Control panel	Change control panel
Fan does not turn	Fan	The fan is damaged, replace the fan
	Micro Switch	Determine whether the micro switch is closed
	Line	Check the line
	Control panel	Replace the control board
The device is not powered	Power supply	Check power supply connects well
	Power wire	Check whether power wire has obvious damage
	Fuse	Check if the fuse is good
	Switch	Whether there is a problem with the power switch and terminal wiring is off or loose

The display does not light up	Transformer	Check whether the transformer works normally
	Control panel	Change it
	Connection winding displacement	Check whether the connection cable is in good contact
	Display screen	Check if the display is good
	Control panel	Replace the circuit board

Table 3

NOTE :

(1) The above trouble shooting methods should be done by qualified electricians under safe conditions(cut off power supply). Other components should not be removed.Risk caused by failing to follow those instructions would be responsible by the user.

(2) The trouble shooting and repair of this equipment only could be undertaken by trained and recognized technicians.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com