

Operation Manual



BONC 200 series

Natural Convection Oven

Thank you for Choosing Biolab products. Please read the “Operating Instructions” and “Warranty” before operating this unit to assure proper operation.

Index

| | |
|--|----|
| I. Summary..... | 03 |
| II. Structure features..... | 03 |
| III. Product structure diagram and parameters..... | 04 |
| IV. Working conditions..... | 06 |
| V. Attentions..... | 06 |
| VI. Operation instruction..... | 07 |
| VII. Meter operation instruction..... | 08 |
| VIII. Fault analysis..... | 12 |

01 Summary

BONC 200 Oven is widely used for drying, baking, melting, sterilizing and curing in labs of industrial enterprises, scientific research institutions, and health and medicine units etc.

02 Structure features

1. The body of the oven is made of high quality cold-rolled steel plate, finished with electrostatic powder spraying in pleasant design and durable using.
2. FW working room is made of stainless steel and foursquare semicircle transits. The heating mantle is of knock-down construction and easy to be cleaned.
3. PID digital and intelligent temperature controller has the functions of setting temperature, setting time to display in duel digital screen, over temperature alarming and time-setting.
4. Heater consists of two shifts: high temperature and low temperature. The latter one not only assures temperature reliability, but also saves the energy.
5. Dependent over temperature alarming equipment is synchronous with temperature control meter of auto switch and over-heat alarming.
6. Loose & Tight adjustable lock can adjust gate lock and ensure well leak tightness.
7. Now forced air ensures temperature evenness and prevents tiny and powdery substances blowing away.

03 Product structure diagram and parameters

i. Main technical parameters


| Model | | BONC-201 | BONC-202 | BONC-203 | BONC-204 |
|---------------|---------------------------------|--|-------------|-------------|--------------|
| Cycle Mode | | Natural convection | | | |
| Function | Temp. Range | RT+10-300℃ | | | |
| | Temp. Resolution Ratio | 0.1℃ | | | |
| | Temp. Motion | ±1℃ | | | |
| | Temp. Uniformity | ±3.5% | | | |
| Structure | Inner Chamber | Mirror Stainless Steel | | | |
| | Outer Shell | Cold rolling steel electrostatic spraying exterior | | | |
| | Insulation layer | Compound silicate heat preservation board | | | |
| | Heater | Stainless steel heater | | | |
| | Power rating | 0.8kW | 1.2kW | 1.6kW | 2.5kW |
| | Exhaust hole | φ28mm top (with function of test hole) | | | |
| Controller | Temp. control mode | Two temperature section PID intelligent | | | |
| | Temp. setting mode | Touch button setting | | | |
| | Temp. display mode | display on LED | | | |
| | Timer | 0-9999min (with timing wait function) | | | |
| | Operation function | Fixed temperature operation, timing function, auto stop. | | | |
| | Additional function | Mechanical independent temperature limiter, sensor deviation correction, temperature overshoot self-tuning, internal parameter locking, power-off parameter memory | | | |
| | Sensor | pt100 | | | |
| Specification | Inner Chamber size (W*L*H) (mm) | 310*310*310 | 350*350*350 | 400*360*450 | 500*450*550 |
| | Exterior size (W*L*H) (mm) | 460*510*695 | 500*550*735 | 550*550*840 | 636*680*915 |
| | Packing size (W*L*H) (mm) | 550*585*845 | 590*625*885 | 640*635*985 | 740*725*1085 |
| | Volume | 30L | 45L | 65L | 125L |

| | | | | | |
|---------------|-------------------------|-----------------|-----------------|-------------|------------------|
| | Shelf number | 6 | 7 | 9 | 13 |
| | Load per rack | 15kg | | | |
| | Shelf space | 35mm | | | |
| | (50/60HZ)Current rating | AC220V/ 3.6A | AC220V/ 5.5A | AC220V/7.2A | AC220V/ 10.5A |
| | NW/GW (kg) | 33/37 | 37/43 | 44/49 | 60/66 |
| Access ory | Shelf | 2 | | | |
| | Shelf frame | 4 | | | |

ii. Structure diagram



iii. Temperature profile



NOTE: The heating time is different between each model.

04 Working Conditions

The drying oven work under the following conditions:

1. Temperature ranges between 5~40°C;
2. Relative humidity less than 85% RH;
3. Power: voltage 220-240v, frequency 50-60Hz;
4. No succession and corrosive gas surround the oven.

05 Attentions

1. Install the outer ground protection to ensure safety of machine and experiment; supply power as the machine nameplate required.
2. This equipment is forbid to use in inflammable and explosive, poisonous and strong corrosive experiments.
3. Make sure horizontal installation.
4. Non-professionals are not allowed to disassemble and repair this machine.
5. Pay attention to the setting temperature when dealing with inflammable matters.
6. Make sure dry the resin container, if the temperature is setting too high by accident, the container would be dissolved and then fall on the heater, which will

cause fire.

7. Overfilled of sample will lead to overheat of working room under part, which will dissolve the inflammable material and cause fire.

8. While the machine is working, don't touch the device top, as well as observation window and exhaust port to keep away from high-temperature burns.

9. Read the instruction book before operation.

06 Operation Instruction

1. Put the material needs drying into container (advice: size of drying material should not over 2/3 of the shelf); then close the container door and switch power, and next switch on the blower.

2. Heating

Set the temperature as needs (find details in meter instruction), then the temperature starts to rise; when temperature inside working room reaches the set point, the indication light will go out, after constant temperature for 30min, the working room goes into constant temperature state.

3. Working time:

Decide the drying time according to humidity of sample.


Note: for example, if the sample humidity is big, the sample on each layer should not be too thick to ensure intensive drying of sample.

4. After finishing drying, turn off power, and then bring the sample out.

5. Keep the drying oven clean, wipe the container sealing rubber strip by soft cloth and clear the dirt out; avoid cleaning it by chemical solution to prevent chemical reaction damage on sealing rubber strip.

6. If the oven is unused for a long time, daub neutral grease or Vaseline on galvanized parts to prevent corrosion; cover the oven with plastic dust cap, and store it in the dry room to keep the electric device against wet.

07 Meter Operation Instruction



Indicator light function

- 1) AT: It flickers during self-tuning, it is not bright in any other state.
- 2) OUT: It is lit when heating output.
- 3) TIME: It is lit when time is set, it flickers in the process of timing.
- 4) ALM: It is lit when there is a temperature alarm.
- 5) SW: It is invalid.

Button function

- 1) [SET]: In normal state, press this button to enter the setting state.
- 2) [←/AT]: "SHIFT" button. In the setting state, click this button to shift the set value.
In normal state, press this button for 6 seconds to enter the auto-tuning selection state.
- 3) [→/RES]: "DEC" button. In the setting state, click this button to reduce the set value.
If you keep pressing this button, the value will reduce continuously. In the normal state, when the timer ends, press this button for 3 seconds, the controller will restart to work.
- 4) [▲]: "INC" button. In the setting state, click this button to increase the set value. If you keep pressing this button, the value will increase continuously.

1. Operation and using

1-1. When the controller is switched on, All displays light up for 2 seconds, display windows show the version number and controller model for 2 seconds, then it starts running.

1-2. Temperature and Time Setting

1) Without Timing Function :

In the normal state , press the “SET” button to enter the temperature setting state, windows display the prompt “SP” and the temperature set point value. Using the “SHIFT”、“DEC” and “INC” buttons, user can edit the temperature set value. Press the “SET” button again, the controller will return to its normal state, the setting value will be saved automatically.

2) With Timing Function :

In the normal state , press the “SET” button to enter the temperature setting state, windows display the prompt “SP” and the temperature set point value. Re-press the “SET” button to enter the time setting state, windows display the prompt “ST” and the time set point value. Press the “SET” button again, the controller will return to its normal state, the set values will be saved automatically.

When the time is set to "0", it indicates the timer is inoperative, the controller will run continuously, the under window will display the temperature set point value. If there is time set, the under window will display the running time, its decimal point and the “TIME” indicator are lit, when the timer starts, its decimal point and the “TIME” indicator flickers. When the timer ends, the under window will display the “End” prompt, the buzzer will sound for 5 minutes, it can be muted by pressing any button, press the “DEC” button for 3 seconds, the controller will restart to work .

1-3. If the upper window show the prompt “---”, it indicates that the temperature sensor has faults or temperature exceeds the measuring range or the controller itself is faulty, the controller will cut off the heat output automatically, the buzzer will sounds continuously, “ALM” indicator is lit, Please check the temperature sensor and its wiring carefully.

1-4. When over temperature alarm, the buzzer beeps continuously, "ALM" indicator is lit, the heat output is cut off. If the over temperature alarm is caused by the change of the temperature setting value, "ALM" indicator is lit, but the buzzer does not beep.

1-5. When the buzzer sounds, press any key to mute.

2. Auto-tuning

In the normal state, press the “SHIFT” button for 6 seconds, the controller will enter the auto-tuning selection state, the upper window displays the prompt “AT”, the under window displays "0", change "0" to "1" by pressing the “INC” button, then press the “SET” button, the controller will run the auto-tuning program, the

“AT” indicator flickers. After auto-tuning end, the indicator stops flickering, PID parameter value is saved automatically. In the auto-tuning process, press the “SHIFT” button for another 6 seconds, the controller will stop the auto-tuning program.

During the Auto-tuning process, if over temperature alarm, the buzzer does not beep, “ALM” indicator is not lit, the heat output will be cut off, the “SET” button is invalid, the under window always displays temperature set point value.

3. Internal parameters settings

In the normal state, press the “SET” button for 3 seconds, windows will display the prompt “Lc” and the password value. Adjust the password to the required value, then press the “SET” button again, it will enter the internal parameters setting state. Press the “SET” button for another 3 seconds, it will return to the normal state, the set value will be saved automatically.

Parameter table 1

| Prompt | Name | Function description | (Setting range) Factory value |
|------------|----------------------------------|---|----------------------------------|
| Lc | Password key | When “Lc=3”, enter the next parameters. | 0 |
| ALH | Over-temp alarm | If “PV>SV+ALH”, the ALM indicator turns on. The buzzer sounds and the heat output turn off. | (0~100.0°C) 20.0 |
| P | Proportional band | Adjustment of proportional function. | (0~300.0°C) 35.0 |
| I | Integration time | Adjustment of integration function. | (1~2000S) 300 |
| D | Differential time | Adjustment of differential function. | (0~1000S) 200 |
| T | Control cycle | The temperature control cycle. | (1~60S) |
| Pb | Temperature deviation correction | It is usually used to correct errors in low temperature measurement. | (-50.0~50.0°C) 0 |
| PL | Temperature slope correction | It is usually used to correct errors in high temperature measurement. | (-999~999) 0 |

Parameter table 2

| Prompt | Name | Function description | (Setting range) Factory value |
|------------|------------------------------|--|----------------------------------|
| Lc- | Password key | When "Lc=9", enter the next parameters. | 0 |
| doT | Temperature decimal point | 0: No decimal point display 1: With decimal point display | (0~1) 1 |
| ndT | Timer mode | 0: No timer function. 1: Start timing when the temp reaches the set value. 2: Start timing as soon as the controller starts working. | (0~2) 1 |
| Hn | Timer unit | 0: Minute. 1: Hour. | (0~1) 0 |
| SPd | Timer parameter | If "ndT=1", Start timing when "SV-SPd≤PV≤SV+SPd" | (0.1~50.0°C) 0.5 |
| EH | Timer end mode | 0: Continue to control the temperature 1: Stop temperature control | (0~1) 0 |
| oPn | Door parameter | Automatic judge door opening. 0: invalid; 0: valid | (0~1) 0 |
| nP | Power percentage | Percentage of max heating power output. | (0~100%) 100 |
| Co | Heating prohibited deviation | When "PV≥SV+Co", heating output will be cut off | (0~50.0°C) 50.0 |
| SPH | Max set value | The maximum temperature set point value. | (0~400°C) 300.0 |


Parameter table 3

| Prompt | Name | Function description | (Setting range) Factory value |
|------------|---------------|---|----------------------------------|
| Lc | Password key | When "Lc=567", enter the next parameters. | 0 |
| rST | Factory reset | 0: cancel; 1: confirm | (0~1) 0 |

6. Wiring

Represents the charged output, should be directly connected to the load.

Represents the switch output without charge.



08 Meter Operation Instruction

| Phenomena | Causation | Treatment Method |
|---|--|--|
| 1.No power supply | 1.Plug is poor contact or line broke | 1. Connect the plug and line. |
| | 2. Fuse protector is broke. | 2. Change the fuse protector. |
| 2. No temperature rising inside container | 1. Low setting temperature | 1. Readjust and set tempe. |
| | 2. Heater is broke. | 2. Change the heater |
| | 3.Temperature controller is broke | 3. Change the temperature controller |
| | 4. Temperature sensor is loose. | 4. Screw up the sensor nut. |
| | 5. Temperature sensor is broke | 5. Change the temperature sensor. |
| 3. No temperature | 1. Setting temp. of Detached tem. limiter is low | 1. Readjust the temp. 30°C above setting temp. |

| | | |
|--|---|---|
| rising alarm | 2. Detached temperature limiter sensor is broke. | 2. Change the detached temperature limiter sensor |
| 4. Temperature cannot reach the setting point. | 1. Exhaust port is fully opened | 1. Shut off the exhaust port. |
| | 2. The container is overfilled, hot air cannot convert. | 2. Decrease amount of sample to improve convection condition. |
| 5. The fan does not work. | The fan motor is broke | Stop work and check electric capacity and motor |
| 6.Displaying---- --- | The sensor is broke | Change the sensor |
| 7.Display STOP | Time-up | Press the program key for 3s to start |



BIOLAB SCIENTIFIC LTD.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8 Canada

Email: contact@biolabscientific.com Tel: +1 707 533 1445

Website: www.biolabscientific.com