



## PATHOLOGY WORKSTATION BAIP-301

## PATHOLOGY WORKSTATION BAIP-301

Workstations are used during dissections, autopsies, necropsies, and other pathological procedures. It includes grossing/trimming stations and dissection tables. It may also include a built-in sink and ventilation system. Used in Hospital pathology department, Pathology laboratory. Also known as Pathological sampling bench.

## BAIP-301 PATHOLOGY WORKSTATION



- Odor-remove function
- Single water sink, two-person operations
- Exhaust air system: High-power suction device, low noise
- Optional water heater is available, realizing hot & cold water adjustment
- The base cabinet is inward-concave type, convenient for sitting operation
- Optional pulverizer can completely smash the waste during test, to avoid blocking drainage pipes
- Automatic flushing device, slightly right-side inclined worktable, easy for waste discharge and cleaning

## SPECIFICATIONS

Model	BAIP-301
Work Surface Height	850 mm
Exhaust Airflow Volume	500~600m <sup>3</sup> /h
Main Body Material	304 stainless steel
Work Table Material	304 stainless steel
Package Size (WxDxH)	1160x890x2170 mm
Gross Weight	250 kg
Working Area Size (WxDxH)	895x430x900 mm
External Size (WxDxH)	1000x750x1950 mm
Consumption	250 W
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 60Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)