



DOUBLE BASKET TISSUE PROCESSOR BHTP-306

DOUBLE BASKET TISSUE PROCESSOR BHTP-306

Tissue processor prepares tissue samples for sectioning and microscopic examination by fixing, staining, dehydrating or decalcifying them. It is mostly single unit devices which can accommodate a variety of processing techniques therefore improving the efficiency of tissue processing.

Used in Clinical and Research Histopathology.

Also known as Tissue Embedding Cassettes.

BHTP-306 DOUBLE BASKET TISSUE PROCESSOR



This machine can run both single hanging basket and double basket, one dehydrator is equivalent to two hanging baskets. Large LCD touch screen, all English menu, simple and clear operation. Adopt intelligent program control technology to protect equipment and specimens. Intelligent separate protection specimen in case of failure. When the dehydrator fails before the anhydrous alcohol cylinder, the tissue is protected in the fixed cylinder. when the dehydrator fails before the paraffin cylinder and after the anhydrous alcohol cylinder, the tissue is protected in the anhydrous alcohol cylinder; when the hanging basket fails in the paraffin cylinder dehydrator, the tissue is protected in the paraffin cylinder. Adopt new water heating mode, low pressure contactless temperature control technology, more safe. Reagent cylinder can be fully heated to avoid the effect of dehydration due to the difference of temperature in winter and summer.

SPECIFICATIONS

Model	BHTP-306
Volume of individual reagent containers	2000 ml
Number of reagent containers	14
Hanging basket filtrate	Time can be adjusted arbitrarily
Number of dehydrated tissues	More than 200
Processing time in a single container	0 ~ 99 hours adjustable
Power supply	AC 220 V \pm 10 % 50 Hz
Power	\leq 1500 W



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com