



TISSUE BAKING PROCESSOR BHTP-607

TISSUE BAKING PROCESSOR BHTP-607

Embedding is a fixation process during which tissues or the specimens of interest are entrapped in the embedding media. The choice of the media depends on type of tissue, type of microtome and types of microscope to examine the tissue.

Used in Pathology, Histology.

Also known as Tissue Embedding Station.

BHTP-607 TISSUE BAKING PROCESSOR

Heat conduct fast, durable. Self memory of temperature. Display presetting temperature & current temperature. Corrosion-resistance of the bath.

SPECIFICATIONS

Model	BHTP-607
Range of setting temperature	Ambient to 90 °C (constant adjusting).
Temperature precise error	+ 2 °C.
Allowance ambient temperature	0 - 40 °C
Working voltage	220 V / 50 Hz
Power	250 W
Dimensions overall	450x400x200 mm (WxDxH)
Weight	7 kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com