



MANUAL OR TABLE BOPT-203

MANUAL OR TABLE BOPT-203

Manual Operation table are suitable for various operation tables conducted by medical unit. The OR table and accessories are ergonomic and can be used across multiple disciplines. Various actions are controlled by hydraulic pneumatic system. 360° rotation helps every surgical requirement.

Used in Medical unit, Surgery, Operation theatre, cardiovascular, pediatric, gynecology, abdominal, gall bladder, plastic surgery and more.

Also known as Manual surgical table, Manual OR Table, Manual operating table.

BOPT-203 MANUAL OR TABLE



Various actions are controlled by hydraulic pneumatic system

All movements of the table are manipulated using the operating handle under the head plate

Imported Y type sealing ring is durable

The table surface adopts quality carbon plastic plate which is easily for conducting X-ray examination, equipped with coincidence rail

Both sides of the bed rail for the high quality 304 stainless steel

The base adopts imported carbon cloth which can be easily cleaned and sterilized

The mattress is wrapped with imported artificial leather. High-density memory sponge is formed in one piece, with no seam

Foot treadle brake device which is easily for moving operating table

SPECIFICATIONS

| | |
|-----------------------|--------------------|
| Model | BOPT-203 |
| Table Length | 2000±50 mm |
| Table Width | 500±20 mm |
| Height Adjustment | (700~950) ±50 mm |
| Left tilt | ≥15° |
| Right tilt | ≥15° |
| Trendelenburg | ≥20° |
| Reverse trendelenburg | ≥15° |
| Head board upwards | ≥30° |
| Head board downwards | ≥90° |
| Backboard upwards | ≥70° |
| Backboard downwards | ≥15° |
| Waist board lifting | ≥70° |
| Leg board downwards | ≥90° |
| Leg board outwards | ≥180° |
| Main Voltage AC | 220V±22V, 50Hz±1Hz |
| Input Power VA | 150±10% |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com