




## FLASH POINT AND FIRE POINT TESTER BPTL-205

## FLASH POINT AND FIRE POINT TESTER BPTL-205

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

### BPTL-205 FLASH POINT AND FIRE POINT TESTER



Use heat resistant glass interlayer on the heater, avoid oil splash and fire breaking out. Automatic stop if there isn't operation after heating 20 minutes. Automatic air cooling when temperature is too high.

### SPECIFICATIONS

Model	BPTL-205
Air cooling	Automatic when temperature is too high
Stop	Automatic if there isn't operation after heating 20 minutes
Heating	Can only be started if press switch 3 seconds, avoiding touch by mistake.
Safety	Use heat resistant glass interlayer on the heater, avoid oil splash and fire breaking out



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)