



## THERMO SHAKER INCUBATOR BSTH-205

# THERMO SHAKER INCUBATOR BSTH-205

Used in Incubation, Cultivation, Re-suspending Pellets, Denaturin, PCR preparations, ELISA Assays, Colorimetric Assays, Enzyme, Protein Analysis.

Also known as Thermo Shaker, Lab Shaker, Thermo Shaker Incubator, Laboratory Thermo Shaker Incubator.

## BSTH-205 THERMO SHAKER INCUBATOR



Optimized 2D mixing designed for simultaneous mixing and vortexing applications  
Pre-programmed soft keys and reliable universal holder  
Automatic imbalance detection and anti-spill technology  
Mixing in a 3-in-1 format: Plates (up to 96-well and 384-well formats), Micro-Tubes, Vortexing numerous tube formats  
Custom blocks available for different experiment needs

## SPECIFICATIONS

|                        |                     |
|------------------------|---------------------|
| Model                  | BSTH-205            |
| Speed Range            | 300~2000 rpm        |
| Temperature            | RT+5-100°C          |
| Temperature Accuracy   | 0.5°C               |
| Temperature Resolution | 0.1°C               |
| Vortex Shaking Speed   | None                |
| Speed Increment        | 10 rpm              |
| Shaking Orbit          | 3 mm                |
| Time of Heating        | ≤12 min             |
| Time of Cooling        | No cooling function |
| Timing Range           | 1min ~ 99h59 min    |
| Motion                 | Orbital             |
| Overall Dimension      | 230x170x165 mm      |
| Weight                 | 5.5 kg              |
| Power                  | 150 W               |
| Power Supply           | 230 V, 60 Hz        |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)