

# PRODUCT CATALOG



## **TUBE ROLLER BTUR-101**





### **TUBE ROLLER BTUR-101**

The versatile, end-over-end tube roller and Rotators are designed for mixing samples in both a vertical and horizontal plane. They can be placed in an environment temperature ranges from 5 to 40° and the rotation speed is adjustable, from 0 to 99 rpm. Tube Rollers, Rotators and Rocker are designed to deliver salient performance, accuracy, and reproducible results for all your application needs.

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Tube Roller, Roller, Laboratory Roller.

#### **BTUR-101 TUBE ROLLER**



Rocking and rolling motions for complete mixing

Speed range of 0-70rpm

Six roller design with small footprint and easy clean

DC brushless motor provides reliable and quiet running, maintenance free

Designed for continuous operation

Can be used in cold rooms or in incubators

#### **SPECIFICATIONS**

Model	BTUR-101
Speed Range	0-70 rpm
Temperature	4-40°C
Roller Lenght	343 mm
Maximum Diameter Tubes	3 inches
Amplitude	24 mm
Maximum Load Capacity	5 kg
Operation	Continuous, No timer
Dimension	533Lx279Wx127H mm
Motor Type	DC brushless
Weight	5.5 kg
Power Supply	100-220V, 50/60Hz



#### Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada Email: contact@biolabscientific.com | Website: www.biolabscientific.com