



## CHILDREN MEDICAL BED

## CHILDREN MEDICAL BED

Children medical Bed built to the most current standards while also meeting the clinician's needs which allows easy care and disinfection. The bed is equipped with Trendelenburg and reverse Trendelenburg positions. Simple and modern design for safety and protection of the child.

Used in Medical unit, Nursing.

Also known as Children Hospital Bed, Medical bed, Hospital bed, Medical child bed, Pediatric bed, Medical pediatric bed, Hospital pediatric bed, Children medical bed, Children pediatric hospital bed.

### BHBD-104 MANUAL 3 FUNCTION CHILDREN MEDICAL BED



ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, pink, with patient name card holder, option: cartoon sticker  
Full length powder coated side rail, drop type with safe lock  
Bed sheet top made from perforated steel, epoxy powder coated and baked finish  
Bed framework made from mild steel ERW rectangular tubes  
Central locking ABS foot pedal with casters dia. 5"  
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

## SPECIFICATIONS

Model	BHBD-104
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	480~720 mm (without mattress)
IV pole	Yes
Dimension	1980Lx980Wx480~720 mm

### BHBD-105 ELECTRIC 3 FUNCTION CHILDREN MEDICAL BED



ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, pink, with patient name card holder, option: cartoon sticker  
Full length powder coated side rail, drop type with safe lock  
Bed sheet top made from perforated steel, epoxy powder coated and baked finish  
Bed framework made from mild steel ERW rectangular tubes  
Central locking S.S foot pedal with casters dia. 5"  
With 3 motors from Taiwan  
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

## SPECIFICATIONS

Model	BHBD-105
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	480~720 mm (without mattress)
IV pole	Yes
Dimension	1980Lx980Wx480~720 mm

## BHBD-107 CHILDREN BED WITH TWO FUNCTION



## SPECIFICATIONS

Model	BHBD-107
Load Bearing	150 kg
Dimension	1770x900 mm
Packing Dimension	wooden carton: 1880x980x1450 mm (assembled)
carton	
1880x1020x310 mm (detached)	Net Weight



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)