



MANUAL HOSPITAL BED

MANUAL HOSPITAL BED

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance

Used in Recuperating, Nursing.

Also known as Medical bed, hospital bed, Hospital cot, Electric ICU bed.

BHBD-509 TWO CRANK MANUAL BED



Steel framework, top class ABS & PP plastic

Collapsible Aluminium alloy protective railing; Detachable PP head & foot board

4 pcs universal castor, dia. 125 mm, Individual brake

A shoe rack under the bed board

SPECIFICATIONS

Model	BHBD-509
Function	2
Back-rest	0° - 75°
Knee Gatch	0° - 45°
Height (from bed board to ground)	470 mm
Castor	360° swivel castors (Dia. 150 mm), with individual brake.
Material	Powder coated steel bed frame, top class ABS & PP plastic
Bed board dimension	1910x830 mm
Standard accessory	IV pole
Net Weight	65 kg
Patient surface	1910Lx830W mm
Height adjustment	470 mm
IV pole	Yes
Dimension	2080Lx980Wx470 mm

BHBD-510 3 CRANK MANUAL BED



Steel framework, top class ABS & PP plastic

PP head & foot board and bed board, detachable and easy to clean 6 columns collapsible Aluminium alloy protective railing

4 pcs universal castor, dia. 125 mm, Individual brake

SPECIFICATIONS

Model	BHBD-510
Function	3
Back-rest	0° - 65°
Knee Gatch	0° - 45°
HI - LOW	750 mm - 575 mm
Castor	360° swivel castors (Dia. 150 mm), with individual brake.
Material	Powder coated steel bed frame, top class ABS & PP plastic
Bed board dimension	1930x830 mm
Standard accessory	IV pole
Overall dimension	2160Lx970Wx575-750 mm
Patient surface	1930x830 mm
Height adjustment	750-575 mm
IV pole	Yes
Dimension	2160x960x530- 730 mm

BHBD-511 4 CRANK MANUAL BED



Steel framework, top class ABS & PP plastic

PP material head / foot board, detachable and easy-clean

4 pcs drop-down PP protective railing, controlled by gas spring

Central brake castors, dia. 150 mm. Brake pedal locate below the foot board

SPECIFICATIONS

Model	BHBD-511
Function	5
Back-rest	0°-60°
Knee Gatch	0° - 30°
Trendelenburg	0 - 7°
Reverse trendelenburg	0 - 7°
HI - LOW	665 mm - 505 mm
Castor	360° swivel castors (Dia. 150 mm), with central brake under the foot board.
Material	Powder coated steel bed frame, top class ABS & PP plastic
Bed board dimension	1930x990 mm
Standard accessory	IV pole
Net Weight	117 kg
Overall dimension	2120x1055x550-665 mm
Height adjustment	755-485 mm
IV pole	Yes
Dimension	2200Lx1050Wx485-755 mm

BHBD-512 FLAT MEDICAL BED



ABS head & foot boards can be removable. with bumper, 8kg/pair, color :blue, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

4 pcs Individual lockable Taiwan Secure castors at 4 corners, diam 5"

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Shoe shelf, mosquito net rack, dinner table

SPECIFICATIONS

Model	BHBD-512
Back-rest lifting angle	0°-75°
Foot-rest lifting angle	0°-35°
IV pole	Yes
Dimension	2160Lx980Wx500H mm

BHBD-515 MANUAL 2 FUNCTION MEDICAL BED



ABS head & foot boards can be removable. with bumper, 8kg/pair, color option: brown, blue, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock.

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

With casters dia. 5" with individual brake

Manual 2 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Shoe shelf, mosquito net rack, dinner table

SPECIFICATIONS

Model	BHBD-515
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2150Lx980Wx500H mm

BHBD-516 MANUAL 3 FUNCTION MEDICAL BED



ABS head & foot boards can be removable. with bumper, 7kg/pair, color blue, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

With casters dia. 5" with individual brake

Manual 3 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Shoe shelf, mosquito net rack, dinner table

SPECIFICATIONS

Model	BHBD-516
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	430~720 mm
IV pole	Yes
Dimension	2120Lx980Wx430~720 mm

BHBD-517 MANUAL 2 FUNCTION MEDICAL BED



ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder

2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes

Central locking ABS foot pedal with casters dia. 5"

Manual 2 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

Optional Accessories: Shoe shelf, mosquito net rack, dinner table

SPECIFICATIONS

Model	BHBD-517
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2160Lx980Wx500 mm

BHBD-519 MANUAL 2 FUNCTION MEDICAL BED



ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder

ABS side rail, 4-piece, tuck away type with safe lock. with angle display

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes.

Central locking both side foot pedal with casters dia. 5"

Manual 2 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

SPECIFICATIONS

Model	BHBD-519
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2220Lx1060Wx500 mm

BHBD-520 MANUAL 3 FUNCTION MEDICAL BED



ABS head & foot boards can be removable.

with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder

ABS side rail, 4-piece, tuck away type with safe lock. with angle display

Bed sheet top made from perforated steel, epoxy powder coated and baked finish

Bed framework made from mild steel ERW rectangular tubes.

Central locking both side foot pedal with casters dia. 5".

Manual 3 functions adjusted by hidden type hand crank

Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes

SPECIFICATIONS

Model	BHBD-520
Function	3
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
Height adjustment	450~740 mm
IV pole	Yes
Dimension	2220Lx1060Wx450~740mm

500 MANUAL HOSPITAL BED



SPECIFICATIONS

Model	BHBD-522	BHBD-525
Load Bearing	240 kg	
Patient surface	2150x950 mm	1950x830 mm
Height adjustment	500 mm (from floor)	-
Back-rest lifting angle	-	0°-75°
Foot-rest lifting angle	-	0°-35°
Dimension	-	2200x1020x500 mm



BHBD-522



BHBD-525

BHBD-524 MANUAL HOSPITAL BED



Material: ABS

4 pcs of European style guardrail, upside down structure

Luxury central control castors, brake the bed with one stamp

With 4 bumpers

SPECIFICATIONS

Model	BHBD-524
Back-rest lifting angle	0°-75°
Foot-rest lifting angle	0°-35°
Dimension	2280x1120x420-780 mm

BHBD-526 MANUAL BED FOUR CRANK 5 FUNCTION



The main bed and auxiliary bed are connected with each other, auxiliary table is easy for dismantling. The auxiliary table with four wheels

The main bed adjusts various movements for the gears and auxiliary bed is ratched adjustment

Various Accessories are complete

SPECIFICATIONS

Model	BHBD-526
Table Length	1900+50 mm
Table width	600+20 mm
Height Adjustment	750+50 mm
Back board upwards	≥50°
Backboard downward	≥15°
Seat board downwards	≥8°
Seat board upwards	≥30°
Leg board upwards	≥30°
Leg board downwards	≥30°
Standard accessory	A pair of arm support, a pair of leg support, a pair of arm rods, a pair of foot pedal



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com