



## STRECHER

## STRECHER

Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Recuperating, Nursing.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

### BHBD-810 PATIENT STRETCHER



## SPECIFICATIONS

|           |                 |
|-----------|-----------------|
| Model     | BHBD-810        |
| Dimension | 1860x590x700 mm |

### BHBD-812 HOSPITAL TRANSFER STRETCHER TROLLEY



It is a type of manual transfer stretcher. Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down;

Aluminium-Alloy protective rails, can be folded and hidden under the surface on no use;

Central-lock castors, diameter150mm. Matched with direction castor, help the operator controlling the direction easily;

Rotate the crand and handle on bedboard, patient on the mattress can be moved parallellly to another bed or flat table

## SPECIFICATIONS

|                   |            |
|-------------------|------------|
| Model             | BHBD-812   |
| Height adjustment | 920-610 mm |
| Net Weight        | 89 kg      |

## BHBD-815 MANUAL TRANSFER STRETCHER



Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down Bright arrow symbols direct the users how to operate clearly

Central-lock castors, diameter 150mm. Matched with direction castor under the chassis, which help users easy control the stretcher direction

Standard accessories: Mattress, I.V.Pole

### SPECIFICATIONS

|                                    |  |
|------------------------------------|--|
| Model                              | BHBD-815   |
| Back-rest                          | 0 - 85°  |
| Height range (bed board to ground) | 862 - 566 mm                                       |
| Standard accessory                 | Mattress, IV pole                                  |
| Bed Board                          | 1906x610 mm  |
| Castor                             | Central lockable 360° swivel castors (Dia. 150 mm) |
| Overall dimension                  | 2013x700 mm  |
| Net Weight                         | 82 kg  |

## BHBD-816 MANUAL TRANSFER STRETCHER



Single crank for height adjustment Turn by clockwise, stretcher will be up; Turn by anti- clockwise, stretcher will be down

Bright arrow symbols direct the user how to operate clearly

Central-lock castors with steering, diameter 150mm

Concave plastic chassis cover, for storage use

Rotary rollover PP protective railing is blew molding by one-shot plastic forming, 3 different positions of the side rail supply different use

PP bed board produced by blowing molding, easy for clean

Standard accessories: Mattress, I.V. pole

### SPECIFICATIONS

|                                    |  |
|------------------------------------|--|
| Model                              | BHBD-816   |
| Back-rest                          | 0 - 89°  |
| Height range (bed board to ground) | 791 - 509 mm                                       |
| Standard accessory                 | Mattress, IV pole                                  |
| Bed Board                          | 1970x685 mm  |
| Castor                             | Central lockable 360° swivel castors (Dia. 150 mm) |

|                   |               |
|-------------------|---------------|
| Overall dimension | 1970x685 mm   |
| Net Weight        | 67 kg         |
| Dimension         | 1960Lx680W mm |

## BHBD-817 CONNECTING TRANSFER STRETCHER



It is a manual Connecting Transfer Stretcher, can be used to defend cross-infection  
 Cart A inside the operation room, Cart B outside, adjusting them to the same height,  
 then push Cart B to Cart A, two carts will be connected, following the orbit, bedboard C  
 can be pushed from one cart to another, and VICE VERSA  
 Bright arrow symbols direct the users how to operate clearly  
 Standard accessories: Mattress, I.V. pole

### SPECIFICATIONS

|                                    |  |
|------------------------------------|--|
| Model                              | BHBD-817   |
| Back-rest                          | 0 - 85°  |
| Height range (bed board to ground) | 660 - 910 mm                                       |
| Standard accessory                 | Mattress, IV pole                                  |
| Bed Board                          | 1906x610 mm  |
| Castor                             | Central lockable 360° swivel castors (Dia. 150 mm) |
| Overall dimension                  | 3870x840 mm  |
| Net Weight                         | 139 kg   |

## BHBD-818 HYDRAULIC TRANSFER STRETCHER



It is a type of hydraulic transfer stretcher  
 Five functions available by pressing the foot-pedals  
 Bright arrow symbols direct the users how to operate clearly  
 Central-lock castors, diameter 150mm. Matched with direction castor under the  
 chassis, which help users easily control the stretcher direction  
 Standard accessories: Mattress, I.V.Pole

### SPECIFICATIONS

|           |          |
|-----------|----------|
| Model     | BHBD-818 |
| Back-rest | 0 - 85°  |

|                                    |  |
|------------------------------------|--|
| Height range (bed board to ground) | 813 - 590 mm                                       |
| Standard accessory                 | Mattress, IV pole                                  |
| Bed Board                          | 1906x610 mm  |
| Castor                             | Central lockable 360° swivel castors (Dia. 150 mm) |
| Trendelenburg                      | 0 - 12°  |
| Reverse trendelenburg              | 0 - 12°  |
| Overall dimension                  | 1906x743 mm  |
| Net Weight                         | 82 kg  |

## BHBD-819 HYDRAULIC TRANSFER STRETCHER



It is a type of hydraulic transfer stretcher, 5 functions available by hydraulic and gas spring and crank system together

At the head-end, there is a pair of push- handles, can be folded up and down

Bright arrow symbols direct the users how to operate clearly

Central-lock castors, diameter 150mm. Matched with direction castor, help the operator controlling the direction easily, Bed board for X-Ray Examination

## SPECIFICATIONS

|                                    |  |
|------------------------------------|--|
| Model                              | BHBD-819   |
| Back-rest                          | 0 - 85°  |
| Height range (bed board to ground) | 860 - 630 mm                                       |
| Knee Gatch                         | 0° - 50°   |
| Trendelenburg                      | 0 - 13°  |
| Reverse trendelenburg              | 0 - 13°  |
| Standard accessory                 | IV pole  |
| Bed Board                          | 1932x630 mm  |
| Castor                             | Central lockable 360° swivel castors (Dia. 150 mm) |
| Overall dimension                  | 2050x745 mm  |
| Net Weight                         | 105 kg   |
| Height adjustment                  | 835-600 mm   |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)