



VACUUM CONCENTRATOR CENTRIFUGE

VACUUM CONCENTRATOR CENTRIFUGE BJZ-2601

CENTRIFUGAL VACUUM CONCENTRATOR



Built-in vacuum pump system.
Centrifugal concentration technology.
Compact and space-saving design.
Cost-effective laboratory solution.
Ideal for routine sample preparation.
Microcomputer control system, 7-inch LCD touch screen.
Delayed vacuum pump start function to avoid sample mixing.
Running parameters can be changed at will during operation.
Stainless steel corrosion-resistant inner liner, transparent acrylic lid.
A highly efficient magnetic drive system protects the airtightness of the cavity.
Chemical corrosion resistant plunger pump, standard low vacuum system.

SPECIFICATIONS

Model	BJZ-2601
Max. Speed	2000rpm
Max. RCF	430xg
Max. Capacity	4x50ml
Temperature setting range	room temperature +5°C to 80°C
Temperature accuracy	±1°C
Display method	LCD display
Overall machine noise	≤65dB (A)
Net weight	44kg
Ultimate vacuum	-92kpa
vacuum pump	built-in
Vacuum pump parameters	2L/min
Timing range	1 min-99h
Supported power supply	AC 220+22V 50Hz 16A
Dimension	580*360*330 (mm)
Alt Name	Centrifugal Vacuum Concentrator

ACCESSORIES FOR PURCHASE

No	name	Product Model	Rotor image
1	Fixed-angle rotor	66x1.5ml	
2	Fixed-angle rotor	24x1.5ml+24x5ml	

3	Fixed-angle rotor	12x10ml	
4	Fixed-angle rotor	4x50ml/4x100ml	
5	ELISA plate rotor/Microplate rotor 5	2x96-hole plate	

FEATURES

Built-in vacuum pump system.
 Centrifugal concentration technology.
 Compact and space-saving design.
 Cost-effective laboratory solution.
 Ideal for routine sample preparation.
 Microcomputer control system, 7-inch LCD touch screen.
 Delayed vacuum pump start function to avoid sample mixing.
 Running parameters can be changed at will during operation.
 Stainless steel corrosion-resistant inner liner, transparent acrylic lid.
 A highly efficient magnetic drive system protects the airtightness of the cavity.
 Chemical corrosion resistant plunger pump, standard low vacuum system.
 A variety of rotor sizes are available for selection, suitable for sample concentration of different capacities.
 The combined use of cold traps yields even better results.

APPLICATIONS

Suitable for concentration of protein, drugs, metabolites and enzymes. Used for solvent removal of DNA, RNA and nucleoside synthetic samples.
 DNA RNA sample concentration and purification.
 Protein sample preparation and purification.
 Laboratory sample processing and concentration.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
 Email: info@biolabscientific.com | Website: www.biolabscientific.com