



DENTAL AUTOCLAVE BADT-702

DENTAL AUTOCLAVE BADT-702

AUTOCLAVE (STEAM STERILIZER) CLASS B

Unique traceability system identifies load in different clinical and medical grade applications. Effectively sterilizes miscellaneous solid, hollow or porous (wrapped/unwrapped) materials and meets the requirement of today's demanding clinical environment. Automatic sterilization along with safety operating features and drying function makes it a perfect solution for effective and rapid steam sterilization. Innovated to comply with international standards, reliable solution for sterilization needs of reusable equipment.



- Fast and simple operation, Open water tank
- Open storage tank, waste water tank, easy to clean
- Concealment type interface
- Drain valve, output interface conceal
- Each interface front cover is open, easy to operate
- Chamber material is Stainless Steel (#304)
- Test the pot body tightness
- Built-in printer: The disinfection process date and data automatically printed
- USB: Built-in USB interface, End of the work, the information stored in the U disk automatically
- Automatic identification of fault protection system, Safe and fast judgment fault

SPECIFICATIONS

| | |
|---------------------------|--------------------------------------|
| Model | BADT-702 |
| Capacity | 23 L |
| Sterilization Temperature | 121°C, 134°C |
| Vacuum degree | -0.09Mpa |
| Steam leakage test | Yes |
| Residual humidity | <0.2% |
| 121°C Solid | Yes |
| 121°C Porous | Yes |
| 121°C Hollow | Yes |
| 134°C Porous | Yes |
| 134°C Hollow | Yes |
| B&D Test | Yes |
| Vacuum Test | Yes |
| Penetration test | Yes |
| Chamber Dimension | Ø249x450mm |
| Overall Dimension (LxWxH) | 471x505x640 mm |
| Display | 7" full-color touch LCD screen |
| Door | Automatic locking mechanism |
| Weight | 65 kg |
| Power | 2000 W |
| Power Supply | AC220V/110V±10% |
| 50/60Hz±1Hz | Printer |
| Alt Name | Autoclave (Steam Sterilizer) Class B |

FEATURES

- Fast and simple operation, Open water tank
- Open storage tank, waste water tank, easy to clean
- Concealment type interface
- Drain valve, output interface conceal

Each interface front cover is open, easy to operate
Chamber material is Stainless Steel (#304)
Test the pot body tightness
Built-in printer: The disinfection process date and data automatically printed
USB: Built-in USB interface, End of the work, the information stored in the U disk automatically
Automatic identification of fault protection system, Safe and fast judgment fault
Semi-automatic fault identifying system
3 pulsating vacuum
Suitable for sterilizing wrap, unwrap, solid, porous, hollow device
Use imported & advanced 16-bit microprocessor, easy to operate
The jet type of steam generator, ensure the efficient sterilization
With alarm system for waste water tank, avoid the waste water to enter into the sterilizing circulation
Advanced self-test system ensure you get the working data freely
The steam generator cleaning procedures, Clean the machine lines automatically every 30 times to avoid jams
Open water tank, Avoid growing chloride toxins
Alarm: Storage tank water shortage, waste water tank filled with water
Powerful cooling system, Effective condensing steam automatic electronic security door lock.

APPLICATIONS

Medical, Dental, Laboratory, Pharmaceutical, Research.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com