



## DISPENSING BOOTH BAIP-105

# DISPENSING BOOTH BAIP-105

## DISPENSING BOOTH (SAMPLING AND WEIGHING BOOTH)

Dispensing booth is a kind of partial purifying equipment for filling, refilling, weighing and sampling of raw material and compounds. It is provided with HEPA filter, which prevents the airborne dusts by down draught technique. The unit eliminates powder contamination to protect the operator and the surrounding environment. Dispensing booth is also called sampling booth or weighing booth.



1. Customized design is welcome.
  2. Unique designed air duct effectively controls the noise.
  3. Smooth transitions of wall and ground can eliminate the blind angle.
  4. Differential pressure gauge is equipped to real-time monitor the filters.
  5. Intelligent control mode and alarm system ensure the reliability of running.
  6. Automatic changeable frequency system monitors the air velocity to ensure its stability in the work area.
  7. Uniform flow design of main working area can protect the operators and prevent environmental pollution and cross contamination of products.
  8. Dispensing Booth(Sampling or Weighing Booth) has primary filters, medium efficiency filters and HEPA filters to keep air cleanliness of work area.
- Reference Picture

## SPECIFICATIONS

Model	BAIP-105
Clean Level	ISO 5 (Class 100), Class A
HEPA Filters	99.995% efficiency at 0.3µm
Air volume	Supply air volume: ≤7500m³/h, Exhaust air volume: ≤2250m³/h, Adjustable
Material	Fully stainless steel
Air Velocity	0.3~0.6m/s, adjustable
Pressure Gauge	3 pcs
Noise	≤75dB
UV Lamp	30W*2 Emission of 253.7 nanometers, with UV timer
Illuminating Lamp	16W*2&14W*4 LED lamp, Illumination: ≥500lux
Power Supply	AC220V±10%, 50/60Hz; 110V±10%, 60Hz
Consumption	1700W
Sockets	2pcs
Work Area Size(W*D*H)	2400*1800*2000mm
External Size(W*D*H)	2500*2400*2400mm
Gross Weight	840kg
Package Size(W*D*H)	2640*960*2380mm 2640*780*680mm 2640*560*2250mm
Alt Name	Dispensing Booth (Sampling and Weighing Booth)

## FEATURES

1. Customized design is welcome.
2. Unique designed air duct effectively controls the noise.
3. Smooth transitions of wall and ground can eliminate the blind angle.
4. Differential pressure gauge is equipped to real-time monitor the filters.
5. Intelligent control mode and alarm system ensure the reliability of running.
6. Automatic changeable frequency system monitors the air velocity to ensure its stability in the work area.
7. Uniform flow design of main working area can protect the operators and prevent environmental pollution and cross contamination of products.
8. Dispensing Booth(Sampling or Weighing Booth) has primary filters, medium efficiency filters and HEPA filters to keep air cleanliness of work area.





**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)