



LABORATORY VERTICAL AUTOCLAVE BAVT-3002

LABORATORY VERTICAL AUTOCLAVE BAVT-3002

VERTICAL AUTOCLAVE



- Hand wheel translation and quick opening structure
- High quality stainless steel material
- With drying function
- Microcomputer automatic control, arbitrary setting sterilization parameters
- Safety interlocking device
- Digital display and Touch buttons
- Standard test interface
- Self-expanding seal
- With automatic protection function over temperature protection; over pressure self discharge protection ; low water level protection; anti dry burning
- With stainless steel bucket or basket
- Buzzer reminder after sterilization, automatic stop

SPECIFICATIONS

| | |
|---|---------------------------|
| Model | BAVT-3002 |
| Capacity | 50 L |
| Power | 3.5 kW |
| Voltage | AC 220V. 50HZ |
| Design pressure | 0.25 Mpa |
| Design temperature | 139°C |
| Rated working pressure | 0.22 Mpa |
| Rated working temperature | 134°C |
| Sterilization temperature setting range | 116°C~134°C |
| Sterilization time setting range | 4 ~ 120 min |
| Drying time setting range | 0 ~ 240 min |
| Chamber volume | Ø 350x525 |
| Bucket dimension | Ø 330x460 |
| Basket dimension | Ø 320x240x2 |
| Outer dimension | 630x500x1180 |
| Packing size | 800x600x1230 |
| Net weight | 112 Kg |
| Gross Weight | 130 Kg |
| Alt Name | Vertical Steam Sterilizer |

FEATURES

- Hand wheel translation and quick opening structure
- High quality stainless steel material
- With drying function
- Microcomputer automatic control, arbitrary setting sterilization parameters
- Safety interlocking device
- Digital display and Touch buttons
- Standard test interface
- Self-expanding seal
- With automatic protection function over temperature protection; over pressure self discharge protection ; low water level protection; anti dry burning
- With stainless steel bucket or basket
- Buzzer reminder after sterilization, automatic stop

Automatic discharge of cold air, automatic steam exhaust after sterilization, Internal steam cycle



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com