



LABORATORY VERTICAL AUTOCLAVE BAVT-604-B

LABORATORY VERTICAL AUTOCLAVE BAVT-604-B

VERTICAL AUTOCLAVE



For the look, the only difference is the operation board:
Operation board (LED screen)


SPECIFICATIONS

| | |
|--------------------------|---|
| Model | BAVT-604-B |
| Type | Economical type |
| Capacity (L) | 100 |
| Dimension (LxWxH, mm) | 620x680x1200 |
| Chamber (Dia x H) | 400x809 |
| Chamber material | SUS304 |
| Power requirement | 220V±10%, 32A, 50/60Hz |
| Working environment | 5-45°C, RH 10-85% |
| Rate power supply (W) | 4600 |
| Sterilization temp. | 105-135°C |
| Temp. display accuracy | 0.1°C |
| Sterilization time range | 1-999min |
| Melting temp. | 60-100°C |
| Melting time range | 1-999min |
| Warming temp. | 45-60°C |
| Warming time range | 1-999min |
| Storage | up to 20 programs |
| Sterilization mode | 1.Solid 2.Agar with warming 3.Self defined |
| Exhaust temp. range | DWS: 73-127°C (Solid mode) |
| Exhaust mode | Automatic internal discharging with embedded steam collector |
| Pressure gauge range | -0.1-0.5Mpa |
| Pressure | Safety valve set pressure: 2.9bar |
| Controller | "Smart II" fast speed microcomputer controller |
| Safety device | Self-induction interlocking, lid closing checking, over temp protection, dry scorch protection, over pressure, safety valve, over current, short circuit, electric leakage, scald protection, auto failure inspection |
| Baskets (Dia x H) | (360x220)x3 |
| Alt Name | Vertical Autoclave |



FEATURES

For the look, the only difference is the operation board:
Operation board (LED screen)



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com