



LABORATORY LOW SPEED CENTRIFUGE BCFL-401

LABORATORY LOW SPEED CENTRIFUGE BCFL-401

FLOOR STANDING LARGE CAPACITY CENTRIFUGE



1. Widely used in the field of hospitals and laboratory.
2. Brushless frequency motor, in great torque, free maintenance, no powder pollution, quick in speed up and down.
3. Digital display which indicates the program, speed, time, RCF, temperature.
4. Micro computer control, there are 10 kinds of program and 10 kinds of acceleration and deceleration for your choice.
5. Electric lid lock, compact design, super speed and imbalance protection.
6. The centrifuge body is made of high-quality steel, safe and reliable.

SPECIFICATIONS

Model	BCFL-401
Max. Speed	5000 rpm
Max. RCF	4730 x g
Max. Capacity	4 x 800 ml
Timer	1~9h59min
RPM/RCF Convert	Yes
Noise (dB)	≤ 58
Temperature Range	/
Acc/Dec	10 Kinds
Speed Accuracy	±20 r/min
Temperature Accuracy	/
Voltage(V/Hz)	AC 220V/110V 50HZ/60HZ
Size (L x W x H mm)	720 x 540 x 830 mm
Net Weight(Kg)	108 KG
Certificates	CE, ISO & Calibration report are available
Alt Name	Floor Standing Large Capacity Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max Speed (r/min)	Volume (ml)	Max RCF (xg)
1	Swing Rotor	4000	4x800ml	3450
2	Swing Rotor	4000	4x750ml	3530
3	Swing Rotor	4000	4x500ml	3380
4	Swing Rotor	4000	4x4x96(well)	2940
5	Microplate rotor	4000	2x4x96(well)	3210
6	Swing rotor	5000	4x1x50ml	4730
7	Swing rotor	5000	4x1x100ml	4730
8	Swing rotor	4000	4x2x100ml	3020
9	Swing rotor	4000	4x2x50ml	3200
10	Swing rotor	4000	4x4x15ml	3200

11	Swing rotor	4000	4x6x15ml	3200
12	Swing rotor	4000	4x8x15ml	3200
13	Swing rotor	4000	4x12x5ml vacuum tube	2480
14	Swing rotor	4000	4x12x7ml vacuum tube	2760
15	Swing rotor	4000	4x12x10ml vacuum tube	2880
16	Swing rotor	4000	4x8x10/7/5ml vacuum tube; 4x6x10/7/5ml vacuum tube; 4x4x10/7/5ml vacuum tube	2790
17	Swing rotor	4000	4x30x7ml vacuum tube	3140
18	Swing rotor	4000	4x30x5ml vacuum tube	3140
19	Swing rotor	4000	4x40x7ml vacuum tube	3270
20	Swing rotor	4000	4x18x10ml vacuum tube	3140
21	Swing rotor	5000	24x15ml	3500
22	Swing rotor	5000	30x15ml	3830
23	Swing rotor	5000	12x50ml	3860
24	Swing rotor	5000	24x50ml	2970
25	Swing rotor	4000	12x100ml	2970



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com