



PASS BOX BEM1CL1 (BAPB-101)

PASS BOX BEM1CL1

PASS BOX



SPECIFICATIONS

Model	BEM1CL1
Old Model	BAPB-101
Type	Air Shower Pass Box
External Size(W*D*H)mm	660*570*1000
Internal Size(W*D*H)mm	500*500*500
Material	304Stainless steel
Interlock	Electronical
UV Lamp	15W*1
HEPA Filter	Efficiency≥99.995% at 0.3μm
Clean Air Speed	0.36~0.54m/s.
Power Supply	220V, 50Hz; 110V, 50/60Hz
Gross Weight (Kg)	80
Packing Size(W*D*H) mm	780*700*1250
Alt Name	Pass Box



APPLICATIONS

Pass Box is used to transfer goods between rooms of different clean levels of clean room. It can reduce times of opening door, prevent the air convection between different rooms and minimize the pollution extent.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com