



VISIBLE SPECTROPHOTOMETER BEQ-1701

VISIBLE SPECTROPHOTOMETER BEQ-1701

VIS SPECTROPHOTOMETER



7-Inch Color Touchscreen: Clear, intuitive operation for routine measurements.
30 Concentration Curves: Supports curve storage, linear regression, coefficient methods and direct concentration readout.

Automatic T/A & Zero Adjustment: Automatic T/A conversion plus 0%T/100%T adjustment reduces manual setup.

5-50 mm Cuvette Support: Accommodates different optical path lengths for various sample requirements.

USB Interface: Enables convenient data connection and processing.

Automatic Light Gate: Protects the photoelectric sensor without requiring a blackbody.

SPECIFICATIONS

Model	BEQ-1701
Wavelength Range	325-1100nm
Spectral Bandwidth	2nm
Wavelength Accuracy	±2nm
Photometry	Split Beam
Monochromator Type	Czerny-Turner
Focal Length	160mm
Grating	1200 lines/mm
Detector	Silicon Photocell
Wavelength Setting	Manual Turn Knob
Wavelength Repeatability	≤1nm
Stray Light	≤0.1%T(at 360nm NaNo ₂)
Photometric Range	0.0-200.0% T-0.097-4.000 A0.000-9999 C
Photometric Accuracy	±0.5% T±0.004 Abs (0-0.5A)±0.008 Abs (0.5-1A)
Photometric Repeatability	≤0.2% T0.002 Abs (0-0.5A)0.004 Abs (0.5-1A)
Noise	100%(T)noise≤0.3%(T) , 0%(T)noise≤0.2%(T)
Power	AC220V±22V 50Hz±1Hz, 120W
Packaging Size	580mm X 460mm X 345mm 0.1M3 16.5kg
Alt Name	VIS Spectrophotometer

ITEMS INCLUDED

No	Name	Quantity
1	Operation manual	1pc
2	UV WIN8 Spectrum data processing software	1pc
3	Glass cuvette 10mm	4pcs
4	Power cable	1pc
5	Fuse	2pcs
6	Cuvette holder 10mm	1pc

ACCESSORIES FOR PURCHASE

No	Name	Size
1	Cuvette holder	50mm
2	Cuvette	20mm 30mm 50mm

APPLICATIONS

Education & Teaching Laboratories

Widely used in teaching laboratories for demonstrating Beer-Lambert Law, calibration curves, and basic qualitative and quantitative analysis of colored solutions.

Chemical & Industrial Process Control

Measuring absorbance-based color intensity and tracking reaction kinetics in dyes, coatings, inks, and electroplating solutions for routine industrial quality control.

Pharmaceutical & Biotech Quality Control

Used for quantitative analysis of APIs and finished products with visible absorbance, dissolution testing, and optional basic DNA/protein quantification and reaction kinetics monitoring via PC software.

Water & Environmental Analysis

Used for routine colorimetric determination of ammonia nitrogen, total phosphorus, Cr, nitrite, and volatile phenols in environmental water quality monitoring.

Food & Agricultural Quality Control

For colorimetric detection of nitrite, nitrate, phosphate, peroxide value, and pigments in food and agricultural products for quality control compliance.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com