



SINGLE BEAM SCANNING UV VISIBLE SPECTROPHOTOMETER BEQ-3201

SINGLE BEAM UV VISIBLE SPECTROPHOTOMETER BEQ-3201

SINGLE BEAM UV/VIS SPECTROPHOTOMETER



Wide UV-Vis Range (190-1100 nm) for versatile quantitative and qualitative analysis.

Stable Single-Beam Optical System with holographic grating for accurate, low-noise measurements.

High Accuracy & Low Stray Light ensure reliable and repeatable laboratory results. User-Friendly Operation & Data Output with LCD display, USB/RS232 connectivity, and optional PC software.

High Photometric Stability & Low Noise improve analytical reliability.

PC & Printer Connectivity via USB and RS232 for easy data management.

Wide Wavelength Range (190-1100 nm) for UV and visible quantitative analysis.

Multiple Testing Modes including absorbance, transmittance, and quantitative ana

SPECIFICATIONS

Model	BEQ-3201
Wavelength Range	190-1100nm
Spectral Bandwidth	2.0nm or 4.0nm optional
Wavelength Accuracy	±0.8nm
Optical System	Single Beam, Blazed Holographic Grating (1200 lines/mm)
Wavelength Repeatability	≤0.2nm
Wavelength Setting	Auto, Resolution 0.1nm
Photometric Range	0-200%T, -0.3 - 3A, 0-9999C
Photometric Accuracy	±0.002A (0-0.5A), ±0.004A (0.5-1A), ±0.3%T (0-100%T)
Photometric Repeatability	≤0.001A (0-0.5A), ≤0.002A (0.5-1A), ≤0.2%T (0-100%T)
Stray Light	≤0.1%T(220/360nm)
Scan Speed	High, Medium, Low. Max.1000nm/minute
Baseline Flatness	±0.003A
Stability	±0.002A/h (500nm,0A)
Noise	≤0.2%T/3min (250/500nm,0%T); ≤0.5%T/3min (250/500nm,100%T)
Sample Compartment	Accommodate 5-100mm Pathlength Cuvette
Detector	Silicon Photodiode
Lamps	Tungsten Lamp & Deuterium Lamp (Pre-aligned)
Display	Graphic LCD (128x64 Dots)
Keypad	11-key Membrane Keypad
Output Port	USB Port & RS232 Port
Printer	Mini Serial Printer; PC Printer
PC Software	Optional PC Scanning Software
Power Requirements	AC 90-250V,50/60Hz
Dimension	431x318x195mm
Net Weight / Gross Weight	9kgs / 13kgs
Alt Name	Single Beam UV/VIS Spectrophotometer

ITEMS INCLUDED

No	Name	Quantity	Unit
1	10mm glass cuvette	4	pcs
2	10 mm Quartz Cuvette	2	pcs
3	Power line	1	pcs

ACCESSORIES FOR PURCHASE

No	Name	Quantity
1	Software card	1pcs
2	Software brochure	1pcs
3	Software dog	1pcs
4	USB cable	1pcs
5	PC scanning software	
6	Tungsten Lamp	
7	Deuterium lamp	
8	Micro cuvette cell holder (one-position)	
9	Manual 4-positions cell holder	
10	Automatic 4/8 positions cell holder	

FEATURES

Wide UV-Vis Range (190-1100 nm) for versatile quantitative and qualitative analysis.

Stable Single-Beam Optical System with holographic grating for accurate, low-noise measurements.

High Accuracy & Low Stray Light ensure reliable and repeatable laboratory results.

User-Friendly Operation & Data Output with LCD display, USB/RS232 connectivity, and optional PC software.

High Photometric Stability & Low Noise improve analytical reliability.

PC & Printer Connectivity via USB and RS232 for easy data management.

Wide Wavelength Range (190-1100 nm) for UV and visible quantitative analysis.

Multiple Testing Modes including absorbance, transmittance, and quantitative analysis.

APPLICATIONS

Education & Teaching Laboratories

Widely used in teaching laboratories for demonstrating Beer-Lambert Law, calibration curves, and basic qualitative and quantitative analysis of colored solutions.

Chemical & Industrial Process Control

Measuring absorbance-based color intensity and tracking reaction kinetics in dyes, coatings, inks, and electroplating solutions for routine industrial quality control.

Pharmaceutical & Biotech Quality Control

Used for quantitative analysis of APIs and finished products with visible absorbance, dissolution testing, and optional basic DNA/protein quantification and reaction kinetics monitoring via PC software.

Water & Environmental Analysis

Used for routine colorimetric determination of ammonia nitrogen, total phosphorus, Cr, nitrite, and volatile phenols in environmental water quality monitoring.

Food & Agricultural Quality Control

For colorimetric detection of nitrite, nitrate, phosphate, peroxide value, and pigments in food and agricultural products for quality control compliance.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com