



## ULTRASONIC CLEANER BER-5806

# ULTRASONIC CLEANER BER-5806

ULTRASONIC CLEANER DIGITAL MODEL, WITH TIMER AND HEATER



1. Use ordinary tap water, alcohol or solvent as cleaning;
2. High power transducer cleaning effect significantly, usable macroscopic observation;
3. After cleaning objects shine light;
4. LCD control boot time;
5. The stainless steel shell, bladder and cover, more upscale;
6. Waterproof properties is greatly improved, more safety products last;
7. constant temperature system with heating;
8. The stainless steel ultrasonic cleaning basket (optional).

## SPECIFICATIONS

Model	BER-5806
Frequency	40KHz
Timer	0-30 min
Heater	RT-80°C
Transducer(PCS)	6
Capacity	15L
Ultrasonic Power(W)	360
Heating Power(W)	450
Tank Size (L*W*H MM)	330*300*150
Unit Size (L*W*H MM)	530*325*285
Packing Size (L*W*H MM)	460*410*340
G.W./N.W	9.1/10.3
Alt Name	Ultrasonic Cleaner



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)