



ULTRASONIC CLEANER BER1CG3

ULTRASONIC CLEANER BER1CG3

ULTRASONIC CLEANER DIGITAL MODEL, WITH TIMER AND HEATER



1. Use ordinary tap water, alcohol or solvent as cleaning;
2. High power transducer cleaning effect significantly, usable macroscopic observation;
3. After cleaning objects shine light;
4. LCD control boot time;
5. The stainless steel shell, bladder and cover, more upscale;
6. Waterproof properties is greatly improved, more safety products last;
7. constant temperature system with heating;
8. The stainless steel ultrasonic cleaning basket (optional).

SPECIFICATIONS

Model	BER1CG3
Frequency	40KHz
Timer	0-30 min
Heater	RT-80°C
Transducer(PCS)	3
Capacity	4.5L
Ultrasonic Power(W)	180
Heating Power(W)	200
Tank Size (L*W*H MM)	300*150*100
Unit Size (L*W*H MM)	325*180*220
Packing Size (L*W*H MM)	410*255*310
G.W./N.W	4.6/5.5
Alt Name	Ultrasonic Cleaner



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com