





## **ULTRASONIC CLEANER BER1CG8**





## **ULTRASONIC CLEANER BER1CG8**

ULTRASONIC CLEANER DIGITAL MODEL, WITH TIMER AND HEATER



- 1. Use ordinary tap water, alcohol or solvent as cleaning;
- 2. High power transducer cleaning effect significantly, usable macroscopic observation;
- 3. After cleaning objects shine light;
- 4. LCD control boot time;
- 5. The stainless steel shell, bladder and cover, more upscale;
- 6. Waterproof properties is greatly improved, more safety products last;
- 7. constant temperature system with heating;
- 8. The stainless steel ultrasonic cleaning basket (optional).

## **SPECIFICATIONS**

Model	BER1CG8
Frequency	40KHz
Timer	0-30 min
Heater	RT-80°C
Transducer(PCS)	10
Capacity	30L
Ultrasonic Power(W)	600
Heating Power(W)	600
Tank Size (L*W*H MM)	500*300*200
Unit Size (L*W*H MM)	530*325*325
Packing Size (L*W*H MM)	635*405*390
G.W./N.W	14.4/16
Alt Name	Ultrasonic Cleaner









## Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada Email: info@biolabscientific.com | Website: www.biolabscientific.com