



## VORTEX MIXER BFA1BI1 (BTUR-403)

# VORTEX MIXER BFA1BI1

## ROTARY MIXING SOLUTION

Rotary Mixing Solution is mainly used for tissue culture, sample extraction, sedimentation rate determination. It is also widely applied in biochemistry, organism, cellular, immune, chemistry, flora and fauna analysis, clinical experiments, etc.




Humanized designed operation keyboard. LED display time and speed.  
Several optional turntable supplied, suit for the tubes from 1.5ml to 50ml.  
The turntables' gradient can be adjusted from 0 to 90°.  
Timing function. Time range from 0 to 9999min, LED display remaining time. Beep alarm when program time ends.

## SPECIFICATIONS

Model	BFA1BI1
Old Model	BTUR-403
Speed Range	10rpm~80rpm
Gradient Angle Range	0°~90°
Time Range	1min~9999min
Voltage	AC 100~230V, 50/60Hz
Power	35W
Fuse	250V, 1A, $\Phi$ 5x20
Dimension	W.206 x D.260 x H.300mm
Net Weight	4.5kgs
Alt Name	Rotary Mixing Solution

## ACCESSORIES FOR PURCHASE

No	Name	Description
1	Platform	Conical tube 50mlX10
2	Platform	Conical tube 15mlX18
3	Platform	1.5ml centrifuge tubeX64
4	Tube holder	50ml centrifuge tubeX12 15ml centrifuge tubeX12 1.5ml centrifuge tubeX30
5	Platform	50ml centrifuge tubeX10 15ml centrifuge tubeX10 1.5ml centrifuge tubeX30



1



2



3



4



5



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)