



VORTEX MIXER BFA1BI3 (BTUR-305)

VORTEX MIXER BFA1BI3

ROTARY MIXING SOLUTION

Rotary Mixing Solution is mainly used for tissue culture, sample extraction, sedimentation rate determination. It is also widely applied in biochemistry, organism, cellular, immune, chemistry, flora and fauna analysis, clinical experiments, etc.



It can be used in humid environments such as refrigerators and refrigerated cabinets.


Several optional turntables are supplied, suitable for the tubes from 1.5ml to 50ml. One button to start/stop, easy to operate, can run for a long time, running smoothly without noise.

SPECIFICATIONS

Model	BFA1BI3
Old Model	BTUR-305
Speed Range	20rpm
Gradient Angle Range	45° (fixed angle)
Time Range	/
Voltage	AC 220V/AC 110V, 50/60Hz
Power	14W
Fuse	250V, 1A, Φ 5x20
Dimension	W.260 x D.300 x H.233mm
Net Weight	1.8kgs
Alt Name	Rotary Mixing Solution

ACCESSORIES FOR PURCHASE

No	Name	Description
1	Platform	Conical tube 50mlX10
2	Platform	Conical tube 15mlX18
3	Platform	1.5ml centrifuge tubeX64
4	Tube holder	50ml centrifuge tubeX12 15ml centrifuge tubeX12 1.5ml centrifuge tubeX30
5	Platform	50ml centrifuge tubeX10 15ml centrifuge tubeX10 1.5ml centrifuge tubeX30



1



2



3



4



5



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com