



BENCHTOP HIGH SPEED CENTRIFUGE BFH-10501

HIGH SPEED CENTRIFUGE BFH-10501

TABLE TOP HIGH AND LOW SPEED UNIVERSAL CENTRIFUGE



1. Widely used for laboratory, chemistry and hospital.
2. 10 kinds of acceleration and deceleration.
3. More than 10 rotors for your choice.
4. High speed, low noise, compact design and small vibration
5. User-friendly LCD touch screen , display speed, time and temperature parameters. With 100 programs can be stored.
6. Brushless frequency motor which in great torque, no powder pollution, free maintenance, quick in speed up and down
7. Imbalance sensor and over speed sensor, the centrifuge body is made of high-quality steel, very safe.

SPECIFICATIONS

Model	BFH-10501
Max. Speed	18500rpm
Max. RCF	23990xg
Max. Capacity	4 x 250 ml
Time Range	0-99 h 59 min 59 sec
RPM/RCF Convert	Yes
Noise (dB)	≤54
Acc/Dcc	10 Kinds
Speed Accuracy	±20 r/min
Voltage (V/Hz)	AC 220V/110V, 50Hz/60Hz
Size (W x D x H)	570 x 400 x 320 mm
Net Weight (Kg)	40 kg
Certificates	CE, ISO & Calibration report are available
Alt Name	Table Top High and Low Speed Universal Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed (rpm)	Max. Volume (ml)	Max. RCF (xg)
1	Swing rotor	4000	4 x 250 ml (round cup)	2880
2	Swing rotor	4000	4 x 250 ml (square cup)	2830
3	Swing rotor	4000	4 x 13 x 10 ml	2830
4	Swing rotor	4000	4 x 16 x 7 ml	2860
5	Swing rotor	4000	4 x 16 x 5 ml	2470
6	Microplate rotor	4000	2 x 3 x 96 well	2270
7	Swing rotor	5000	4 x 1 x 50 ml	4730
8	Swing rotor	5000	4 x 1 x 100 ml	4730
9	Swing rotor	4000	4 x 2 x 100 ml	3000
10	Swing rotor	4000	4 x 2 x 50 ml	3020
11	Swing rotor	4000	4 x 4 x 15 ml	3020
12	Swing rotor	4000	4 x 6 x 15 ml	3020
13	Swing rotor	4000	4 x 4 x 10 ml	2770
14	Swing rotor	4000	4 x 4 x 7 ml vacuum tube adapter	2360
15	Swing rotor	4000	4 x 6 x 10 ml	2770
16	Swing rotor	4000	4 x 6 x 7 ml vacuum tube adapter	2360
17	Swing rotor	4000	4 x 8 x 10 ml	2770

18	Swing rotor	4000	4 x 8 x 7 ml vacuum tube adapter	2360
19	Swing rotor	4000	4 x 12 x 7 ml vacuum tube	2480
20	Swing rotor	4000	4 x 12 x 5 ml vacuum tube	2760
21	Swing rotor	4000	4 x 10 x 10 ml vacuum tube	2880
22	Fixed Rotor	14000	4 x 8PCR	12070
23	Fixed Rotor	13000	6 x 8PCR	16080
24	Fixed Rotor	14000	8 x 8PCR	13390
25	Fixed Rotor	13000	12 x 8PCR	17220
26	Fixed Rotor	14000	40 x 0.5 ml	19960
27	Fixed Rotor	18500	12 x 1.5 ml	23990
28	Fixed Rotor	15000	24 x 1.5 ml	20600
29	Fixed Rotor	13500	30 x 1.5 ml	19340
30	Fixed Rotor	11000	48 x 1.5 ml	12840
31	Fixed Rotor	15000	16 x 5 ml	19350
32	Fixed Rotor	14000	6 x 10 ml	16460
33	Fixed Rotor	13000	12 x 10 ml	17030
34	Fixed Rotor	11000	16 x 10 ml	13950
35	Fixed Rotor	11000	8 x 15 ml	12730
36	Fixed Rotor	9000	12 x 15 ml	9590
37	Fixed Rotor	9000	24 x 15 ml	3500
38	Fixed Rotor	7000	30 x 15 ml	3830
39	Fixed Rotor	12000	6 x 30 ml	14000
40	Fixed Rotor	12000	6 x 50 ml (conical)	16050
41	Fixed Rotor	12000	6 x 50 ml (round)	15960
42	Fixed Rotor	11000	4 x 85 ml	13560
43	Fixed Rotor	10000	4 x 100 ml	10310
44	Fixed Rotor	8000	6 x 100 ml	7280
45	Fixed Rotor	4000	12 x 100 ml	2970
46	Fixed Rotor	12000	24 capillaries	15800
47	Swing Rotor	15000	4 x 5 ml	19920
48	Vertical rotor	15000	16 x 5 ml	16450
49	Vertical rotor	14000	8 x 30 ml	19750



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com