



BENCHTOP HIGH SPEED CENTRIFUGE BFH-10701

HIGH SPEED CENTRIFUGE BFH-10701

TABLE TOP HIGH AND LOW SPEED UNIVERSAL CENTRIFUGE



1. Widely used for laboratory, chemistry and hospital.
2. 10 kinds of acceleration and deceleration
3. More than 10 rotors for your choice. 100 programs can be stored.
4. High speed, low noise, compact design and small vibration
5. User-friendly LCD display touch screen provides speed, time and temperature parameters.
6. Brushless frequency motor which in great torque, no powder pollution, free maintenance, quick in speed up and down
7. Imbalance sensor and over speed sensor, the centrifuge body is made of high-quality steel, very safe.
8. Rotor automatic identification is available for choice.

SPECIFICATIONS

Model	BFH-10701
Max. Speed	21000rpm
Max. RCF	30910xg
Max. Capacity	4 x 250 ml
Time Range	0-99 h 59 min 59 sec
RPM/RCF Convert	Yes
Noise (dB)	≤55
Acc/Dcc	10 Kinds
Speed Accuracy	±20 r/min
Voltage (V/Hz)	AC 220V/110V, 50Hz/60Hz
Size (W x D x H)	570 x 400 x 320 mm
Net Weight (Kg)	45 kg
Certificates	CE, ISO & Calibration report are available
Alt Name	Table Top High and Low Speed Universal Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max Speed (rpm)	Max Volume (ml)	Max. RCF (xg)
1	Swing rotor	4000	4 x 250 ml	2880
2	Swing rotor	4000	4 x 16 x 7 ml	2860
3	Swing rotor	4000	4 x 16 x 5 ml	2470
4	Swing rotor	4000	4 x 13 x 10 ml	2830
5	Microplate rotor	4000	2 x 3 x 96 Well	2310
6	Swing rotor	5000	4 x 1 x 50 ml	4730
7	Swing rotor	5000	4 x 1 x 100 ml	4730
8	Swing rotor	4000	4 x 2 x 100 ml	3020
9	Swing rotor	4000	4 x 2 x 50 ml	3020
10	Swing rotor	4000	4 x 4 x 15 ml	3020
11	Swing rotor	4000	4 x 6 x 15 ml	3020
12	Swing rotor	4000	4 x 4 x 10/7 ml	2770
13	Swing rotor	4000	4 x 4 x 5 ml vacuum tube adapter	2360
14	Swing rotor	4000	4 x 6 x 10/7 ml	2770
15	Swing rotor	4000	4 x 6 x 5 ml vacuum tube adapter	2360
16	Swing rotor	4000	4 x 8 x 10/7 ml	2770

17	Swing rotor	4000	4 x 8 x 5 ml vacuum tube adapter	2360
18	Swing rotor	4000	4 x 12 x 5 ml vacuum tube	2480
19	Swing rotor	4000	4 x 10 x 10 ml vacuum tube	2880
20	Fixed Rotor	16000	4 x 8PCR	15760
21	Fixed Rotor	15000	6 x 8PCR	21420
22	Fixed Rotor	16000	8 x 8PCR	17480
23	Fixed Rotor	15000	12 x 8PCR	22930
24	Fixed Rotor	15000	40 x 0.5 ml	22920
25	Fixed Rotor	21000	12 x 1.5 ml	30910
26	Fixed Rotor	16000	24 x 1.5 ml	23440
27	Fixed Rotor	14000	30 x 1.5 ml	20800
28	Fixed Rotor	13000	48 x 1.5 ml	17930
29	Fixed Rotor	16000	16 x 5 ml	22020
30	Fixed Rotor	16000	6 x 10 ml	21500
31	Fixed Rotor	15000	12 x 10 ml	22680
32	Fixed Rotor	13000	16 x 10 ml	19490
33	Fixed Rotor	13000	8 x 15 ml	17790
34	Fixed Rotor	10000	12 x 15 ml	11840
35	Fixed Rotor	5000	24 x 15 ml	3500
36	Fixed Rotor	5000	30 x 15 ml	3830
37	Fixed Rotor	14000	6 x 30 ml	19060
38	Fixed Rotor	13000	6 x 50 ml (falcon)	18840
39	Fixed Rotor	13000	6 x 50 ml (round)	18730
40	Fixed Rotor	5000	12 x 50 ml	3860
41	Fixed Rotor	13000	4 x 85 ml	18940
42	Fixed Rotor	12000	4 x 100 ml	14850
43	Fixed Rotor	10000	6 x 100 ml	11380
44	Fixed Rotor	12000	24 capillaries	15800
45	Swing Rotor	15000	4 x 5 ml	19920
46	Vertical Rotor	16000	16 x 5 ml	16450
47	Vertical Rotor	14000	8 x 30 ml	19750



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com