



DOUBLE BEAM SCANNING UV VISIBLE SPECTROPHOTOMETER BFP-201

DOUBLE BEAM UV VISIBLE SPECTROPHOTOMETER BFP-201

DOUBLE BEAM UV/VIS SPECTROPHOTOMETER



Wide wavelength range, satisfying requirements of various fields.
Five options for spectral bandwidth selection, 5nm, 4nm, 2nm, 1nm and 0.5nm, made according to customer's need and satisfying the requirements of pharmacopoeia.
Fully automated design, realizing easy measurement.
Optimized optics and large scale integrated circuits design, light source and receiver from world famous manufacturer all add up to high performance and reliability.
Rich measurement methods, wavelength scan, time scan, multi-wavelength determination, multi-order derivative determination, double-wavelength method and triple-wavelength method etc., meet different measurement requirements.
Automatic 10mm 8-cell holder, changeable to automatic 5mm-50mm 4-position cell holder for more choices.

SPECIFICATIONS

Model	BFP-201
Wavelength Range	190-1100nm
Spectral Bandwidth	2nm (5 nm, 4 nm, 1 nm, 0.5 nm optional)
Wavelength Accuracy	±0.3nm
Wavelength Reproducibility	≤0.15 nm
Photometric System	Double beam, auto scan, dual detectors
Photometric Accuracy	±0.3%T (0-100%T), ±0.002A (0-0.5A), ±0.004A (0.5A-1A)
Photometric Reproducibility	≤0.15%T
Working Mode	T, A, C, E
Photometric Range	-0.3-3.5A
Stray Light	≤0.05%T (NaI, 220 nm, NaNO ₂ 360 nm)
Baseline Flatness	±0.002A
Stability	≤0.001A/h (at 500 nm, after warming up)
Noise	±0.001A (at 500 nm, after warming up)
Display	6 inches high light blue LCD
Detector	Silicon photo-diode
Power	AC 220V/50Hz, 110V/60Hz, 180W
Dimensions	630 x 470 x 210 mm
Weight	26 kg
Alt Name	Double Beam UV/VIS Spectrophotometer

FEATURES

Wide wavelength range, satisfying requirements of various fields.
Five options for spectral bandwidth selection, 5nm, 4nm, 2nm, 1nm and 0.5nm, made according to customer's need and satisfying the requirements of pharmacopoeia.
Fully automated design, realizing easy measurement.
Optimized optics and large scale integrated circuits design, light source and receiver from world famous manufacturer all add up to high performance and reliability.
Rich measurement methods, wavelength scan, time scan, multi-wavelength determination, multi-order derivative determination, double-wavelength method and triple-wavelength method etc., meet different measurement requirements.
Automatic 10mm 8-cell holder, changeable to automatic 5mm-50mm 4-position cell holder for more choices.
Data output can be obtained via a printer port.
Parameters and data can be saved in case of power failure for user's convenience.
PC controlled measurement can be achieved via RS-232 interface (USB port) for more accurate and flexible requirements.

APPLICATIONS

This double beam UV/VIS spectrophotometer successfully realizes high accuracy and reliability measurement to meet various application requirements. It can be used extensively for qualitative and quantitative analysis in such fields as biochemical research and industry, pharmaceutical analysis and production, education, environmental protection, food industry, clinical examination, sanitation and antiepidemic etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com