



DRY BATH INCUBATOR BGH1D1 (BDIB-119)

DRY BATH INCUBATOR BGH1D1

DRY BATH INCUBATOR



Adopt innovative replaceable aluminum alloy block design; block automatic identification function.
Stepless adjustable thermo lid design, suitable for 0.2 mL, 0.5 mL, 1.5 mL, 2.0 mL centrifuge tubes, effectively preventing tube cover explosion and evaporation.
Customize the thermo lid switch status.
Easy to Exchange Blocks

SPECIFICATIONS

Model	BGH1D1
Old Model	BDIB-119
Temperature control range	RT. +5 °C~105 °C
Lid temperature range	RT. +5 °C~110 °C
Timer range	0~99 h 59 min / 0~99 min 59 sec
Temp. accuracy for blocks	±0.5 °C
Temp. accuracy for lid	±1 °C
Uniformity within the block	±0.5 °C
Display accuracy	0.1 °C
Heating time (RT. 25 °C)	≤15 min (25 °C-100 °C, @20 °C-30 °C)
Cooling time	/
Dimension (mm)	267x166x167 mm
Weight	2.3 kg
Alt Name	Dry Bath Incubator

ACCESSORIES FOR PURCHASE

No	Name	Series	Block	Description
1	Block	IB10HL	IB10-A	96x0.2 mL PCR tube
2	Block	IB10HL	IB10-B	54x0.5 mL tube
3	Block	IB10HL	IB10-C	35x1.5 mL tube
4	Block	IB10HL	IB10-D	35x2.0 mL tube
5	Block	IB10HL	IB10-E	15x0.5 mL+20x1.5 mL tube
6	Block	IB20HL	IB20-A	96x0.2 mL PCR tube
7	Block	IB20HL	IB20-B	54x0.5 mL tube
8	Block	IB20HL	IB20-C	35x1.5 mL tube
9	Block	IB20HL	IB20-D	35x2.0 mL tube
10	Block	IB20HL	IB20-E	15x0.5 mL+20x1.5 mL tube



1



2



3



4



5

FEATURES

Adopt innovative replaceable aluminum alloy block design; block automatic identification function.

Stepless adjustable thermo lid design, suitable for 0.2 mL, 0.5 mL, 1.5 mL, 2.0 mL centrifuge tubes, effectively preventing tube cover explosion and evaporation.

Customize the thermo lid switch status.

Easy to Exchange Blocks



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com