



DRY BATH INCUBATOR BGH1G4

DRY BATH INCUBATOR BGH1G4

DRY BATH INCUBATOR

Dry bath incubators are ideal for incubation and activation of cultures, enzyme reactions, immunoassays, melting / boiling points and a wide variety of other laboratory procedures.



Precise temperature control, fast and safe heating, temperature deviation calibration function.

Multiple block options, support block customization.

Convenient block replacement, easy to clean and sterilize.

Touch screen design, real-time display of temperature and time curves.

60 programs can be edited and stored, up to 4 sets of temperature points and time can be programmed.

Can be used for simple PCR experiments, with temperature and time cycling functions.

PCR thermo lid technology, dual blocks, with stainless plate.

SPECIFICATIONS

Model	BGH1G4
Temperature control range	RT. +5 °C~100 °C
Timer range	1 min-99 h 59 min
Temp. stability @100 °C	≤0.5 °C
Uniformity: within the block @37 °C / across similar block @37 °C	±0.2 °C / ±0.3 °C
Heating time (25 °C to 105 °C)	≤20 min
Display accuracy	0.1 °C
Block quantity	2
Dimension (LxWxH) mm	420x225x145 mm
Weight	5.5 kg
Alt Name	Dry Bath Incubator

ACCESSORIES FOR PURCHASE

No	Name	Type	Description
1	Block	BH01	96x0.2 mL PCR tube
2	Block	BH02	45x0.5 mL tube
3	Block	BH03	35x1.5 mL tube
4	Block	BH04	35x2.0 mL tube
5	Block	BH05	15x0.5 mL+20x1.5 mL tube
6	Block	BH06	20x1.5 mL+15x2.0 mL tube
7	Block	BH07	32x0.2 mL+22x0.5 mL+9x1.5 mL tube
8	Block	BH08	20x5 mL tube
9	Block	BH09	20x10 mL tube
10	Block	BH10	12x15 mL tube
11	Block	BH11	6x50 mL tube
12	Block	BH12	96x0.2 mL Elisa plate or 4 slides
13	Block	BH13	Solid block (no holes)
14	Block	BH14	384-well PCR plate
15	Block	BH15	40xΦ6 mm

16	Block	BH16	28xΦ10 mm
17	Block	BH17	24xΦ12 mm
18	Block	BH18	24xΦ13 mm
19	Block	BH19	14xΦ15 mm
20	Block	BH20	14xΦ16 mm
21	Block	BH21	12xΦ19 mm
22	Block	BH22	11xΦ20 mm
23	Block	BH23	6xΦ26 mm
24	Block	BH24	6xΦ28 mm
25	Block	BH25	2xΦ40 mm
26	Block	BH26	2x7-12.5x12.5 cuvette
27	Block	L4	Stainless plate



1



2



3



4



5



6



7



8



9



10



11



12



13



14



15



16



17



18



19



20



21



22



23



24



25



26



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com