



## DRYING OVEN BGH1I2 (BOHA-326)

## DRYING OVEN BGH112

### MICROPLATE INCUBATOR

They are commonly used for incubation of microplates and cell culture plates in microplate insect culture, animal and plant cell incubation, ELISA and other experimental processes.



Heating plates above and below the microplates, more even heating.  
Exact correspondence of the designed and actual temperature.  
LCD display. It is easy to set up and use.  
Real-time display of time and temperature, safe and reliable.

## SPECIFICATIONS

Model	BGH112
Old Model	BOHA-326
Temp. control range	RT. +5 °C~70 °C
Timer range	1 min~99 h 59 min
Temp. control accuracy	±0.5 °C
Display accuracy	±0.1 °C
Shaking speed	/
Mixing orbit	/
Heating time (20 to 70 °C)	≤20 min
Capacity	4 microplates / deepwell plates
Dimension (mm)	320x350x185 mm
Weight	8 kg
Alt Name	Microplate Incubator



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)