



THERMO SHAKER INCUBATOR BGI-4403

THERMO SHAKER INCUBATOR BGI-4403

THERMO-SHAKER



Microprocessor controller (with timing function)

- Temperature digital display and time count down display
- The double control segments of temperature and time
- Self-detection and alarm system
- Metal blocks for easy replacement, cleaning and disinfecting
- Over-temperature alarm protection
- Transparent cover to reduce the experiment error
- CD and film switch

SPECIFICATIONS

Model	BGI-4403
Temperature Range	-10°C~100°C
Time for temp. increasing	
Time for temp. decreasing	≤25 minutes (from 20°C to 100°C)
Temperature Accuracy	±0.5°C
Display Resolution	0.1°C
Uniformity of the block temperature	±0.5°C
Maximum Power	200W
Temperature-segment	2
Max. Time setting	5999 min
Standard	E-15X0.5ml+20X1.5ml
Packing Dimension	43*32*30
CBM	0.04
Net Weight	4.70
Gross weight	8
Alt Name	Thermo-shaker



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com