



# THERMO SHAKER INCUBATOR BGI-4501

# THERMO SHAKER INCUBATOR BGI-4501

THERMO-SHAKER



Microprocessor controller (with timing function)

- LCD display
- Various optional blocks, easy for changing and cleaning
- Available customizing blocks
- Microprocessor controller for temperature, time and shaking speed
- Excellent temperature accuracy and uniformity
- Simultaneously display setting and actual temperature, time and shaking speed
- Temperature calibrating function
- Over temperature protector
- Beep-signal and auto stop shaking after program completion

## SPECIFICATIONS

Model	BGI-4501
Temperature Range	RT+5°C ~ 100°C
Temperature Uniformity	≤±0.3°C
Temperature Stability	±0.2°C
Display Resolution	0.1°C
Temperature Increasing Time	< 20min (20°C ~100°C)
Shaking Speed	200-1500rpm
Amplitude	2mm
Timing Range	0-5999min
Electrical Requirement	AC220V 2.5A, 50/60Hz. Max, 150VA
Exterior Dimension(WxDxH)mm	300x212x180
Net weight	≤5.0kg
Remark	Standard Block:A-96X0.2ml B-54X0.5ml C-35X1.5ml D-35X2.0ml E-15X0.5ml+20x1.5ml
Alt Name	Thermo-shaker



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)