

## PRODUCT CATALOG



# STRECHER BGW1B1 (BHBD-817)





#### STRECHER BGW1B1

#### CONNECTING TRANSFER STRETCHER



- 1. Ideal for operation rooms and used to prevent cross-infection.
- 2. Single stretcher A is for sterile area (in the operation room), single stretcher B is for non-sterile area (outside the operation room). Adjust them to same heights (arrows at same level), push B to A, 2 stretchers will be connected and locked automatically. Then, the patient together with bed board C can slide from stretcher B to stretcher A safely along the rails, and VICE VERSA.
- 3. Central lockable 360° swivel castors (Dia. 150mm). Retractable 5th wheel provides effortless directional movement and cornering.

#### **SPECIFICATIONS**

Model	BGW1B1
Old Model	BHBD-817
Overall dimension (connected)	3870*840 mm (L*W)
Height range (bed board C to ground)	660 - 910 mm
Bed board C dimension	1906*610 mm (L*W)
Backrest	0 - 85°
Protective railings	Damping PP protective railings, lifting controlled by gas spring. Seamless connection with stretcher/operation table
Net weight	139 KG
Standard accessories	Mattress, I.V. Pole
Alt Name	Connecting Transfer Stretcher

#### **FEATURES**

- 1. Ideal for operation rooms and used to prevent cross-infection.
- 2. Single stretcher A is for sterile area (in the operation room), single stretcher B is for non-sterile area (outside the operation room). Adjust them to same heights (arrows at same level), push B to A, 2 stretchers will be connected and locked automatically. Then, the patient together with bed board C can slide from stretcher B to stretcher A safely along the rails, and VICE VERSA.
- 3. Central lockable 360° swivel castors (Dia. 150mm). Retractable 5th wheel provides effortless directional movement and cornering.











### Biolab Scientific Ltd.