



STRECHER BGW1D1 (BHBD-816)

STRECHER BGW1D1

MANUAL TRANSFER STRETCHER



1. Turn the crank to adjust height. Turn clockwise, the bed board will go up. Turn anticlockwise, the bed board will go down.
2. Clear arrow symbols to instruct the user to operate.
3. Central lockable 360° swivel castors (Dia. 150mm). Retractable 5th wheel provides effortless directional movement and cornering.
4. Integrated utility tray to store the patient's belongings and medical supplies.

SPECIFICATIONS

Model	BGW1D1
Old Model	BHBD-816
Overall dimension	1970*685 mm (L*W)
Height range (bed board to ground)	791 - 509 mm
Bed board dimension	1970*685 mm (L*W)
Backrest	0 - 89°
Side rails	Multifunctional rotary PP side rails, act as transfer board, fixable at horizontal level for IV/admin use
Bed board	Easy-to-clean PP bed board, integrally blow-molded
Net Weight	67 KG
Standard accessories	Mattress, I.V. Pole
Alt Name	Manual Transfer Stretcher

FEATURES

1. Turn the crank to adjust height. Turn clockwise, the bed board will go up. Turn anticlockwise, the bed board will go down.
2. Clear arrow symbols to instruct the user to operate.
3. Central lockable 360° swivel castors (Dia. 150mm). Retractable 5th wheel provides effortless directional movement and cornering.
4. Integrated utility tray to store the patient's belongings and medical supplies.





Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com