



TISSUE PROCESSOR BHTP-304

TISSUE PROCESSOR BHTP-304

SPIN TISSUE PROCESSOR (VACUUM)



1. Adopt Japan Mitsubishi PLC to control the whole working process, easy to operate, stable and reliable work
2. The starting cylinder can be set and has the function of customizing the cylinder position
3. Have optional function that can provide power -off protection (UPS).
4. Adopting agitation (vacuum) mode for dehydration.
5. Large LCD screen showing every step of processing and easy for monitoring, 160*90mm (7inch)
6. Timing : The configurable delay time is 30 days
7. Time setting range: It can be turned on any day or every day, the working time of each cylinder can be up to 999 minutes, the minimum is 1 minute for the 1st and 12th cylinders, and the minimum is 0 for the rest

SPECIFICATIONS

Model	BHTP-304
Carousel type	12 stations (9 reagent, 3 wax)
Basket loading	Two basket loading
Capacity	2.3 L reagent, 1.8 L wax parafin
Temperature of wax container	60-99°C
Material & capacity	Metal baskets up to 100 cassettes
Stirring interval	1 min every 5 min
Battery backup	4 hrs (optional)
Ergonomic control panel	Foil-protected keyboard & LCD
Infiltration programmable	Separate per station
Editing during run	Yes
Programmable infiltration	30 days per station, delayed 30 days
Programs selectable	10 selectable
Vacuum pressure	Max 0.05 MPa (0.5 bar)
Standard tissue basket size	Φ130x110 mm
Basket number & capacity	1 (2 optional), max 90 cassettes
Instrument dimension	Φ780 x H 570-690 mm
Nominal voltage	AC220V 50Hz / AC110V 50Hz
Nominal frequency	50Hz
Power	400-500 W
Net Weight	70 kg
Instrument Dimensions	780 x 570 x 690 mm
Package Info	
--Package	Wooden Case
--Package Dimensions	880 x 880 x 800 mm (machine), 620 x 410 x 530 mm (tanks)
--Gross Weight	110 kg (150 kg with bubble goods)
Alt Name	Spin Tissue Processor (Vacuum)

ITEMS INCLUDED

No	Name	Unit
1	Reagent bottle	9 pcs
2	Paraffin wax tank	3 pcs

3	Paraffin wax tank connection line	3 pcs
4	Tissue basket	1 pc
5	Tissue cassette	20 pcs
6	Hexagonal Wrench	1 pc
7	Screwdriver	1 pc
8	Fuse	2 pcs (5A)
9	Power line	1 line

ACCESSORIES FOR PURCHASE

No	Name
1	Extra paraffin tank
2	Extra reagent tank

FEATURES

1. Adopt Japan Mitsubishi PLC to control the whole working process, easy to operate, stable and reliable work
2. The starting cylinder can be set and has the function of customizing the cylinder position
3. Have optional function that can provide power -off protection (UPS).
4. Adopting agitation (vacuum) mode for dehydration.
5. Large LCD screen showing every step of processing and easy for monitoring, 160*90mm (7inch)
6. Timing : The configurable delay time is 30 days
7. Time setting range: It can be turned on any day or every day, the working time of each cylinder can be up to 999 minutes, the minimum is 1 minute for the 1st and 12th cylinders, and the minimum is 0 for the rest
8. Constant temperature controlling system adopted on the paraffin jar to ensure high precision of temperature controller . Internal surface oxidation treatment.
9. There are 10 programs selectable
10. It can be equipped with power supply according to user requirements, and can operate for 4-6 hours for other functions besides wax bath after power failure to ensure the safety of the organization
11. Anti-collision stopping device for hanging basket.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com