



## COOLED INCUBATOR BICL-403

# COOLED INCUBATOR BICL-403

COOLING INCUBATOR



Microprocessor control

- Polished stainless-steel chamber, semicircular arcs at corners for easy cleaning, and the space between the shelves in the chamber is adjustable.
- Latest PID controller
- R134a refrigerant, imported compressor
- Independent temperature-limiting alarm system ensures experiments run safely.(Option)
- Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature.(Option)
- There is a 25mm instruction connection hole on the left side of the chamber for easy testing operation and temperature measurement.

## SPECIFICATIONS

Model	BICL-403
Temperature Range	-40°C~65°C
Display Resolution	0.1°C
Temperature Stability	High Temp $\pm 0.5^{\circ}\text{C}$ / Low Temp $\pm 1.0^{\circ}\text{C}$
Electrical Requirement	380V 50Hz
Ambient Temperature	+5°C~35°C
Power Consumption	6100W
Interior Dimension (WxDxH)mm	600x500x830
Exterior Dimension (WxDxH)mm	800x1160x1850
Shelves	3(pcs)
Timing Range	0~9999min
Packing Dimension	205*93*130
CBM	2.48
Net Weight	282.00
Gross weight	313
Alt Name	Cooling Incubator

## FEATURES

Microprocessor control

- Polished stainless-steel chamber, semicircular arcs at corners for easy cleaning, and the space between the shelves in the chamber is adjustable.
- Latest PID controller
- R134a refrigerant, imported compressor
- Independent temperature-limiting alarm system ensures experiments run safely.(Option)
- Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature.(Option)
- There is a 25mm instruction connection hole on the left side of the chamber for easy testing operation and temperature measurement.

Option:

- Printer
- Independent temperature-limiting Alarm system
- RS485 connector

## APPLICATIONS

Serve for preserve culture medium, serum, medicine as well as microorganism training and environmental testing etc.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)