



VISIBLE SPECTROPHOTOMETER BJN-1101

VISIBLE SPECTROPHOTOMETER BJN-1101

VISIBLE SPECTROPHOTOMETER



Controlled by a single-chip microcomputer with digital display.
Supports manual zero calibration and 100% transmittance adjustment.
F factor can be directly input to build a standard curve for direct concentration reading.
Standard samples and their corresponding concentration values can be directly entered for concentration measurement.
Comes with a standard parallel printer interface for printing measurement data.
Built-in USB communication interface.

SPECIFICATIONS

Model	BJN-1101
Wavelength Range	325-1000nm
Spectral Bandwidth	4nm
Wavelength Accuracy	±2nm
Wavelength Repeatability	1.0 nm
Photometric Accuracy	±0.5% T
Photometric Repeatability	0.2% T
Stray Light	≤0.1% T
Photometric Mode	T, A, C
Wavelength Setting	Manual
Photometric Display Range	-0.3-3 A
Display Mode	LED display
Detector	Import Silicon Photodiode
Light Source	Tungsten Halogen Lamp
Power Requirement	AC 220V/50Hz, 110V/60Hz
Power	120 W
Dimensions (W x D x H)	470 x 325 x 180 mm
Weight	8 kg
Alt Name	Visible Spectrophotometer

ITEMS INCLUDED

No	Name	Quantity
1	1 cm glass cuvette	4 pcs
2	Dust cover	1 pc
3	Power line	1 pc
4	Instruction manual	1 pc



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com