



NANO SPECTROPHOTOMETER BJT-4401

NANO SPECTROPHOTOMETER BJN-4401

MICRO-VOLUME SPECTROPHOTOMETER



Microscale sample testing with a minimum sample volume of 0.5 μ L.
Wide detection range, up to 100 times that of a traditional spectrophotometer.
Most samples can be tested directly without dilution.
No warm-up required, enabling immediate testing.
Full wavelength coverage of 190-1100 nm with 1 nm accuracy; automatic wavelength scanning from 190-850 nm.
Compact and portable design, suitable for on-site testing.
More accurate and flexible testing is achieved through PC control.

SPECIFICATIONS

Model	BJN-4401
Wavelength Range	190-850nm
Spectral Bandwidth	
Wavelength Accuracy	1nm
Minimum Sample Size	0.5 μ L
Optical Distance	1.0 / 0.2 / 0.05 mm
Cuvette Optical Distance	10.0 / 5.0 / 2.0 / 1.0 mm
Wavelength Resolution	≤ 0.3 nm (FWHM at Hg 253.7 nm)
Absorbance Accuracy	2% (0.76 at 257 nm)
Absorbance Resolution	0.002 Abs (at 1 mm optical distance)
Absorbance Range	0.005-300 A (10 mm equivalent)
Detection Range	Lower limit: 2 ng/ μ L (dsDNA); Upper limit: 15,000 ng/ μ L (dsDNA)
Sample Pedestal Material	304 Stainless Steel and Quartz
Measurement Time	<5 s
Detector	3648 CCD
Light Source	Xenon Flash Lamp
Absorbance Test Range	0.01-0.5 Abs: ± 0.01 Abs; 0.5-80 Abs: $\pm 2\%$ (can extend to 300 Abs)
Dimensions (W x D x H)	200 x 130 x 136 mm
Weight	≤ 2.5 kg
Operating Power	24 W
Voltage	AC 100-240 V, 50-60 Hz
Display Language	English and Chinese
Alt Name	Micro-Volume spectrophotometer



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com