



DUCTLESS FUME HOOD BJU-3003

DUCTLESS FUME HOOD BJU-3003

DUCTED FUME HOOD

Fume Hood is designed as a negative pressure ventilation system that ensures the safety of operators and the laboratory environment. It works by drawing room air inward through the front opening of the hood and continuously ventilating it with an extraction blower. This prevents the spread of contaminated fumes, vapors, toxic gases, aerosols, and corrosive substances that may be generated during experiments. Furthermore, an optional active carbon filter can purify the polluted air before it is exhausted through the duct, thereby protecting the outdoor environment as well.



SPECIFICATIONS

Model	BJU-3003
Internal Dimensions (W*D*H)	1320*670*730
External Dimensions (W*D*H)	1540*800*2200
Maximum Opening	520mm
Work Surface Height	850mm
Display	Microprocessor controller with LCD display, display the operating status of each function
Airflow Velocity	0.3~0.8 m/s, speed adjustable
Exterior	Electro-Galvanized steel with antimicrobial powder coating
Interior	Compact laminate board with good acid and alkali resistance function
Work Surface	Solid physicochemical board, resistant to acid and alkali.
Front Window	Manual (FMH-DA series), Motorized (FMH-DA(II) series)
CabinetConstruction	5mm toughened glass; Height adjustable.
UV Lamp	30Wx1
ElectricalSystem	With timing function, Interlock with front sliding sash.
LED Lamp	9.6Wx1
Blower	Built-in centrifugal blower; Speed adjustable.
Consumption	500W
Power Supply	AC220V±10%, 50/60HZ; 110V±10%, 60Hz
Illumination	≥900Lux
Sound Emission	≤68dB(A)
Standard Accessory	1.Base Cabinet; 2.UV Lampx1; 3.LED Lampx1; 4.Waterproof Socketx2; 5.Water Sink; 6.Water Tap; 7.Gas Tap; 8.Gas & Water Remote Control Valve; 9.10m PVC aluminum foil composite exhaust duct, Pipe strapx2(Duct is retractable and can be cut freely)
Optional Accessory	1.Active carbon filter/HEPA filter; 2.Air Velocity Sensor.
Gross Weight (kg)	350
Shipping Size (WxDxH) (Main Body)	1690x990x1690
Base Cabinet	1690x990x900

Shipping Volume (m ³)	4.3
Alt Name	Ducted Fume Hood



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
 Email: info@biolabscientific.com | Website: www.biolabscientific.com