





# PRODUCT CATALOG



## BENCHTOP SHAKING INCUBATOR BJW-301

# BENCHTOP SHAKING INCUBATOR BJW-301

CELL CULTURE ROLLER



Stackable design: 1-3 units can be stacked to maximize usable capacity with minimal floor space.

Portable, all-in-one operation: integrated display, control panel, and unit body for convenient, space-saving use with easy parameter setting.

Parameter adjustment via keypad, allowing quick setting of speed, time, and other functions.

Accurate frequency display with long-term stable shaking performance.

## SPECIFICATIONS

Model	BJW-301
No-load shaking speed	2~50rpm
Shaking speed accuracy	±1rpm
Display mode	Non-touch LCD screen
Rolling bottle capacity	1000mlx3 or 2000mlx2
Timing range	0-999.9hours
Powersupply	AC220±10%, 50~60Hz
Overall Dimensions(L*W*H)	500mmx315mmx410mm
Net weight	30kg
Noise level	Below 65dB
Alt Name	Cell Culturing Roller

## MORE INFO

High-performance motor and belt drive ensure smooth, durable, and stable operation.

Compact, space-saving design with a streamlined stainless-steel aluminum exterior-lightweight, durable, and corrosion-resistant.

Provides steady, uniform mixing for cell suspensions.

Data protection features: automatic power-off in abnormal conditions and power-failure recovery to prevent data loss.

Premium exterior finish with professional anti-rust coating for enhanced durability and appearance.

Thoughtful safety details: anti-collision and anti-scratch design to protect users and culture flasks during operation.



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)