



BENCHTOP HIGH SPEED CENTRIFUGE BKA-1401

HIGH SPEED CENTRIFUGE BKA-1401

BENCHTOP HIGH CAPACITY CENTRIFUGE



1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.

SPECIFICATIONS

Model	BKA-1401
Max. Speed	6000rpm
Max. RCF	5080xg
Max. Capacity	4x500ml
Noise(dBA)	≤58dB(A)
Time range	1-99minutes
Speed accuracy	±20r/min
Dimension(mm)(W*D*H)	490*640*400mm
Net weight(without rotor)	60KGS
Power supply	AC220-240V,50/60HZ,10A,0.75KW
Alt Name	Benchtop High Capacity Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max.Speed (rpm)	Max.Capacity	Max. RCF (xg)	Note
1	Swing rotor	4000	4x500ml	3520	To order
2	Swing rotor	4000	96x2ml(75mm blood	2975	To order
3	Swing rotor	4000	96x5ml(100mm blood tube)	3405	To order
4	Swing rotor	4000	80x2ml(75mm blood	3260	To order
5	Swing rotor	4000	80x5ml(100mm blood	3620	To order
6	Swing rotor	4000	24x50ml	3405	To order
7	Swing rotor	4000	48x15ml	3000	To order
8	Swing rotor	4000	4x250ml(conical)	3640	To order

FEATURES

1. Microprocessor control with AC Brushless Motor
2. LCD Display the parameters of speed ,timer,RCF,rotors information,Program(Memory)information, Acceleration& Deceleration information and Error code information.
3. Steel centrifuge body and stainless steel chamber are autoclavable and resistant detergents and disinfection liquids.
4. Gas hinge for lid drop protection.
5. Safety Electrical lid interlock,door can not open when rotor is running and machine won't run when lid open and error code appears.
6. Emergency lid lock release when power failure.
7. Window on lid for speed calibration and noise measurement.
8. Continuous operation, working continuously according to the set time.
9. 10 kinds of Program(Memory) and 10 kinds of Acceleration&Deceleration for choice.If ask for more Programs please specially

mention when place the order.

10. Automatic troubleshooting with indication of error code and alarm on the screen.

11. Wide choice of Rotors and adapters.

12. With CE and ISO9001:2008, ISO13485 certificate.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com