



PRP CENTRIFUGE BKB-4401

PRP CENTRIFUGE BKB-4401

PRP CENTRIFUGE



The Plasma Centrifuge with microprocessor control, maintenance-free brush-less motor, and low noise are widely used in qualitative analysis of blood serum, plasma, and urea in hospitals, blood centers, and biochemistry.

RCF value can be calculated automatically. A separate short spin key for quick centrifugation.

Save up to 25 user-defined programs. Special programs can be saved for PRP/PRF/CGF extraction.

Soft-touch lid closure for ergonomic lid locking. The PRP centrifuge for sale cannot work if the door is open and the door cannot be opened when it is working.

9 acceleration and 10 deceleration rates, to achieve the best centrifugation effect.

SPECIFICATIONS

Model	BKB-4401
Max. Speed	4000rpm
Max. RCF	2760xg
Max. Capacity	4x50ml
Speed Accuracy	±30r/min
Time Setting Range	1min to 99min
Noise	≤65dB(A)
Power Supply	AC220V±22V 50/60Hz
Power Consumption	250W
Dimension (WxDxH)	380x500x250mm
Packing Dimension (WxDxH)	465x660x350mm
Weight	19kg
Alt Name	PRP Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Max. Speed	Max. RCF	Max. Capacity
1	Angle Rotor	4000	2200	4x50ml
2	Adapter	4000	2200	Upper(30ml Syringe)
3				Lower(30ml Syringe)
4				Upper(20ml Syringe)
5				Lower(20ml Syringe)

FEATURES

The Plasma Centrifuge with microprocessor control, maintenance-free brush-less motor, and low noise are widely used in qualitative analysis of blood serum, plasma, and urea in hospitals, blood centers, and biochemistry.

RCF value can be calculated automatically. A separate short spin key for quick centrifugation.

Save up to 25 user-defined programs. Special programs can be saved for PRP/PRF/CGF extraction.

Soft-touch lid closure for ergonomic lid locking. The PRP centrifuge for sale cannot work if the door is open and the door cannot be opened when it is working.

9 acceleration and 10 deceleration rates, to achieve the best centrifugation effect.

A special adapter for the syringe and stainless stopper reduces the possibility of pollution in the process of injection transplantation.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com