



BENCHTOP HIGH SPEED CENTRIFUGE BKE-1501

HIGH SPEED CENTRIFUGE BKE-1501

HIGH SPEED TABLE CENTRIFUGE



Intelligent control system

With five groups of optional programs and can store 99 program groups. Electric lid lock:

Long press the start/short-time key to realize instantaneous centrifugation.

Up to 10 level acceleration/ deceleration with Linear drive.

Excellent kernel performance

HD large LCD touch screen and the parameters can be revised during operation without shutdown.

Microcomputer control ,AC frequency conversion brush-less motor with fast speed and stable operation.

Two modes can be set: start timing/arrival revolution timing.

Careful and safe

SPECIFICATIONS

Model	BKE-1501
Max. Speed	16500rpm
Max. RCF	21615xg
Max. Capacity	12x5ml
Timer	1s~99h59min59s
Revolutions/min	±10r/min
Voltage	AC 220±22V 50Hz 10A or custom
Power	350W
Noise Level	≤ 62dB(A)
Size	270x390x240 (mm)
Net Weight	21kg
Alt Name	High Speed Table Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Capacity	Max Speed(r/min)	Max RCF(xg)
1	Angle rotor	16x1.5/2.0ml	16500rpm	20880xg
2		12x5ml	16000rpm	17344xg
3		24x1.5/2.0ml	15000rpm	21615xg
4		4x8x0.2ml PCR Row	15000rpm	14866xg
5		8x8x0.2ml PCR Row	14000rpm	15032xg
6		Capillary rotor 24 hole	12000rpm	13845xg

FEATURES

Intelligent control system

With five groups of optional programs and can store 99 program groups. Electric lid lock:

Long press the start/short-time key to realize instantaneous centrifugation.

Up to 10 level acceleration/ deceleration with Linear drive.

Excellent kernel performance

HD large LCD touch screen and the parameters can be revised during operation without shutdown.

Microcomputer control ,AC frequency conversion brush-less motor with fast speed and stable operation.

Two modes can be set: start timing/arrival revolution timing.

Careful and safe

Over speed,over heat,unbalance protections with auto lid lock.

Steel Structure,centrifuge chamber made of stainless steel with Environmentally friendly anti-corrosion coating treatment.

Intimate detail design

An intimate tool rack is provided for accompanying tools to guard the instrument at all times.

The front of the fuselage is designed with an operation SOP guide map, which makes it easier and more standard for laboratory personnel to use the instrument.

It is equipped with a rotor cover that can be quickly locked for easy use.

Special noise reduction system ;



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com