



FLOOR STANDING HIGH SPEED REFRIGERATED CENTRIFUGE BKE-3701

HIGH SPEED REFRIGERATED CENTRIFUGE BKE-3701

HIGH SPEED LARGE CAPACITY REFRIGERATED CENTRIFUGE



Superior features

Controlled by computer, AC frequency variable motor drive, toner contamination free and long life. Lower noise and less vibration, equipped with the most recent super light large capacity rotor to enhance efficiency significantly.

Memory-based programmed operation: ten programmed operations.

Easy to operate

LCD/LED display, programmable.

Optimal effect

User able to set up acceleration/deceleration timing at his will, automatic calculation of RCF value, heating and refrigerating dual control, optimal effect on centrifuge.

Safe and reliable

SPECIFICATIONS

Model	BKE-3701
Max. Speed	11500rpm
Max. RCF	18300xg
Max. Capacity	6x1000ml
timing range	1min~23h59min
temperature range	-20°C~40°C
temperature accuracy	±2.0°C
RPM accuracy	±20r/min
/	/
power supply	AC 220±22V,50Hz,30A
total power	4500W
noise level	< 65dB(A)
centrifuge chamber diameter	600mm
dimensions (length x width x height)	830x720x1250(mm)
packaging dimensions (length x width x height)	930x830x1350(mm)
Net Weight	230kg
Gross Weight	300kg
Alt Name	High Speed Large Capacity Refrigerated Centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Rotor Capacity	Maximum Speed	RCF
1	Angle rotor	6x250ml	10000rpm	18300xg
2		6x300ml	10000rpm	18300xg
3		6x500ml	8000rpm	11620xg
4		6x1000ml	7200rpm	12166xg
5	Horizontal Rotor	4x1000ml	4500rpm	5000xg
6		6x1000ml	4200rpm	5100xg

FEATURES

Superior features

Controlled by computer, AC frequency variable motor drive, toner contamination free and long life. Lower noise and less vibration, equipped with the most recent super light large capacity rotor to enhance efficiency significantly.

Memory-based programmed operation: ten programmed operations.

Easy to operate

LCD/LED display, programmable.

Optimal effect

User able to set up acceleration/deceleration timing at his will, automatic calculation of RCF value, heating and refrigerating dual control, optimal effect on centrifuge.

Safe and reliable

Equipped with over speed, overheat, imbalance and door lid auto lock protection systems to ensure safety of machine and users



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com