



FLOOR STANDING HIGH SPEED REFRIGERATED CENTRIFUGE BKE-3801

HIGH SPEED REFRIGERATED CENTRIFUGE BKE-3801

FLOOR HIGH SPEED LOW TEMPERATURE CENTRIFUGE



The series is stable but not heavy, excellent compatibility, large capacity with light rotor.

Microcomputer controlled, AC frequency variable motor drive, without toner power pollution, and longer lifetime. Lower noise and less vibration, equipped with multiple types of rotors for users to choose.

Auto recognition of rotor and speed limit control, to make centrifuge safer.

Imported compressor set, refrigeration and heating double temperature control, unique heat radiation system, to ensure that the cooling effect of the instrument to achieve the best.

Able to configure acceleration/deceleration mode, provide configuration of RPM and RCF, and convert RPM/RCF using one key.

Equipped with over speed, overheat, imbalance and door lid auto lock protection

SPECIFICATIONS

Model	BKE-3801
Max. Speed	25000rpm
Max. RCF	64800xg
Max. Capacity	4x1000ml
timing range	1min~23h59min
Temp Range	-20°C~40°C
Temp Accuracy	±2.0°C
RPM accuracy	±20r/min
/	/
power supply	AC 220±22V,50Hz,35A
total power	4500W
Noise level	≤ 65dB(A)
centrifuge chamber diameter	600mm
dimensions (length x width x height)	830x720x1250(mm)
packaging dimensions (length x width x height)	930x830x1350(mm)
Net Weight	320kg
Gross Weight	355kg
Alt Name	Floor High speed Low temperature centrifuge

ACCESSORIES FOR PURCHASE

No	Name	Rotor Capacity	Speed	RCF
1	Angle Rotor	24x1.5/2.2ml	25000rpm	64800xg
2		18x10ml	23000rpm	60460xg
3		8x50ml	21000rpm	52356xg
4		6x250ml	14000rpm	30074xg
5		8x100ml	13000rpm	19830xg
6		6x500ml	10000rpm	17700xg
7		4x1000ml	8000rpm	12040xg

FEATURES

The series is stable but not heavy, excellent compatibility, large capacity with light rotor.

Microcomputer controlled, AC frequency variable motor drive, without toner power pollution, and longer lifetime. Lower noise and less vibration, equipped with multiple types of rotors for users to choose.

Auto recognition of rotor and speed limit control, to make centrifuge safer.

Imported compressor set, refrigeration and heating double temperature control, unique heat radiation system, to ensure that the cooling effect of the instrument to achieve the best.

Able to configure acceleration/deceleration mode, provide configuration of RPM and RCF, and convert RPM/RCF using one key.

Equipped with over speed, overheat, imbalance and door lid auto lock protection systems to ensure safety of machine and users

Compact size, equipped with large torque variable-frequency motor to facilitate acceleration and deceleration and enhance efficiency. Higher performance to cost ratio.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com